

Activation and Expansion of Human T Cells With Immobilized CD3 and Soluble CD28 Antibodies

In a G-Rex[®] Bioreactor

Introduction

Adoptive T cell therapy is a rapidly growing segment of cell and gene therapy. This protocol outlines how to activate and expand purified human T cells using immobilized CD3 and soluble CD28 antibodies, all within a G-Rex bioreactor. T cells are stimulated with GMP CD3 and CD28 antibodies for 2 days followed by 7 days of expansion in Bio-Techne's xeno-free GMP T cell media. Robust CD3+ T cell expansion and high cell viabilities are achieved in the G-Rex.

Abbreviations:

BSC: Biological Safety Cabinet
 GMP: Good Manufacturing Practice
 PBMCs: Peripheral Blood Mononuclear Cells

TABLE // 01

Materials Required

| Material | Catalog Number |
|--|-------------------------|
| GMP T Cell Media | CCM038-GMP |
| Human IL-7 * | BT-007-GMP / BT-007-AFL |
| Human IL-15 * | BT-015-GMP / BT-015-AFL |
| Human CD3 GMP Antibody | MAB11411-GMP |
| Human CD28 GMP Antibody | MAB11412-GMP |
| G-Rex [®] 6 Well Plate M-Series | 80240M |
| Human AB Serum | |
| RectroNectin | Multiple vendors |
| 15 and 50 mL Centrifuge Tubes | |
| 1X PBS | |
| Cell Counter | |

*Note that IL-7 and IL-15 are available Animal-Free and GMP in lyophilized or liquid formulations

General Guidelines

- When handling bio-hazardous materials such as human cells, safe laboratory procedures should be followed.
- Maintain sterile technique by performing work in a BSC, wearing gloves, and using nuclease-free reagents and sterile pipette filter tips.
- All reagents should be stored according to the manufacturer's recommendations.
- If using the recommended Bio-Techne media and G-Rex bioreactor, cells will be confluent after 9 days of culture. If using a different media or reaction vessel, follow manufacturer's protocol to determine optimal culture time.

TABLE // 02

Timeline

| Vessel Coating | Activation | Transposition | Expansion |
|---|---|--|--|
| Day -1 <ul style="list-style-type: none"> • Coat vessel with RetroNectin and CD3 antibody (overnight) | Day 0 <ul style="list-style-type: none"> • Prepare media • Isolate/thaw T cells • Activate T Cells • Incubate (2 days) | Day 2 <ul style="list-style-type: none"> • Prepare media to dilute the culture • Mix and count cells • Incubate (7 days) | Day 6 (Optional) <ul style="list-style-type: none"> • Check cell growth Day 9 <ul style="list-style-type: none"> • Harvest cells for desired application |

Protocol

Vessel Preparation

Day -1 of Total Culture Time

Immobilization of CD3 Antibody to the Vessel Surface

1. Dilute CD3 antibody and RetroNectin in sterile 1X PBS in a 50 mL conical tube.
 - a. Recommended final concentration of CD3 antibody is 1 µg/mL and RetroNectin is 40 µg/mL.
 - b. These concentrations can be optimized for the end-user's workflow.
2. Add 2 mL of diluted CD3 antibody and RetroNectin to each well of G-Rex.
 - a. For different size G-Rex bioreactors, add 0.2 mL of diluted CD3 antibody and RetroNectin per cm² of G-Rex well.
3. Incubate G-Rex bioreactor overnight at 4 °C.

Activation

Day 0 of Total Culture Time

Prepare Media and Vessel

1. Transfer G-Rex bioreactor from 4 °C to room temperature to acclimate.

- Mix and sterile filter complete media.
 - Human GMP T Cell Media
 - 5% Human AB Serum
 - 10 ng/mL IL-7
 - 10 ng/mL IL-15
- Remove diluted CD3 antibody and RetroNectin from G-Rex bioreactor.
- Wash G-Rex bioreactor with sterile 1X PBS.
- Add at least 0.2 mL sterile 1X PBS per cm² of G-Rex well to prevent coated vessel surface from drying out.
 - Remove 1X PBS prior to adding T cells.

Prepare Purified T Cells

- Either:
 - Isolate T cells from human PBMCs following desired protocol or
 - Thaw purified T cells and wash with complete media.
- Aspirate supernatant.
- Resuspend cell pellet in complete media.
- Count resuspended cells.

Activate T Cells

- Seed 5 x 10⁶ purified T cells per well in G-Rex bioreactor.
 - Refer to the G-Rex Plating Reference below for seeding cell number of other G-Rex vessel sizes.
- Add CD28 antibody to final concentration of 1 µg/mL.
 - This concentration can be optimized for the end-user's workflow.
- Fill G-Rex bioreactor to activation volume of 10 mL with pre-warmed complete media (1 mL/cm²).
 - Refer to G-Rex Plating Reference to determine activation volume for different G-Rex bioreactor sizes.
- Transfer G-Rex bioreactor to a humidified incubator (37 °C, 5% CO₂) and incubate for 2 days.

Increase Culture Volume

Day 2 of Total Culture Time

Prepare Media

- Pre-warm desired volume of complete media needed to dilute cell culture.

Fill G-Rex Bioreactor

- Gently mix the activated cell complexes to break them apart, approximately 10-15 times.
 - If desired, cell counts can be performed at this step.
 - Note: We have observed disrupting cell complexes at this time improves cell expansion.
- Fill G-Rex bioreactor to maximum fill volume of 40 mL with pre-warmed complete media.
 - Refer to the G-Rex Plating Reference for fill volumes of different G-Rex vessel sizes.
- Place G-Rex bioreactor into humidified incubator at 37 °C, 5% CO₂.

Expansion

Optional: Day 6 of Total Culture Time

Prepare Media and Perform ½ Media Exchange

1. Pre-warm volume of complete media needed to complete a ½ media exchange in a G-Rex 6 well bioreactor.
 - a. Note: If using G-Rex 6M, 100M, or 500M, cells and media do not need handling and can be left untouched until day 9 of total culture time.
2. Transfer G-Rex bioreactor to BSC.
3. Carefully volume reduce wells of G-Rex to ~2 mL/cm².
4. Add pre-warmed complete media to wells to fill G-Rex bioreactor.
5. Proceed to check cell growth, or place G-Rex bioreactor into incubator at 37 °C, 5% CO₂.

Check Cell Growth

1. Cell counts and other analysis can be performed at this time.
2. Mix cells and sample each well for cell counts and desired analysis.
3. Place G-Rex bioreactor into incubator at 37 °C, 5% CO₂ until day 9 of total culture.

Day 9 of Total Culture Time

Harvest Cells

1. Reduce the volume of the cell culture media in the G-Rex bioreactor to ~1-2 mL/cm².
2. Mix the cells and sample each well for final cell counts and desired flow cytometry applications for phenotype characterization.
3. Cryopreserve remaining cells or use directly for desired applications.

TABLE // 03


G-Rex Plating Reference for Purified T Cells

| G-Rex Format | cm ² | Cell # for Activation | Activation Capacity | Maximum Capacity | Confluency |
|----------------------------|---------------------|-----------------------|---------------------|------------------|---------------------------|
| G-Rex 24 Well Plate | 2 cm ² | 1 x 10 ⁶ | 2 mL | 8 mL | 60-80 x 10 ⁶ |
| G-Rex 6 Well Plate | 10 cm ² | 5 x 10 ⁶ | 10 mL | 40 mL | 350-400 x 10 ⁶ |
| G-Rex 6M Well Plate | 10 cm ² | 5 x 10 ⁶ | 10 mL | 100 mL | 350-400 x 10 ⁶ |
| G-Rex 100M | 100 cm ² | 50 x 10 ⁶ | 100 mL | 1,000 mL | 3.5-4 x 10 ⁹ |
| G-Rex 500M | 500 cm ² | 250 x 10 ⁶ | 500 mL | 5,000 mL | 15-20 x 10 ⁹ |


Representative Data

FIGURE // 01


A. Expression of Activation Markers on Day 2



B. T Cell Count at Harvest



C. Cell Viability at Harvest



Immobilized CD3 and soluble CD28 GMP antibodies promote T cell expansion as well as competitor activation reagents. Purified human T cells from three independent donors were activated with GMP CD3 (Catalog # [MAB11411-GMP](#)) and CD28 (Catalog # [MAB11412-GMP](#)) antibodies for 2 days and expanded in GMP Human T Cell Media (Catalog # [CCM038-GMP](#)) supplemented with 5% human AB serum and 10 ng/mL of GMP IL-7 (Catalog # [BT-007-GMP](#)) and IL-15 (Catalog # [BT-015-GMP](#)) for 7 days in a 6M G-Rex bioreactor. After 48 hours of stimulation, levels of activation markers were analyzed by flow. Upregulation of CD25 and CD69 markers in CD4+ and CD8+ T cells indicate that cells were successfully activated with immobilized CD3 and soluble CD28 antibodies after two days of stimulation (A). Total cell numbers (B) and cell viabilities (C) were measured on day 9 of culture. Data is average \pm SD of n=3 donors.

For research use or manufacturing purposes only. Trademarks and registered trademarks are the property of their respective owners. 7679843487_1124

biotechne // Global Developer, Manufacturer, and Supplier of High-Quality Reagents, Analytical Instruments, and Precision Diagnostics.

INCLUDES R&D Systems™ Novus Biologicals™ Tocris Bioscience™ ProteinSimple™ ACD™ ExosomeDx™ Asuragen® Lunaphore™