PRODUCT INFORMATION & MANUAL

Human IL-1 β /IL-1F2 Valukine™ ELISA

Catalog Number: VAL101

For the quantitative determination of natural and recombinant
human Interleukin (IL)-1 β /IL-1F2 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202404.6

TABLE OF CONTENTS

I. BACKGROUND2

II. OVERVIEW 3

III. ADVANTAGES4

IV. EXPERIMENT 7

V. KIT COMPONENTS AND STORAGE 8

VI. PREPARATION 10

VII. ASSAY PROCEDURE 12

VIII. REFERENCES 13

I. BACKGROUND

The Interleukin-1 (IL-1) family of proteins consists of the classic members IL-1 α , IL-1 β , and IL-1ra, plus IL-18, IL-33 and IL-1F5-F10. IL-1 α and IL-1 β bind to the same cell surface receptors and share biological functions (1). IL-1 is not produced by unstimulated cells of healthy individuals with the exception of skin keratinocytes, some epithelial cells, and certain cells of the central nervous system. However, in response to inflammatory agents, infections, or microbial endotoxins, a dramatic increase in the production of IL-1 by macrophages and various other cell types is seen. IL-1 β plays a central role in immune and inflammatory responses, bone remodeling, fever, carbohydrate metabolism, and GH/IGF-I physiology. Inappropriate or prolonged production of IL-1 has been implicated in a variety of pathological conditions including sepsis, rheumatoid arthritis, inflammatory bowel disease, acute and chronic myelogenous leukemia, insulin-dependent diabetes mellitus, atherosclerosis, neuronal injury, and aging-related diseases (2-5).

IL-1 α and IL-1 β are structurally related polypeptides that show approximately 25% homology at the amino acid level. Both are synthesized as 31 kDa precursors that are subsequently cleaved into mature proteins of approximately 17.5 kDa (6, 7). Cleavage of the IL-1 β precursor by Caspase-1/ICE is a key step in the inflammatory response (2, 8). Neither IL-1 α nor IL-1 β contains a typical hydrophobic signal peptide (9-11), but evidence suggests that these factors can be secreted by non-classical pathways (12, 13). A portion of unprocessed IL-1 α can be presented on the cell membrane and may retain biological activity (14). The precursor form of IL-1 β , unlike the IL-1 α precursor, shows little or no biological activity in comparison to the processed form (13, 15). Both unprocessed and mature forms of IL-1 β are exported from the cell.

IL-1 α and IL-1 β exert their effects through immunoglobulin superfamily receptors that additionally bind IL-1ra. The 80 kDa transmembrane type I receptor (IL-1 RI) is expressed on T cells, fibroblasts, keratinocytes, endothelial cells, synovial lining cells, chondrocytes, and hepatocytes (16, 17). The 68 kDa transmembrane type II receptor (IL-1 RII) is expressed on B cells, neutrophils, and bone marrow cells (18). The two IL-1 receptor types show approximately 28% homology in their extracellular domains but differ significantly in that the type II receptor has a cytoplasmic domain of only 29 amino acids (aa), whereas the type I receptor has a 213 aa cytoplasmic domain. IL-1 RII does not appear to signal in response to IL-1 and may function as a decoy receptor that attenuates IL-1 function (19). The IL-1 receptor accessory protein (IL-1 RAcP) associates with IL-1 RI and is required for IL-1 RI signal transduction (20). IL-1ra is a secreted molecule that functions as a competitive inhibitor of IL-1 (21, 22). Soluble forms of both IL-1 RI and IL-1 RII have been detected in human plasma, synovial fluids, and the conditioned media of several human cell lines (23, 24). In addition, IL-1 binding proteins that resemble soluble IL-1 RII are encoded by vaccinia and cowpox viruses (25).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human IL-1 β has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human IL-1 β present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for human IL-1 β is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of IL-1 β bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ♦ **FOR RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ♦ This kit is suitable for cell culture supernate, human serum and human plasma.
- ♦ The kit should not be used beyond the expiration date on the kit label.
- ♦ Do not mix or substitute reagents with those from other lots or sources.
- ♦ If samples generate values higher than the highest standard, dilute the samples with Diluent (1 $\times$) and repeat the assay.
- ♦ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	20.7	65.0	141	21.4	66.9	137
Standard Deviation	1.20	4.0	5.3	2.3	6.1	9.3
CV%	5.8	6.2	3.8	10.5	9.2	6.8

B. RECOVERY

The recovery of human IL-1 β spiked to different levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range%
Cell culture media (n=4)	97	85 - 107
Serum (n=4)	89	84-92
Plasma (n=4)	94	93 - 97

C. SENSITIVITY

The minimum detectable dose (MDD) of IL-1 β is typically less than 1.0 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant human IL-1 β produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, samples were containing or spiked with high concentrations of human IL-1 β in various matrices and diluted with Diluent (1 $\times$) to produce samples with values within the dynamic range of the assay.

Dilution		Cell culture media (n=4)	Serum (n=4)	Plasma (n=4)
1:2	Average % of Expected	102	99	96
	Range (%)	98 - 107	98 - 101	92-100
1:4	Average % of Expected	103	102	101
	Range (%)	101 - 110	101 - 103	91-105
1:8	Average % of Expected	106	102	105
	Range (%)	103 - 110	101 - 102	103-110
1:16	Average % of Expected	107	105	104
	Range (%)	105 - 113	104 - 106	97-109

F. SAMPLE VALUES

Cell Culture Supernates - Human peripheral blood mononuclear cells (1 $\times$ 10⁶ cells/mL) were cultured in RPMI supplemented with 10% fetal calf serum, 50 μ M β -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate. Cells were stimulated as described below. Aliquots of the cell culture supernate was removed and assayed for levels of natural IL-1 β .

Stimulant	Day 1 (pg/mL)	Day 3 (pg/mL)	Day 5 (pg/mL)
10 μ g/mL PHA	2185	2004	2383
10 μ g/mL PHA + 10 ng/mL rhIL-2	1938	1973	2839
50 ng/mL PMA	1767	1027	1159
50 ng/mL LPS	4158	2145	1308

Human serum - Four serum samples were evaluated for the presence of human IL-1 β in this assay. All samples measured below the lowest standard, 3.9 pg/mL.

Human plasma - Four human plasma samples were evaluated for the presence of human IL-1 β in this assay. All samples measured ranged from 20.0 to 47.1 pg/mL with an average of 32.3 pg/mL.

G. SPECIFICITY

This assay recognizes both natural and recombinant human IL-1 β . The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhIL-1 β control were assayed for interference. No significant cross-reactivity or interference was observed.


Recombinant human			Recombinant mouse	
EGF	IL-2	IL-9	IL-1 α	IL-13
G-CSF	IL-2 sR α	IL-10	IL-1 β	TNF- α
GM-CSF	IL-3	IL-11	IL-3	
sgp130	IL-3 sR α	IL-12	IL-4	
GRO α	IL-4	IL-13	IL-5	
GRO β	IL-4 sR	IL-18	IL-5 sR α	
IFN- γ	IL-5	PDGF-AA	IL-6	
IGF-I	IL-5 sR β	PDGF-AB	IL-7	
IGF-II	IL-6	PDGF-BB	IL-9	
IL-1a	IL-6sR	VEGF	IL-10	
IL-1ra	IL-7			

Recombinant human IL-1 sRI and IL-1 sRII do not cross-react in this assay. However, interference was observed at concentrations greater than 10,000 pg/mL.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

pg/mL	OD	Average Corrected	
0	0.044 0.047	0.046	-
3.9	0.081 0.084	0.083	0.037
7.8	0.120 0.121	0.121	0.075
15.6	0.197 0.212	0.205	0.159
31.3	0.346 0.350	0.348	0.303
62.5	0.655 0.656	0.656	0.610
125	1.202 1.131	1.167	1.121
250	2.097 2.155	2.126	2.081

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human IL-1 β Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human IL-1 β	1 plate
Human IL-1 β Conjugate	Antibody against human IL-1 β conjugated to horseradish peroxidase	1 vial
Human IL-1 β Standard	Recombinant human IL-1 β in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Calibrator Diluent (5 $\times$)/RD5P	A 5 $\times$ concentrated buffered protein base used to dilute standard and samples.	1 vial
Wash Buffer Concentrate (25 $\times$)	A 25 $\times$ concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/ TMB Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Conjugate	
	TMB Substrate	
	Standard	Aliquot and store for up to 1 month at -20°C in a manual defrost freezer. * Avoid repeated freeze-thaw cycles.
	Calibrator Diluent (5×)/RD5P	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ♦ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ♦ Pipettes and pipette tips.
- ♦ Deionized or distilled water.
- ♦ Squirt bottle, manifold dispenser, or automated microplate washer.
- ♦ 500 mL graduated cylinder.

D. PRECAUTION

- ♦ The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.
- ♦ IL-1 β is detectable in saliva. Take the necessary precautions to prevent contamination of the kit reagents.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Diluent (1 $\times$).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Diluent (1 $\times$).

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Diluent (1 $\times$).

B. SAMPLE PREPARATION

Human serum samples recommend a 2-fold dilution. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Diluent (1 $\times$). Optimal dilutions should be determined by the end user.

Human plasma samples recommend a 2-fold dilution. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Diluent (1 $\times$). Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Note: Bring all reagents to room temperature before use.


Wash Buffer (1 $\times$) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 $\times$).

Diluent (1 $\times$) - Use deionized or distilled water to prepare Diluent (1 $\times$).

Human IL-1 β Standard - Refer to the vial label for the reconstitution volume*. This reconstitution produces a stock solution of 250 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 500 μ L of Diluent (1 $\times$) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 250 pg/mL standard serves as the high standard. The Diluent (1 $\times$) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of standard and prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 200 μ L of human IL-1 β Conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
6. Repeat the aspiration/wash as in step 4.
7. Add 200 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
8. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.
10. **CALCULATION OF RESULTS.**

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human IL-1 β concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Sims, J.E and D.E Smith (2010) Nat. Rev. Immuno. 10:89.
2. Martinon, F. and J. Tschopp (2007) Cell Death Differ. 14:10.
3. Isoda, K. and F. Ohsuzu (2006) J. Atheroscler. Thromb. 13:21.
4. Allan, S.M. et al. (2005) Nat. Rev. Immunol. 5:629.
5. Kornman, K.S. (2006) Am. J. Clin. Nutr. 83:475S.
6. Giri, J.G. et al. (1985) J. Immunol. 134:343.
7. Hazuda, D.J. et al. (1988) J. Biol. Chem. 265:6318.
8. Cerretti, D.P. et al. (1992) Science 256:97.
9. Lomedico, P.T. et al. (1984) Nature 312:458.
10. Auron, P.E. et al. (1987) J. Immunol. 138:1447.
11. March, C.J. et al. (1985) Nature 315:641.
12. Rubartelli, A. et al. (1990) EMBO J. 9:1503.
13. Rubartelli, A. et al. (1993) Cytokine 5:117.
14. Kurt-Jones, E.A. et al. (1985) Proc. Natl. Acad. Sci. 82:1204.
15. Hazuda, D. et al. (1989) J. Biol. Chem. 264:1689.
16. Urdal, D.L. et al. (1988) J. Biol. Chem. 263:2870.
17. Sims, J.E. et al. (1988) Science 241:585.
18. McMahan, C.J. et al. (1991) EMBO J. 10:2821.
19. Slack, J. et al. (1993) J. Biol. Chem. 268:2513.
20. Greenfeder, S.J. et al. (1995) J. Biol. Chem. 270:13757.
21. Eisenberg, S.P. et al. (1990) Nature 343:341.
22. Carter, D.B. et al. (1990) Nature. 344:633.
23. Dayer, J-M. and D. Burger (1994) Eur. Cytokine Netw. 5:563.
24. Svenson, M. et al. (1993) Cytokine 5:427.
25. Sims, J.E and S.K. Dower (1994) Eur. Cytokine Netw. 5:539.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

	A	B	C	D	E	F	G	H
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								

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