PRODUCT INFORMATION & MANUAL

Human IL-8/CXCL8 Valukine™ ELISA

Catalog Number: VAL103

For the quantitative determination of natural and recombinant
human IL-8/CXCL8 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202412.7

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	4
III. ADVANTAGES	5
IV. EXPERIMENT	8
V. KIT COMPONENTS AND STORAGE	9
VI. PREPARATION	11
VII. ASSAY PROCEDURE	13
VIII. REFERENCES	14

I. BACKGROUND

Interleukin-8 (IL-8; also known as GCP-1, NAP-1, and CXCL8) is a heparin-binding, 8-9 kDa member of the alpha, or CXC family of chemokines. Currently, there are 15 human CXC family members that generally range from 8-12 kDa in size. The majority are found on human chromosome 4, all contain a typical three β -sheet/one α -helix structure, and most show an N-terminal GluLeuArg (ELR) tri-peptide motif (1-3). Human IL-8 is synthesized as a 99 amino acid (aa) precursor that contains a 20 aa signal sequence plus a 79 aa mature region (4-6). It circulates as a monomer, homodimer, and heterodimer with CXCL4/PF4. The monomer is considered the most bioactive, while the heterodimer is proposed to potentiate PF4 activity (7-10). Mature human IL-8 shares 65% and 70% aa sequence identity with porcine and canine IL-8, respectively (11, 12). There is no IL-8 gene counterpart in rodents. Multiple isoforms of IL-8 exist that are generated through both alternative splicing and differential proteolytic processing. Alternative splicing occurs at the C-terminus where there is an 11 aa substitution for amino acids 92-99 (13). Proteolytic processing is likely a cell-specific event and results in IL-8 N-terminal truncation. For example, fibroblasts and endothelial cells cleave amino acids 21 and 22, generating IL-8 (1-7), while monocytes and lymphocytes cleave amino acids 21-25, generating IL-8 (6-7). These short forms generally show increased bioactivity, particularly towards the CXCR1 IL-8 receptor (6, 14). Approximately 15% of IL-8 also undergoes citrullination on Arg27 of the precursor. This creates an IL-8 with an increased half-life and ability to induce leukocytosis (15, 16). A wide variety of cells secrete IL-8, and these include monocytes and neutrophils (17), fibroblasts and keratinocytes (18), mast cells (19), visceral smooth muscle cells (20), dendritic cells (21), type II Greater alveolar cells (22), and endothelial cells (23).

There are two G-protein-coupled receptors for IL-8, termed CXCR1/IL-8RA and CXCR2/IL-8RB, that share 77% aa sequence identity (24). CXCR1 is 45-50 kDa in size and used almost exclusively by IL-8. CXCR2 is 35-40 kDa in size and used by almost all CXC chemokines (25, 26). Both CXCR1 and CXCR2 constitutively homodimerize, and this appears to be their functional configuration. When expressed on the same cell, heterodimerization will also occur, but this configuration undergoes disassembly following IL-8 binding (27). CXCR2 responds to low concentrations of IL-8 and is principally associated with chemotaxis and MMP-9 release. CXCR1, by contrast, responds to high concentrations of IL-8 and is associated with respiratory burst and

phospholipase D2 activation (26). Thus, and in the case of neutrophils, CXCR2 is suggested to initiate PMN migration to an inflammatory site followed by CXCR1 mediated priming for anti-microbial activity.

IL-8 is perhaps best known for its proinflammatory effects on immune cells. In essence, IL-8 is secreted by multiple cell types exposed to inflammatory stimuli. For monocytes/macrophages, microbial exposure results in IL-8 release. This is followed by CXCR2-mediated chemotaxis of PMNs to sites of antigenic challenge, and their subsequent activation/priming for anti-microbial activity. This activity is complemented by IL-8's ability to potentiate M-CSF action in the bone marrow that results in the release and maturation of granulocytes (14). IL-8 is also reported to have an angiogenic effect on tumor-associated endothelial cells. Here, tumor-derived IL-8 can act in a paracrine manner to activate CXCR1 and R2 on endothelial cells. Both receptors have been associated with PI3-K/Akt and RasGTP signaling, resulting in cell survival and proliferation. In addition, IL-8 positively regulates VEGF R2 and EGF R, two TKRs that mediate cell growth and migration (28).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human IL-8 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human IL-8 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for human IL-8 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human IL-8 bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ♦ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ♦ This kit is suitable for cell culture supernate, human serum and human plasma.
- ♦ The kit should not be used beyond the expiration date on the kit label.
- ♦ Do not mix or substitute reagents with those from other lots or sources.
- ♦ If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- ♦ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	166	414	708	164	410	751
Standard Deviation	6.3	18.3	33.2	9.4	23.4	33.4
CV%	3.8	4.4	4.7	5.7	5.7	4.4

B. RECOVERY

The recovery of Human IL-8 spiked to levels throughout the range of the assay in various matrices was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture media (n=4)	108	96 - 118
Serum (n=4)	83	79 - 88
Plasma (n=4)	96	84 - 101

C. SENSITIVITY

The minimum detectable dose (MDD) of human IL-8 is typically less than 7.8 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E. coli*-expressed recombinant human IL-8 produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, samples were containing or spiked with high concentrations of Human IL-8 in various matrices and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution		Cell culture media (n=4)	Serum (n=4)	Plasma (n=4)
1:2	Average % of Expected	114	103	92
	Range (%)	105-125	101-105	90-94
1:4	Average % of Expected	114	106	93
	Range (%)	105-122	105-108	91-95
1:8	Average % of Expected	112	114	92
	Range (%)	104-124	112-116	86-98
1:16	Average % of Expected	108	111	95
	Range (%)	102-107	107-115	90-102

F. SAMPLE VALUES

Cell Culture Supernates - Human peripheral blood mononuclear cells (1×10^6 cells/mL) were cultured in RPMI supplemented with 10% fetal calf serum, 50 μ M -mercaptoethanol, 2 mM L-glutamine, 100 U/mL penicillin, and 100 μ g/mL streptomycin sulfate and stimulated for 3 days with 10 μ g/mL PHA. An aliquot of the cell culture supernate was removed, assayed for levels of natural human IL-8, and measured 403,000 pg/mL.

Human serum - Four serum samples were evaluated for the presence of human IL-8 in this assay. All samples measured ranged from 92.8 to 100.6 pg/mL with an average of 96.9 pg/mL.

Human plasma - Four human plasma samples were evaluated for the presence of human IL-8 in this assay. All samples measured ranged from 137 to 226 pg/mL with an average of 171 pg/mL.

G. SPECIFICITY


This assay recognizes both natural and recombinant human IL-8. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range rhIL-8 control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human:	Recombinant mouse:
GRO α	KC
GRO β	MIP-1 α
GRO γ	MIP-1 β
I-309	
IP-10	
MCP-1	
MCP-2	
MCP-3	
MIP-1 α	
MIP-1 β	

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human IL-8 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human IL-8	1 plate
Human IL-8 Conjugate	Solution of antibody against human IL-8 conjugated to horseradish peroxidase	1 vial
Human IL-8 Standard	Recombinant human IL-8 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume	1 vial
Calibrator Diluent (5×)/RD5P	A 5× concentrated buffered protein base used to dilute standard and samples	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/ TMB Substrate Solution	2 vials
Stop Solution	2 N sulfuric acid	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8°C.*
	Stop Solution	
	Conjugate	
	TMB Substrate	
	Standard	Aliquot and store for up to 1 month at -20°C in a manual defrost freezer.* Avoid repeated freeze-thaw cycles.
	Calibrator Diluent (5×)/RD5P	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ♦ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ♦ Pipettes and pipette tips.
- ♦ Deionized or distilled water.
- ♦ Squirt bottle, manifold dispenser, or automated microplate washer.
- ♦ 500 mL graduated cylinder.

D. PRECAUTION

The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 $\times$).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 $\times$).

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000 x g within 30 minutes of collection. Assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1 $\times$).

B. SAMPLE PREPARATION

Human serum samples recommend a 2-fold dilution. A suggested 2-fold dilution is 100 μL of sample + 100 μL of Calibrator Diluent (1 $\times$). Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Bring all reagents to room temperature before use.

Note: High concentrations of human IL-8 are found in human saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.


Wash Buffer (1 $\times$) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25 $\times$) into deionized or distilled water to prepare 500 mL of Wash Buffer (1 $\times$).

Calibrator Diluent (1 $\times$) - Use deionized or distilled water to prepare Calibrator Diluent (1 $\times$).

Human IL-8 Standard - Refer to the vial label for the reconstitution volume*. This reconstitution produces a stock solution of 2000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 500 μL of Calibrator Diluent (1 $\times$) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 1000 pg/mL standard serves as the high standard. The Calibrator Diluent (1 $\times$) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

Note: *High concentrations of human IL-8 are found in human saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of standard and prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 μ L of human IL-8 Conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
6. Repeat the aspiration/wash as in step 4.
7. Add 200 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
8. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
9. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

10. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human IL-8 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Rosenkilde, M.M. and T.W. Schwartz (2004) *APMIS* 112:481.
2. Colobran, R. et al. (2007) *Clin. Exp. Immunol.* 148:208.
3. Lazennec, G. and A. Richmond (2010) *Trends Mol. Med.* 16:134.
4. Schmid, J. and C. Weissmann (1987) *J. Immunol.* 139:250.
5. Matsushima, K. et al. (1988) *J. Exp. Med.* 167:1883.
6. Van den Steen, P.E. et al. (2000) *Blood* 96:2673.
7. Rajarathnam, K. et al. (2006) *Biochemistry* 45:7882.
8. Burrows, S.D. et al. (1994) *Biochemistry* 33:12741.
9. Nesmelova, I.V. et al. (2005) *J. Biol. Chem.* 280:4948.
10. Nesmelova, I.V. et al. (2008) *J. Biol. Chem.* 283:24155.
11. Goodman, K.B. et al. (1992) *Biochemistry* 31:10483.
12. Ishikawa, J. et al. (1993) *Gene* 131:305.
13. Jang, J. and B. Kim (1998) GenBank Accession #:AAK00048.
14. Mortier, A. et al. (2008) *Pharmacol. Ther.* 120:197.
15. Loos, T. et al. (2009) *Haematologica* 94:1346.
16. Proost, P. et al. (2008) *J. Exp. Med.* 205:2085.
17. Smedman, C. et al. (2009) *J. Immunol. Methods* 346:1.
18. Lim, C.P. et al. (2009) *J. Invest. Dermatol.* 129:851.
19. Moller, A. et al. (1993) *J. Immunol.* 151:3261.
20. Vanaudenaerde, B.M. et al. (2003) *J. Heart Lung Transplant* 22:1280.
21. Sandor, N. et al. (2009) *Mol. Immunol.* 47:438.
22. Brasier, A.R. et al. (1998) *J. Biol. Chem.* 273:3551.
23. Watchorn, T.M. et al. (2002) *Am. J. Physiol. Endocrinol. Metab.* 282:E763.
24. Lee, J. et al. (1992) *J. Biol. Chem.* 267:16283.
25. Nasser, M.W. et al. (2009) *J. Immunol.* 183:3425.
26. Stillie, R. et al. (2009) *J. Leukoc. Biol.* 86:529.
27. Monoz, L.M. et al. (2009) *J. Immunol.* 183:7337.
28. Waugh, D.J.J. and C. Wilson (2008) *Clin. Cancer Res.* 14:6735.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

12								
11								
10								
9								
8								
7								
6								
5								
4								
3								
2								
1								
	A	B	C	D	E	F	G	H

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