PRODUCT INFORMATION & MANUAL

Human IL-12 p70 Valukine™ ELISA

Catalog Number: VAL115

For the quantitative determination of natural and recombinant human
Interleukin (IL)-12 p70 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202503.4

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	3
III. ADVANTAGES	4
IV. EXPERIMENT	6
V. KIT COMPONENTS AND STORAGE	7
VI. PREPARATION	9
VII. ASSAY PROCEDURE	11
VIII. REFERENCES	12

I. BACKGROUND

Interleukin 12 p70 (IL-12), also known as natural killer cell stimulatory factor (NKSF) or cytotoxic lymphocyte maturation factor (CLMF), is a heterodimeric pleiotropic cytokine made up of a 40 kDa (p40) subunit and a 35 kDa (p35) subunit (1,2-4). Human p40 will circulate either as a monomer, homodimer, or in a complex bound to either p35, forming IL-12, or to p19, forming IL-23 (2-4, 5). Both the homodimeric p40, and IL-23 can bind to the IL-12 R, serving as nonsignaling antagonists (2, 6, 7). Alternatively, the p40 homodimer may also bind to R β 1, activating microglia and macrophages (3, 8). The IL-12 p40 subunit is shared by IL-23, another heterodimeric cytokine that has biological activities similar to, as well as distinct from, IL-12. Cells known to produce IL-12 include macrophages and dendritic cells (9), monocytes (10), Langerhans cells (11), neutrophils (12), keratinocytes (13), plasmacytoid dendritic cells (14), microglia (4), CD8⁺ DC (mouse cells only) (15) and non-germinal center (CD38-CD44⁺) B cells (human cells only) (2, 16). Functionally, IL-12 has been shown to both enhance cytotoxic activity and induce interferon gamma (IFN- γ) production in NK cells, T cells and dendritic epidermal T cells (2, 17-19). IL-12 has also been reported to induce IFN- γ production in macrophages (29). IL-12, in conjunction with the other IL-12 family members IL-23 and IL-27, promotes the development of a CD4⁺ Th1 immune response (3, 4, 20). In response to infection, IL-27 is released initially, promoting a Th0 to Th0/1 transition. IL-12 production follows, generating Th1 effector cells. In combination with IL-18, IL-12 creates Th1 memory cells out of effector cells, and these cells are later activated by IL-23 (3). While mouse IL-12 is active on both human and mouse cells, human IL-12 is not active on mouse cells.

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human IL-12 p70 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any human IL-12 p70 present is bound by the immobilized antibody. After washing away any unbound substances, a biotinylated detection antibody specific for human IL-12 p70 is pipetted into the wells. After washing away any unbound substances, Streptavidin-HRP is added. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of human IL-12 p70 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and human serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	49.9	256	1333.2	53.4	282.9	1484.9
Standard Deviation	2.9	11.6	79.2	4.4	23.4	135.5
CV%	5.7	4.5	5.9	8.2	8.3	9.1

B. RECOVERY

The recovery of human IL-12 p70 spiked to three levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 116.0 to 119.0% with an average of 117.6%.

The recovery of human IL-12 p70 spiked to three levels throughout the range of the assay in human serum was evaluated. The recovery ranged from 97.9 to 105.9% with an average of 101.3%.

C. SENSITIVITY

The minimum detectable dose (MDD) of human IL-12 p70 is typically less than 0.7 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified Sf 21-expressed recombinant human IL-12 p70 produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were spiked with high concentrations of human IL-12 p70 and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	91.8	83.8-98.3
1:4	94.6	81.9-104.4
1:8	100.9	86.6-111.7
1:16	106.9	86.3-117.7

F. SAMPLE VALUES

Human serum - Forty samples from apparently healthy volunteers were evaluated for the presence of human IL-12 p70 in this assay. No medical histories were available for the donors used in this study. All samples measured less than the lowest human IL-12 p70 standard, 31.3 pg/mL.

G. SPECIFICITY


This assay recognizes both natural and recombinant human IL-12 p70. The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity or interference.

Recombinant human:	Recombinant mouse:
IL-12 p35	IL-12 p70
IL-12/IL-23 p40	
IL-23	

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

pg/mL	O.D.		Average	Corrected
0	0.115	0.116	0.116	—
31.3	0.22	0.214	0.217	0.101
62.5	0.315	0.296	0.305	0.189
125	0.496	0.498	0.497	0.381
250	0.812	0.818	0.815	0.699
500	1.278	1.267	1.272	1.156
1000	1.893	1.885	1.889	1.773
2000	2.468	2.439	2.453	2.337

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Human IL-12 p70 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against human IL-12 p70.	1 plate
Human IL-12 p70 Standard	Recombinant human IL-12 p70 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Human IL-12 p70 Detection Antibody	Biotinylated human IL-12 p70 Antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Calibrator Diluent Concentrate (2×)	A 2× concentrated buffered diluent used to dilute standard and samples.	1 vial
Streptavidin-HRP B (40×)	40× Streptavidin conjugated to horseradish peroxidase.	1 vial
Reagent Diluent Concentrate (10×)	A 10× concentrated buffered protein base used to dilute Detection Antibody and HRP.	1 vial
Normal goat serum	Diluted, heat-inactive normal goat serum.	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant.	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution.	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strip.	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Streptavidin-HRP B	May be stored for up to 1 month at 2-8°C.*
	Normal goat serum	
	Wash Buffer (1×)	
	TMB Substrate	
	Stop Solution	
	Standard	Prepare fresh for each assay. Standards may be stored for up to 1 month at -20 °C.*
	Detection Antibody	Aliquot and store for up to 1 month at -20°C in a manual defrost freezer. *
	Calibrator Diluent Concentrate (2×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Reagent diluent Concentrate (10×)	May be stored for up to 1 month at 2-8°C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ♦ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ♦ Pipettes and pipette tips.
- ♦ Deionized or distilled water.
- ♦ Squirt bottle, manifold dispenser, or automated microplate washer.
- ♦ Test tubes for dilution of standards.
- ♦ 100 mL and 500 mL graduated cylinder.

D. PRECAUTION

- ♦ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ♦ The Stop Solution provided with this kit is an acid solution.
- ♦ Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 × g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

B. REAGENT PREPARATION

Bring all reagents to room temperature before use.

Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Reagent Diluent (1×) - Use deionized or distilled water to prepare Reagent Diluent (1×).

Calibrator Diluent (1×) - use deionized or distilled water to prepare Calibrator Diluent (1×).


Detection Antibody (1×) - **Centrifuge briefly before opening. Reconstitution volume refer to vial label to prepare Detection Antibody (100×).** Detection Antibody to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Aliquot and store if needed. Dilute stock solution in Reagent Diluent (1×) with 2% heat inactivated normal goat serum (NGS) (200 µL NGS per 10 mL solution) to Detection Antibody (1×). Prepare at least 15 minutes prior to use.

Streptavidin-HRP B (1×) - **Centrifuge briefly before opening.** Dilute to the working concentration specified on the vial label with Reagent Diluent (1×).

Human IL-12 p70 Standard - **Centrifuge briefly before opening. Refer to the vial label for the reconstitution volume*.** This reconstitution produces a stock solution of 4000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 500 µL of Calibrator Diluent (1×) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 2000 pg/mL standard serves as the high standard. The **Calibrator Diluent (1×)** serves as the zero standard (0 pg/mL).

C. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of Standard and sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 μ L of the Detection Antibody (1 $\times$) diluted in Reagent Diluent (1 $\times$) with NGS, to each well. Cover with a new adhesive strip and **incubate for 2 hours at room temperature.**
6. Repeat the aspiration/wash as in step 4.
7. Add 100 μ L of the working dilution of Streptavidin-HRP B to each well. Cover the plate and **incubate for 20 minutes at room temperature. Protect from light.**
8. Repeat the aspiration/wash as in step 4.
9. Add 100 μ L of TMB Substrate to each well. **Incubate for 20 minutes at room temperature. Protect from light.**
10. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

12. CALCULATION OF RESULTS.

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the human IL-12 p70 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Schoenhaut, D.S. et al. (1992) *J. Immunol.* 148:3433.
2. Hamza, T. et al. (2010) *Int. J. Mol. Sci.* 11:789.
3. Brombacher, F. et al. (2003) *Trends Immunol.* 24:207.
4. Trinchieri, G. (2003) *Nat. Rev. Immunol.* 3:133.
5. Gee, K. et al. (2009) *Inflamm. Allergy Drug Targets* 8:40.
6. Heinzl, F.P. et al. (1997) *J. Immunol.* 158:4381.
7. Sieve, A.N. et al. (2010) *Eur. J. Immunol.* 40:2236.
8. Jana, M. et al. (2009) *Glia* 57:1553.
9. Kato, T. et al. (1997) *Cell. Immunol.* 181:59.
10. Blotta, M.H. et al. (1997) *J. Immunol.* 158:5589.
11. Kang, K. et al. (1996) *J. Immunol.* 156:1402.
12. Romani, L. et al. (1997) *J. Immunol.* 158:5349.
13. Yawalkar, N. et al. (1996) *J. Invest. Dermatol.* 106:80.
14. Krug, A. et al. (2001) *Eur. J. Immunol.* 31:3026.
15. Shortman, K. and W. Heath (2010) *Immunol. Rev.* 234:18.
16. Schultze, J.L. et al. (1999) *J. Exp. Med.* 189:1.
17. Novelli, F. and J.L. Casanova (2004) *Cytokine Growth Factor Rev.* 15:367.
18. Sugaya, M. et al. (1999) *J. Invest. Dermatol.* 113:350.
19. Tominaga, K. et al. (2000) *Int. Immunol.* 12:151.
20. Collison, L.W. and D. Vignali (2008) *Immunol. Rev.* 226:248.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1												
2												
3												
4												
5												
6												
7												
8												
9												
10												
11												
12												
	A	B	C	D	E	F	G	H				

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