PRODUCT INFORMATION & MANUAL

Mouse IL-17 Valukine™ ELISA

Catalog Number: VAL610

For the quantitative determination of natural and recombinant
mouse Interleukin IL-17 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202408.5

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	3
III. ADVANTAGES	4
IV. EXPERIMENT	7
V. KIT COMPONENTS AND STORAGE	8
VI. PREPARATION	10
VII. ASSAY PROCEDURE	12
VIII. REFERENCES	14

I. BACKGROUND

Mouse Interleukin 17 (IL-17; also IL-17A and CTLA-8) is a 21 kDa, variably glycosylated polypeptide that belongs to the IL-17 family of cytokines containing a cysteine-knot fold (1-3). Its sequence was originally isolated from an activated hybridoma created from the fusion of a mouse cytotoxic and rat T cell lymphoma cell line (2-5). It is synthesized as a 158 amino acid (aa) precursor that contains a 25 aa signal sequence and a 15 kDa, 133 aa mature segment (5). In both mouse and human, there is one conserved N-linked glycosylation site that likely contributes 5 kDa to its native molecular weight. IL-17A forms both a 35-38 kDa homodimer, and a 45-48 kDa heterodimer with IL-17F (6, 7). Mature mouse IL-17A is 61% and 89% aa identical to human and rat IL-17A, respectively (4, 5, 8). While rodent and human mature sequences show modest aa sequence identity, human IL-17 is active on both mouse and rat cells (5, 9). Cells known to produce IL-17 are the CD4⁺ Th17 T cells, Paneth cells, GR1⁺CD11b⁺ myeloid suppressor cells, CD27- $\gamma\delta$ T cells, CD1⁺NK1.1⁺ iNKT cells and CD3-CD4⁺ LTI-like cells (3, 5, 6, 10-12).

A high affinity receptor for mouse IL-17 has been reported, and appears to be a heteromultimer of IL-17 RA and IL-17 RC, likely in a 2:1 ratio (1). IL-17 RA is a 130 kDa, type I transmembrane glycoprotein that bears no resemblance to members of the cytokine, TNF or immunoglobulin receptor superfamily (2, 10, 13). IL-17 RC is also a type I transmembrane protein, approximately 90-95 kDa in size, that shares less than 30% aa identity with IL-17 RA (14, 15). Both receptors are needed for IL-17A and IL-17A: F activity. The two receptors appear to form a functional association following ligand binding to IL-17 RA (1, 16).

IL-17 is best known for its participation in the recruitment and survival of neutrophils (3, 10, 17, 18). Its induction was initially described to be the result of antigen stimulation of DC, resulting in IL-23 secretion. In a TCR-independent event, IL-23 induces T cell production of IL-17 (3). Once secreted, IL-17, in the bone marrow would seem to induce stromal/fibroblast expression of both G-CSF and stem cell factor (membrane form), an effect that increases PMN differentiation and production. IL-17 may complement this by directly blocking neutrophil apoptosis, promoting greater circulating PMN numbers (17). In the tissues, IL-17 would also seem to promote neutrophil extravasation, principally through its effects on macrophages and endothelial cells (EC). On macrophages, IL-17 induces TNF- α , IL-1 β and IL-6 production (19). TNF- α and IL-1 β then act on local ECs to induce G-CSF secretion, an effect that is potentiated by IL-17 (20). IL-17 further contributes to PMN influx by inducing EC CXC chemokine release and NO production, which may increase vascular permeability (3, 9). IL-17 effects are not limited to neutrophils. In synovial joints, IL-17 upregulates RANKL expression on osteoblasts. This provides a stimulus for osteoclast formation and subsequent bone resorption (18).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for mouse IL-17 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse IL-17 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for mouse IL-17 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of mouse IL-17 bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Two samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision		Inter-assay Precision		
Sample	1	2	1	2	3
Mean (pg/mL)	40.6	432	41.2	51.9	480
Standard Deviation	3.5	26.2	3.7	4.6	40.9
CV%	8.6	6.1	9.0	8.9	8.5

B. RECOVERY

The recovery of mouse IL-17 spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 82 to 110% with an average of 94%.

The recovery of mouse IL-17 spiked to different levels throughout the range of the assay in mouse serum was evaluated. The recovery ranged from 79.5 to 92.0% with an average of 86.1%.

C. SENSITIVITY

The minimum detectable dose (MDD) of mouse IL-17 is typically less than 1.8 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *E.coli*-expressed recombinant mouse IL-17 produced at R&D Systems®.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of mouse IL-17 and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	107	100 - 111
1:4	109	104 - 119
1:8	107	95 - 120
1:16	102	88 - 120

F. SAMPLE VALUES

Cell Culture Supernates - Two spleen organ tissues from a mouse were homogenized and seeded in 100 mL of RPMI1640 supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 100 µg/mL streptomycin sulfate, and 10 µg/mL Con A for 2 days. The cell culture supernate was assayed for mouse IL-17 and measured 3424 pg/mL.

EL-4 cells (Mouse thymoma) were seeded at 2×10^5 cells/mL and cultured for 4 days in 100 mL of DMEM supplemented with 10% horse serum, 10 µg/mL PHA and 10 ng/mL PMA. The cell culture supernate was assayed for mouse IL-17 and measured 17 pg/mL.

Mouse serum - Four serum samples were evaluated for the presence of mouse IL-17 in this assay. All samples measured ranged from 33.5 to 46.0 pg/mL with an average of 38.4 pg/mL.

G. SPECIFICITY

This assay recognizes both natural and recombinant mouse IL-17. The following factors were prepared at 50 ng/mL and assayed for cross-reactivity. Preparations of the following factors at 50 ng/mL in a mid-range recombinant mouse IL-17 control were assayed for interference. No significant cross-reactivity or interference was observed.


Recombinant Mouse:	
IL-17B (aa 1-180)	IL-17C
IL-17B (aa 21-180)	IL-17D
IL-17BR	IL-17E
IL-17 RC	IL-17F
IL-17 RD	

There is 66% cross-reactivity observed with recombinant mouse IL-17A/F Heterodimer. At concentration 10 ng/mL or greater, mouse IL-17 R interferes in this assay. At concentrations greater than 78 pg/mL, rIL-17 R1 interferes in this assay.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse IL-17 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against mouse IL-17	1 plate
Mouse IL-17 Conjugate	Solution of an antibody against mouse IL-17 conjugated to horseradish peroxidase	1 vial
Mouse IL-17 Standard	Recombinant mouse IL-17 in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Calibrator Diluent(5×) /RD5P	A 5× concentrated buffered protein base used to dilute standard and samples	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant	1 vial
TMB Substrate	TMB ELISA Substrate Solution/ TMB Substrate Solution	1 vial
Stop Solution	Diluted hydrochloric acid solution	1 vial
Plate Sealers	Adhesive strip	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8 °C.*
	Stop Solution	
	Conjugate	
	TMB Substrate	
	Standard	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer.* Avoid repeated freeze-thaw cycles.
	Calibrator Diluent (5×)/RD5P	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

- Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 × g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

B. SAMPLE PREPARATION

Mouse serum samples recommend a 5-fold dilution. A suggested 5-fold dilution is 40 μ L of sample + 160 μ L of Calibrator Diluent (1×). Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Note: *Bring all reagents to room temperature before use.*


Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Calibrator Diluent (1×) - Use deionized or distilled water to prepare Calibrator Diluent (1×).

Mouse IL-17 Standard - Refer to the vial label for reconstitution volume*. This reconstitution produces a stock solution of 700 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 200 μ L of Calibrator Diluent (1×) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted standard 700 pg/mL serves as the high standard. The Calibrator Diluent (1×) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 μ L of Calibrator Diluent (1 $\times$) to each well.
4. Add 50 μ L of Standard and sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 μ L of mouse IL-17 conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
7. Repeat the aspiration/wash as in step 5.
8. Add 100 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
9. Add 100 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse IL-17 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Gaffen, S. (2009) *Nat. Rev. Immunol.* 9:556.
2. Iwakura, Y. et al. (2008) *Immunol. Rev.* 226:57.
3. Kolls, J.K. & A. Linden (2004) *Immunity* 21:467.
4. Yao, Z. et al. (1995) *J. Immunol.* 155:5483.
5. Kennedy, J. et al. (1996) *J. Interf. Cytokine Res.* 16:611.
6. Liang, S.C. et al. (2007) *J. Immunol.* 179:7791.
7. Chang, S. & C.Dong (2007) *Cell Res.* 17:435.
8. Rouvier, E. et al. (1993) *J. Immunol.* 150:5445.
9. Miljkovic, Dj. Et al. (2003) *Cell. Mol. Life Sci.* 60:518.
10. Witowski, J. et al. (2004) *Cell. Mol. Life Sci.* 61:567.
11. Cua, D.J. & C.M. Tato (2010) *Nat. Rev. Immunol.* 10:479.
12. Shin, H.C. et al. (1998) *Cytokine* 10:841.
13. Yao, Z. et al. (1995) *Immunity* 3:811.
14. Haudenschild, D. et al. (2002) *J. Biol. Chem.* 277:4309.
15. Toy, D. et al. (2006) *J. Immunol.* 177:36.
16. Hu, Y. et al. (2010) *J. Immunol.* 184:4307.
17. Schwarzenberger, P. et al. (2000) *J. Immunol.* 164:4783.
18. Yu, J.J. & S.L. Gaffen (2008) *Front. Biosci.* 13:170.
19. Jovanovic, D.V. et al. (1998) *J. Immunol.* 160:3513.
20. Numasaki, M. et al. (2004) *Immunol. Lett.* 95:97.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

	A	B	C	D	E	F	G	H
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								

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