PRODUCT INFORMATION & MANUAL

Mouse/Rat/Porcine/Canine TGF- β 1 Valukine™ ELISA

Catalog Number: VAL611

For the quantitative determination of natural and recombinant
Mouse/Rat/Porcine/Canine TGF- β 1 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202407.9

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	4
III. ADVANTAGES	5
IV. EXPERIMENT	8
V. KIT COMPONENTS AND STORAGE	9
VI. PREPARATION	11
VII. ASSAY PROCEDURE	14
VIII. RERERENCES	16

I. BACKGROUND

Transforming growth factor beta (TGF- β) proteins (including the three closely related mammalian isoforms TGF- β 1, -2 and -3) are pleiotropic cytokines that regulate extracellular matrix production, wound healing, immune functions, cell proliferation and differentiation. They belong to the large TGF- β superfamily, which also includes the activins/inhibins/MIS, bone morphogenetic proteins (BMPs), growth/differentiation factors (GDFs), Lefty1 and 2, and the distantly related GDNF family of neurotrophic factors. All family members show a characteristic cysteine-knot that is formed from multiple intrachain disulfide bonds (1-3).

The mouse, rat, porcine and canine TGF- β 1 cDNAs encode a 390 amino acid (aa) residue precursor that contains a 29 aa signal peptide and a 361 aa pro-protein. The pro-protein for each species is proteolytically processed via a furin-like convertase to generate an N-terminal 249 aa latency-associated peptide (LAP), and a C-terminal 112 aa mature TGF- β 1 (4-9). Both LAP and mature TGF- β 1 exist as disulfide-linked homodimers. After proteolytic cleavage and secretion, the two homodimers remain non-covalently associated as the small latent TGF- β 1 complex. In most cell types, this complex is also covalently linked via LAP to a latent TGF- β binding protein (LTBP). This creates a secreted, large latent complex. The TGF- β 1 present in either the small or the large latent complex is not available for TGF- β receptor binding and activation, and is, therefore, latent. Whereas LAP is both necessary and sufficient to confer latency to TGF- β 1, LTBP facilitates the proper folding and secretion of the small latent complex. LTBP is also a structural component of the extracellular matrix and directs the localization of the latent complex to the extracellular matrix (10-15). To date, four LTBPs that share multiple EGF-like, LTBP- and fibrillin-specific domains have been cloned. Three of the four LTBPs (#1, 3, and 4) have been shown to bind the small latent complexes of all TGF- β isoforms (16). Activation of the latent TGF- β complex is an important step that regulates TGF- β function in vivo. Multiple activation mechanisms have been identified. These involve protease-dependent (plasmin and matrix metalloprotease) and protease-independent (binding of LAP to thrombospondin 1 or a subset of integrins) pathways (5, 15, 17-21). Mature mouse TGF- β 1 shows 100% aa sequence identity to mature rat TGF- β 1 (6, 7), and both rat and mouse TGF- β 1 show 99% aa sequence identity with mature porcine, human, and

canine TGF- β 1 (4, 9, 22). Mature porcine TGF- β 1 shows absolute identity to both mature human and canine TGF- β 1.

The signalling high-affinity receptor for TGF- β 1 is a heteromeric complex consisting of a type I (TGF- β RI) and a type II (TGF- β RII) transmembrane serine/threonine kinase receptor. TGF- β RII is a constitutively active kinase. Upon binding TGF- β 1, it phosphorylates and activates the TGF- β RI. In turn, TGF- β RI phosphorylates and activates Smad proteins that regulate transcription. TGF- β RI, alternatively named activin receptor-like kinase (ALK-5), is present in almost all cell types. Two other type I receptors, ALK-1 and ALK-2 have also been implicated as alternative partners for TGF- β RII in the TGF- β signalling receptor complex. Besides the type I and type II receptors, accessory receptors including the type III receptor, TGF- β RIII (also known as betaglycan), and endoglin, which modulate TGF- β responses, have been identified (2, 23-26).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for TGF- β 1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TGF- β 1 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for TGF- β 1 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate solution (Chromogenic agent) is added to the wells and color develops in proportion to the amount of TGF- β 1 bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse/rat/porcine/canine serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1 $\times$) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Two samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision		Inter-assay Precision		
Sample	1	2	1	2	3
Mean (pg/mL)	88	198	816	97	209
Standard Deviation	3	4.9	19.6	8.1	17.3
CV%	3.4	2.5	2.4	8.4	8.3

B. RECOVERY

The recovery of TGF- β 1 spiked to different levels throughout the range of the assay in activated cell culture supernates was evaluated. The recovery ranged from 83 to 114% with an average of 105%.

The recovery of TGF- β 1 spiked to different levels throughout the range of the assay in activated serum was evaluated. The recovery ranged from 74.5 to 90.2% with an average of 78.9%.

C. SENSITIVITY

The minimum detectable dose (MDD) of TGF- β 1 is typically less than 1.7-15.4 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified *CHO* cell-expressed recombinant TGF- β 1 produced at R&D Systems®.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with

high concentrations of TGF- β 1 and diluted with Calibrator Diluent (1 $\times$) to produce samples with values within the dynamic range of the assay.

Dilution	Mouse		Rat		Porcine	
	Average% of Expected	Range (%)	Average% of Expected	Range (%)	Average% of Expected	Range (%)
1:2	99	95-103	102	101-102	100	98-102
1:4	99	92-104	104	102-107	102	99-104
1:8	97	92-99	106	102-111	114	112-116
1:16	93	80-103	108	102-115	107	95-119

F. SAMPLE VALUES

Cell Culture Supernates - Cell culture supernate samples were evaluated for the presence of TGF- β 1 in this assay.

Sample	Mean (pg/mL)
DMEM + 10% fetal calf serum	933
RPMI + 5% fetal calf serum	493
Mouse EL4/IL-2*	1522
PK15*	1414
Rat splenocytes*	1215

*Values have not been corrected for TGF- β 1 levels in the fetal calf serum.

EL-4 mouse lymphoblast cells (1×10^6 cells/mL) were cultured for 3 days in DMEM supplemented with 10% fetal calf serum and stimulated with 10 μ g/mL PHA and 10 ng/mL PMA. An aliquot of the cell culture supernate was removed and assayed for levels of TGF- β 1.

PK-15 porcine kidney epithelial cells (1×10^6 cells/mL) were grown to confluency in RPMI supplemented with 10% fetal calf serum and stimulated with 100 ng/mL LPS.

Samples were cultured for 24 hours. An aliquot of the cell culture supernate was removed and assayed for levels of TGF- β 1.

Rat splenocyte cells (1×10^6 cells/mL) were cultured for 3 days in DMEM supplemented with 10% fetal calf serum. An aliquot of the cell culture supernate was removed and assayed for levels of TGF- β 1.

Serum - Four mouse serum samples were evaluated for the presence of TGF- β 1 in this assay. All samples measured ranged from 94.5 to 120 ng/mL with an average of 106 ng/mL.

G. SPECIFICITY

This assay recognizes natural and recombinant TGF- β 1. This assay also recognizes human TGF- β 1.

The factors listed below were prepared at 50 ng/mL in Calibrator Diluent (1 $\times$) and assayed for cross reactivity. Preparations of the following factors at 50 ng/mL in a mid-range TGF- β 1 control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant Human		Recombinant Mouse
Activin A	BMPR-IB	BMP-3b
Activin RIA	BMPR-II	TGF- β RI
Activin RIIA	Follistatin288	BMPR-IA
Activin RIIB	Follistatin300	BMPR-IB
BMP-2	Follistatin315	
BMP-3	Inhibin A	Other Recombinant
BMP-3b	Inhibin B	Rat Agrin
BMP-4	LAP (TGF- β 1)	Zebrafish BMP-2a
BMP-5	TGF- α	Amphibian TGF- β 5
BMP-6	TGF- β RI	
BMP-8b	TGF- β RII	Nature protein
BMP-10	TGF- β RIII	Porcine TGF- β 2
BMP-15	TGF- β 2	
BMPR-IA	TGF- β 3	


Cross-reactivity - This kit detects Latent TGF- β 1 complex after acid treatment.

Interference - Significant interference was observed with rmTGF- β RII.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
TGF- β 1 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against TGF- β 1	1 plate
TGF- β 1 Conjugate	Solution of antibody against TGF- β 1 conjugated to horseradish peroxidase with preservatives	1 vial
TGF- β 1 Standard	Recombinant TGF- β 1 in a buffered protein base with preservatives; lyophilized. Refer to the vial label for reconstitution volume	2 vials
Calibrator Diluent (2 $\times$)/RD6-11	A 2 $\times$ concentrated buffered protein base used to dilute standard and samples	2 vials
Assay Diluent/RD1-21	A buffered protein base with preservatives	1 vial
Wash Buffer Concentrate (25 $\times$)	A 25 $\times$ concentrated solution of buffered surfactant with preservatives	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	1 vial
Stop Solution	Diluted hydrochloric acid	1 vial
Plate Sealers	Adhesive strips	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8 °C.*
	Stop Solution	
	Conjugate	
	Assay Diluent/RD1-21	
	TMB Substrate	
	Standard	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer.* Avoid repeated freeze-thaw cycles.
	Calibrator Diluent (2×)/RD6-11	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

- Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at 1000 × g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

B. SAMPLE ACTIVATION PROCEDURE

To activate latent TGF- β 1 to immunoreactive TGF- β 1 detectable by the TGF- β 1 immunoassay, follow the activation procedure below. Assay samples is pH 7.2-7.6 after neutralization. **Use polypropylene test tubes.**

Cell Culture Supernates	Serum
To 100 μ L of cell culture supernate, add 20 μ L of 1 N HCl.	To 40 μ L serum, add 20 μ L of 1 N HCl.
Mix well.	Mix well.
Incubate 10 minutes at room temperature.	Incubate 10 minutes at room temperature.
Neutralize the acidified sample by adding 20 μ L of 1.2 N NaOH/0.5 M HEPES.	Neutralize the acidified sample by adding 20 μ L of 1.2 N NaOH/0.5 M HEPES.
Mix well.	Mix well.
Assay immediately.	Prior to the assay, dilute the activated sample with Calibrator Diluent (1×). See the following for suggested dilutions.
The concentration read off the standard curve must be multiplied by the dilution factor, 1.4.	The concentration read off the standard curve must be multiplied by the appropriate dilution factors*.

***Activated serum samples may be stored for up to 24 hours at 2-8 °C before use. Activated cell culture supernates samples must be assayed immediately after activation. Do not freeze activated samples.**

C. SAMPLE PREPARATION

Cell culture supernate samples recommend at least a 1.4-fold dilution prior to the assay. A suggested 1.4-fold dilution is activation of the sample. Optimal dilutions should be determined by the end user.

*Activated serum samples recommend at least a 60-fold dilution prior to the assay. For example, add 10 µL of serum into a tube with 90 µL Calibrator Diluent (1×) to prepare a 10-fold diluted sample. Mix through and then pipette 20 µL of prepared 10-fold diluted sample into a tube with 100 µL Calibrator Diluent (1×) to prepare a final 60-fold diluted sample. (The final dilution factor, 120.). Optimal dilutions should be determined by the end user.

D. REAGENT PREPARATION

Note: *Bring all reagents to room temperature before use.*


Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Calibrator Diluent (1×) - Use deionized or distilled water to prepare Calibrator Diluent (1×).

TGF-1 Standard - Refer to the vial label for reconstitution volume*. Reconstitute the TGF-β1 Standard with Calibrator Diluent (1×). This reconstitution produces a stock solution of 2000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Use polypropylene tubes. Pipette 200 µL of Calibrator Diluent (1×) into each tube. Use the standard stock solution to produce a 2-fold dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted TGF-β1 Standard serves as the high standard (2000 pg/mL). Calibrator Diluent (1×) serves as the zero standard (0 pg/mL).

E. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- For best results, pipette reagents and samples into the center of each well.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- When using an automated plate washer, adding a 30 second soak period following the addition of wash buffer, and /or rotating the plate 180 degrees between steps may improve assay precision.
- TMB substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB substrate.

VII. ASSAY PROCEDURE

Note: Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 μ L of Assay Diluent to each well.
4. Add 50 μ L of Standard and sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 μ L of TGF- β 1 conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
7. Repeat the aspiration/wash as in step 5.
8. Add 100 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
9. Add 100 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the TGF- β 1 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. RERERENCES

1. Kingsley, D.M. (1994) *Genes Dev.* 8:133.
2. Javelaud, D. and A. Mauviel (2004) *Int. J. Biochem. Cell Biol.* 36:1161.
3. Chang, H. et al. (2002) *Endocr. Rev.* 23:787.
4. Kondaiah, P. et al. (1988) *J. Biol. Chem.* 263:18313.
5. Clark, D.A. and R. Coker (1998) *Int. J. Biochem. Cell Biol.* 30:293.
6. Qian, S.W. et al. (1990) *Nucleic Acids Res.* 18:3059.
7. Derynck, R. et al. (1986) *J. Biol. Chem.* 261:4377.
8. Dubois, C.M. et al. (1995) *J. Biol. Chem.* 270:10618.
9. Manning, A.M. et al. (1995) *Gene* 155:307.
10. Brunner, A.M. et al. (1989) *J. Biol. Chem.* 264:13660.
11. Mittl, P.R.E. et al. (1996) *Protein Sci.* 5:1261.
12. Gray, A.M. and A.J. Mason (1990) *Science* 247:1328.
13. Oklu, R. and R. Hesketh (2000) *Biochem. J.* 352:601.
14. Saharinen, J. et al. (1996) *EMBO J.* 15:245.
15. Mangasser-Stephan, K. and A.M. Gressner (1999) *Cell Tissue Res.* 297:363.
16. Koli, K. et al. (2001) *J. Cell Sci.* 114:2869.
17. Yang, L. et al. (2000) *Wound Rep. Reg.* 8:538.
18. Pedrozo, H.A. et al. (1999) *Endocrinology* 140:5806.
19. Yu, Q. and I. Stamenkovic (2000) *Genes Dev.* 14:163.
20. Maeda, S. et al. (2001) *J. Bone Miner. Res.* 16:1281.
21. Munger, J.S. et al. (1999) *Cell* 96:319.
22. Derynck, R. et al. (1985) *Nature* 316:701.
23. Derynck, R. and X-H. Feng (1997) *Biochim. Biophys. Acta* 1333:F105.
24. Ten Dijke, P. et al. (1996) *Curr. Opin. Cell. Biol.* 8:139.
25. Chai, Y. et al. (2003) *Crit. Rev. Oral Biol. Med.* 14:78.
26. De Caestecker, M. et al. (2004) *Cytokine Growth Factor Rev.* 15:1.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
	A	B	C	D	E	F	G	H

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