PRODUCT INFORMATION & MANUAL

Human/Mouse/Rat/Porcine/Canine TGF- β 1 Valukine™ ELISA

Catalog Number: VAL611

For the quantitative determination of natural and recombinant
Human/Mouse/Rat/Porcine/Canine TGF- β 1 concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202604.10

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	4
III. ADVANTAGES	5
IV. EXPERIMENT	15
V. KIT COMPONENTS AND STORAGE	16
VI. PREPARATION	18
VII. ASSAY PROCEDURE	24
VIII. RERERENCES	26

I. BACKGROUND

Human TGF- β 1 is a 25 kDa, disulfide-linked homodimeric protein involved in a number of key developmental, immunologic, and homeostatic processes (1-4). The molecule is synthesized as a 390 amino acid (aa) precursor that contains a 23 aa signal sequence, a 255 aa pro-region, and a 112 aa mature segment. Processing of the molecule is complex, and it is generally secreted as a latent form (5). Prior to release, the prepro-form is cleaved of its signal sequence, followed by glycosylation of its proregion. The glycosylation process includes the unusual attachment of mannose-6 phosphate residues. This is followed by furin convertase-mediated cleavage of the prohormone, creating an 80 kDa disulfide-linked proregion (termed LAP for latency-associated protein), plus a 25 kDa disulfidelinked mature segment (termed TGF- β 1) (6-8). These two independent disulfide-linked polypeptides associate in a non-covalent interaction that renders TGF- β 1 inactive. Although direct secretion of this 80K:25K complex can occur, it does so inefficiently. To facilitate secretion plus extracellular storage, a third 200 kDa component termed LTBP is covalently-linked to the N-terminus of one of the two LAP polypeptide chains. This promotes secretion and subsequent storage within the extracellular matrix (9, 10). After secretion, TGF- β 1, via LTBP, covalently links to ECM. This complex is later cleaved by proteases and released, exposing mannose residues on LAP. It is postulated that exposed LAP mannose residues now are able to bind to cell surface IGF-II R, where dissociative events disrupt the LAP-TGF- β 1 complex. This results in the release of active, homodimeric TGF- β 1 (7, 10). Mature mouse TGF- β 1 shares 100% aa sequence identity with rat and cotton rat TGF- β 1 (11, 12), 99% aa identity with human, canine, and porcine TGF- β 1 (13, 14, 15), and 97% aa identity with guinea pig TGF- β 1 (16). Relative to mouse TGF- β 2 and β 3, mature mouse TGF- β 1 shares 72% and 78% aa sequence identity, respectively (17, 18).

The traditional high-affinity receptor for TGF- β 1 is a heteromeric complex consisting of transmembrane serine/threonine kinases. Two types are involved; a constitutively phosphorylated, ligand-binding 80 kDa glycoprotein termed TGF- β 1 RII and a signal-transducing, non-ligand-binding 55 kDa glycoprotein termed TGF- β RI/ALK-5 (19-22). It is suggested that TGF- β 1 first binds TGF- β 1 RII, which then initiates a cross-phosphorylation of TGF- β 1 RI, culminating in signal transduction. There is also a third TGF- β receptor termed TGF- β 1 RIII, which can be either the 250 kDa

proteoglycan named betaglycan, or the 180 kDa glycoprotein termed endoglin/CD105 (23, 24). It has been proposed that TGF- β 1 RIII captures TGF- β and "passes" it to TGF- β 1 RII (20). This is perhaps true for betaglycan but not endoglin. Endoglin does not bind TGF- β by itself; only within the context of TGF- β 1 RII ligand binding. Evidence suggests that rather than "passing" on ligand, endoglin may actually enter the receptor complex and modulate TGF- β downstream signaling (25, 26). Finally, and although ALK-5 has traditionally been assumed to be the only type I signaling receptor for TGF- β 1, it is also possible that ALK-1 may serve as a condition-dependent, type I TGF- β receptor (27).

TGF- β 1 has a wide range of activities. During an immune response, TGF- β 1 impacts antibody production by preferentially inducing IgA production in both mouse and human (28). It also regulates dendritic cell chemotaxis by altering the expression of chemokine receptors (29). Finally, it can downmodulate an inflammatory response by dampening macrophage activity and proinflammatory cytokine secretion (30). During wound healing, TGF- β 1 is released from activated platelets. This local source of TGF- β 1 has marked stimulatory effects on fibroblasts, where it induces matrix synthesis; on monocytes, where it induces proinflammatory mediator and growth factor secretion; and on keratinocytes, where it may promote keratinocyte proliferation by downmodulating its own signaling pathway (31). Finally, during development, TGF- β 1 may play a role in endochondral ossification, and its absence results in severely defective yolk sac vasculogenesis and hematopoiesis (32, 33).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for TGF- β 1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TGF- β 1 present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for TGF- β 1 is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate solution (Chromogenic agent) is added to the wells and color develops in proportion to the amount of TGF- β 1 bound in the initial step. The color development is stopped and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and human/mouse/rat/porcine/canine serum, human urine, human/mouse/rat/porcine/canine platelet-poor plasma.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1 \times) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

CELL CULTURE SUPERNATE ASSAY

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	317	683	1271	312	657	1184
Standard Deviation	11.3	42.2	31.1	26.1	46.6	92.3
CV%	3.6	6.2	2.4	8.4	7.4	7.8

SERUM/PLASMA ASSAY

Sample	Intra-assay Precision			Inter-assay Precision		
	1	2	3	1	2	3
Mean (pg/mL)	309	708	1072	303	622	1092
Standard Deviation	15.3	40.0	44.3	21.2	47.2	70.0
CV%	5.0	5.6	4.1	7.0	7.6	6.4

B. RECOVERY

The recovery of TGF- β 1 spiked to levels throughout the range of the assay in activated samples was evaluated.

Sample Type	Average % Recovery	Range (%)
Human urine (n=4)	114	105-123
Human platelet-poor EDTA plasma (n=4)	90	77-107
Human platelet-poor heparin plasma (n=4)	90	83-98
Media + FBS (n=4)	103	80-125
Serum-free media (n=3)	107	96-117
Porcine platelet-poor EDTA plasma (n=4)	87	82-95
Porcine platelet-poor heparin plasma (n=4)	85	76-95

C. SENSITIVITY

Thirty-three assays were evaluated and the minimum detectable dose (MDD) of human/mouse/rat/porcine/canine TGF- β 1 ranged from 0.889-5.50 pg/mL. The mean MDD was 2.38 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against highly purified CHO cell-expressed recombinant TGF- β 1 produced at R&D Systems®.

The NIBSC/WHO TGF- β 1 International Standard 89/514 (Human, rDNA Derived) was evaluated in this kit. The dose response curve of the reference reagent 89/514 parallels the Valukine standard curve. To convert sample values obtained with the Valukine Human/Mouse/Rat/Porcine/Canine TGF- β 1 kit to approximate NIBSC/WHO 89/514 Units, use the equation below:

Cell culture supernate/Urine - NIBSC/WHO 89/514 approximate value (IU/mL) =
 $0.0159 \times \text{Valukine Human/Mouse/Rat/Porcine/Canine TGF-}\beta\text{1 value (pg/mL)}$

Serum/Platelet-poor Plasma - NIBSC/WHO 89/514 approximate value (IU/mL) =
 $0.0171 \times \text{Valukine Human/Mouse/Rat/Porcine/Canine TGF-}\beta\text{1 value (pg/mL)}$

Note: Based on data generated in May 2021.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of TGF- β 1 and diluted with Calibrator Diluent (1 \times) to produce samples with values within the dynamic range of the assay.

Human Samples		Cell culture supernates (n=4)	Serum* (n=4)	Urine (n=4)	Platelet-poor	
					EDTA plasma* (n=4)	Heparin plasma* (n=4)
1:2	Average % of Expected	100	97	93	105	105
	Range (%)	90-114	95-102	92-94	103-107	102-111
1:4	Average % of Expected	93	97	88	110	107
	Range (%)	82-117	94-104	83-91	108-113	103-112
1:8	Average % of Expected	93	95	88	112	109
	Range (%)	76-121	93-100	83-92	111-112	102-116
1:16	Average % of Expected	100	97	92	120	118
	Range (%)	85-124	92-104	86-101	115-124	110-125

*Samples were diluted after activation. See the Sample Activation Procedure.

Mouse Samples		Cell culture supernates (n=4)	Serum* (n=4)	Platelet-poor	
				EDTA plasma* (n=4)	Heparin plasma* (n=4)
1:2	Average % of Expected	97	100	104	105
	Range (%)	91-102	92-105	90-115	103-111
1:4	Average % of Expected	90	100	106	107
	Range (%)	86-96	92-104	95-115	102-113
1:8	Average % of Expected	86	97	101	111
	Range (%)	79-95	91-101	88-110	100-119
1:16	Average % of Expected	85	96	99	109
	Range (%)	79-94	87-101	88-110	93-117

*Samples were diluted after activation. See the Sample Activation Procedure.

Rat Samples		Cell culture supernates (n=2)	Serum* (n=4)	Platelet-poor	
				EDTA plasma* (n=4)	Heparin plasma* (n=4)
1:2	Average % of Expected	107	100	102	104
	Range (%)	107-108	93-104	99-109	101-108
1:4	Average % of Expected	101	102	101	106
	Range (%)	100-103	99-107	96-104	104-110
1:8	Average % of Expected	97	100	103	107
	Range (%)	96-97	96-103	99-109	103-111
1:16	Average % of Expected	97	101	104	106
	Range (%)	96-98	96-108	97-115	97-113

*Samples were diluted after activation. See the Sample Activation Procedure.

Porcine Samples		Cell culture supernates (n=4)	Serum* (n=4)	Platelet-poor	
				EDTA plasma* (n=5)	Heparin plasma* (n=4)
1:2	Average % of Expected	97	103	106	110
	Range (%)	88-106	100-106	102-111	108-111
1:4	Average % of Expected	92	110	112	100
	Range (%)	81-102	97-106	109-111	108-115
1:8	Average % of Expected	94	99	112	117
	Range (%)	81-103	95-104	110-114	108-124
1:16	Average % of Expected	97	104	117	119
	Range (%)	86-106	97-109	115-118	109-124

*Samples were diluted after activation. See the Sample Activation Procedure.

Canine Samples		Cell culture supernates (n=3)	Serum* (n=4)	Platelet-poor	
				EDTA plasma* (n=4)	Heparin plasma* (n=4)
1:2	Average % of Expected	99	100	98	100
	Range (%)	91-104	96-102	91-103	96-105
1:4	Average % of Expected	99	99	94	96
	Range (%)	96-102	95-103	87-98	87-107
1:8	Average % of Expected	100	100	92	98
	Range (%)	93-109	97-103	87-95	88-106
1:16	Average % of Expected	104	100	95	101
	Range (%)	95-120	95-107	86-105	96-106

*Samples were diluted after activation. See the Sample Activation Procedure.

F. SAMPLE VALUES

Serum/Plasma/Urine - Samples were evaluated for the presence of TGF- β 1 in this assay.

Human Samples	Mean (pg/mL)	Range (pg/mL)	Standard Deviation (pg/mL)
Serum (n=31)	51,640	32,091-95,147	11,901
Platelet-poor EDTA plasma (n=31)	2377	1414-4641	729
Platelet-poor heparin plasma (n=31)	2383	1445-3653	561
Urine* (n=10)	65.9	ND-108	—

*Only 40% of the urine samples measured detectable levels (> 31.3 pg/mL).

ND=Non-detectable

Mouse Samples	Mean (pg/mL)	Range (pg/mL)	Standard Deviation (pg/mL)
Serum (n=10)	97,352	78,751-114,130	12,623
Platelet-poor EDTA plasma (n=10)	41,763	16,577-68,505	16,887
Platelet-poor heparin plasma (n=10)	41,101	18,195-66,700	17,419

Rat Samples	Mean (pg/mL)	Range (pg/mL)	Standard Deviation (pg/mL)
Serum (n=10)	68,108	43,721-89,138	14,343
Platelet-poor EDTA plasma (n=5)	5833	4912-7820	1141
Platelet-poor heparin plasma (n=5)	13,971	5274-26,114	8392

Porcine Samples	Mean (pg/mL)	Range (pg/mL)	Standard Deviation (pg/mL)
Serum (n=5)	17,597	13,481-23,914	4146
Platelet-poor EDTA plasma (n=5)	1955	1530-2810	518
Platelet-poor heparin plasma (n=5)	1567	1111-2475	551

Canine Samples	Mean (pg/mL)	Range (pg/mL)	Standard Deviation (pg/mL)
Serum (n=10)	33,455	2193-79,590	19,795
Platelet-poor EDTA plasma** (n=5)	2195	ND-2515	—
Platelet-poor heparin plasma*** (n=5)	3763	ND-4278	—

**Only 60% of the EDTA plasma samples measured detectable levels with the required 40 fold dilution. Two of five samples read just below the standard curve (< 31.3 pg/mL).

***Only 80% of the heparin plasma samples measured detectable levels with the required 40 fold dilution. One of five samples read just below the standard curve (< 31.3 pg/mL).

Cell Culture Supernates - Human peripheral blood mononuclear cells (PBMCs) were separated from whole blood by a density gradient centrifugation method using Ficoll-Paque Plus. CD4⁺ T cells were isolated from PBMCs using the MagCelect™ Human CD4⁺ T cell Isolation Kit (R&D Systems®, Catalog # MAGH102). Cells were seeded at 5 x 10⁵ /mL and cultured using Excellerate™ Human T Cell Expansion Media, Xeno-Free (R&D Systems, Catalog # CCM030). T cells were left untreated or treated with 10 ng/mL GMP recombinant human (rh) IL-7 (R&D Systems, Catalog # 207-GMP), 10 ng/mL GMP rhIL-15 (R&D Systems, Catalog # 247-GMP), and stimulated via their T cell receptor (TCR) and co-stimulatory receptor for 5 days. TCR stimulation was mediated using 25 µL Cloudz™ CD3/28 particles (Cloudz T Cell Activation Kit - CD3/CD28, (R&D Systems, Catalog # CLD001) per mL of culture

media. CD4⁺ T cells were maintained in a 5% CO₂ incubator at 37 °C for 5 days. An aliquot of the cell culture supernates was removed, assayed for TGF-β1, and measured 77.7 pg/mL (untreated) and 471 pg/mL (treated).

Human PBMCs were seeded at 1 x 10⁶ /mL and cultured in RPMI supplemented with 10% fetal bovine serum (FBS), 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin and left untreated. An aliquot of the cell culture supernate was removed, assayed for TGF-β1, and measured 2854 pg/mL.

Mouse EL-4 cells were cultured in DMEM High Glucose supplemented with 10% FBS, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin. Cells were then treated with 10 ng/mL PMA and 10 µg/mL PHA for 24 hours. An aliquot of the cell culture supernate was removed, assayed for TGF-β1, and measured 2969 pg/mL.

Rat spleen was taken from a pregnant Sprague Dawley rat, homogenized, and cultured in DME with 10% FBS, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin at 37 °C and 5% CO₂. Rat splenocytes were treated with 50 ng/mL recombinant rat IL-2 and 5 µg/mL PHA for 3 days. An aliquot of the cell culture supernate was removed, assayed for TGF-β1, and measured 2716 pg/mL.

Porcine PK-15 cells were cultured in MEM supplemented with 10% FBS, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin. Cells were treated with 1 µg/mL LPS for 72 hours. An aliquot of the cell culture supernate was removed, assayed for TGF-β1, and measured 1631 pg/mL.

G. SPECIFICITY

This assay recognizes natural and recombinant TGF-β1. This assay also recognizes human TGF-β1.2.

The factors listed below were prepared at 50 ng/mL in Calibrator Diluent (1×) and assayed for cross reactivity. Preparations of the following factors at 50 ng/mL in a mid-range TGF-β1 control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant human		Recombinant mouse	Other recombinant
Activin A	BMPR-IB	Activin RI	rat Agrin
Activin RIA	BMPR-II	Activin RIIA	zebrafish BMP-2
Activin RIIA	Follistatin ₂₈₈	Activin RIB	amphibian TGF- β 5
BMPR-IB	Follistatin ₃₀₀	Activin RIIB	chinese hamster TGF- β 1-LAP
Activin RIIB	Follistatin ₃₁₅	BMP-3	
Agrin	GDF-5	BMP-3b	Natural protein
BMP-1.1	GDF-7	BMP-4	porcine TGF- β 2
BMP-2	GDF-9	BMP-6	
BMP-3	GDF-11	BMP-7	
BMP-3b	GDF-15	BMP-9	
BMP-4	Inhibin A	BMPR-IA	
BMP-5	Inhibin B	BMPR-IB	
BMP-6	LAP (TGF- β 1)	BMPR-II	
BMP-7	TGF- α	Follistatin ₂₈₈	
BMP-8a	TGF- β RI	Inhibin A	
BMP-8b	TGF- β RIII	TGF- β RI	
BMP-10	TGF- β 3	TGF- β RIII	
BMP-15			
BMPR-IA			

Cross-reactivity - Cross-reactivity was observed with the factors listed below.

Recombinant Protein	Assay Diluent	
	RD1-21	RD1-73
Human BMP-9	<1%	<1%
Mouse BMP-10	<1%	<1%
Human Latent TGF- β 1	<1%	<1%
Activated Human Latent TGF- β 1	15%	15%
Equine TGF- β 1	3%	<1%
Human TGF- β 1.2	14%	12%
Human TGF- β 2	13%	12%

Interference - Significant interference was observed with recombinant human TGF- β RII and recombinant mouse TGF- β RII.


Human or Mouse TGF- β RII (ng/mL)	TGF- β 1 (pg/mL)	% TGF- β 1 Detected in presence of Human TGF- β RII		% TGF- β 1 Detected in presence of Mouse TGF- β RII	
		RD1-21	RD1-73	RD1-21	RD1-73
0	500	100	100	100	100
1.56	500	98	99	90	93
3.13	500	93	95	85	88
6.25	500	82	86	77	79
12.5	500	63	66	69	67
25	500	42	37	60	55
50	500	28	18	52	40

IV. EXPERIMENT

EXAMPLE STANDARD


The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

CELL CULTURE SUPERNATE/URINE ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.027 0.031	0.029	—
31.3	0.070 0.071	0.071	0.042
62.5	0.117 0.119	0.118	0.089
125	0.208 0.214	0.211	0.182
250	0.385 0.407	0.396	0.367
500	0.710 0.720	0.715	0.686
1000	1.243 1.252	1.248	1.219
2000	2.035 2.040	2.038	2.009

SERUM/PLATELET-POOR PLASMA ASSAY

(pg/mL)	O.D.	Average	Corrected
0	0.013 0.014	0.014	—
31.3	0.064 0.065	0.065	0.051
62.5	0.123 0.124	0.124	0.110
125	0.218 0.230	0.224	0.210
250	0.413 0.423	0.418	0.404
500	0.774 0.785	0.780	0.766
1000	1.395 1.460	1.428	1.414
2000	2.249 2.307	2.278	2.264

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
TGF- β 1 Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against TGF- β 1	1 plate
TGF- β 1 Conjugate	Solution of antibody against TGF- β 1 conjugated to horseradish peroxidase with preservatives	1 vial
TGF- β 1 Standard	Recombinant TGF- β 1 in a buffered protein base with preservatives; lyophilized. <i>Refer to the vial label for reconstitution volume</i>	2 vials
Calibrator Diluent (2 \times) RD6-11	A 2 \times concentrated buffered protein base used to dilute standard and samples	1 vial
Assay Diluent RD1-21	A buffered protein base with preservatives. <i>For cell culture supernate/urine samples</i>	1 vial
Assay Diluent RD1-73	A buffered base with preservatives. <i>For serum/platelet-poor plasma samples</i>	1 vial
Wash Buffer Concentrate (25 \times)	A 25 \times concentrated solution of buffered surfactant with preservatives	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution	1 vial
Stop Solution	Diluted hydrochloric acid	1 vial
Plate Sealers	Adhesive strips	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8 °C.*
	Stop Solution	
	Conjugate	
	Assay Diluent RD1-21	
	Assay Diluent RD1-73	
	TMB Substrate	
	Standard	Use a new standard for each assay. Discard after use.
Calibrator Diluent (2×) RD6-11	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.	
Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*	

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- Pipettes and pipette tips.
- Deionized or distilled water.
- Squirt bottle, manifold dispenser, or automated microplate washer.
- 500 mL graduated cylinder.

D. PRECAUTION

- Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Activated samples must be assayed immediately. Do not freeze activated samples.

Cell Culture Supernates - Remove particulates by centrifugation and assay (see activation procedure) immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Note: *Animal serum used in the preparation of cell culture media may contain high levels of latent TGF- β 1. For best results, do not use animal serum for growth of cell cultures when assaying for TGF- β 1 production. If animal serum is used as a supplement in the media, precautions should be taken to prepare the appropriate control and run the control in the immunoassay to determine the baseline concentration of TGF- β 1.*

Human Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature. Centrifugation for 15 minutes at 1000 \times g. Remove serum and assay (see activation procedure) immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Non-human Serum - Allow blood samples to clot for 2 hours at room temperature before centrifuging for 20 minutes at 2000 \times g. Remove serum and assay (see activation procedure) immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Human Urine - Aseptically collect the first urine of the day (mid-stream), voided directly into a sterile container. Centrifuge to remove particulate matter. Assay immediately (see activation procedure) or aliquot and store at ≤ -20 °C. Avoid repeated freeze/thaw cycles.

Note: *Neat unactivated urine samples exhibit a decrease in TGF- β 1 concentration in the first 24 hours of storage (frozen or refrigerated). Care should be taken that samples are assayed under identical storage conditions and durations.*

Platelet-poor Plasma* - Collect plasma on ice using EDTA or heparin as an

anticoagulant. Centrifuge for 20 minutes at 1000 × g within 30 minutes of collection. An additional centrifugation step of the plasma at 10,000 x g for 10 minutes at 2-8 °C is recommended for complete platelet removal. Assay (see activation procedure) immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Note: *Citrate plasma has not been validated for use in this assay.*

***TGF- β 1 is present in platelet granules and is released upon platelet activation. Therefore, to measure circulating levels of TGF- β 1, platelet-poor plasma should be collected for measurement. It should be noted that many protocols for plasma preparation, including procedures recommended by the Clinical and Laboratory Standards Institute (CLSI), result in incomplete removal of platelets from plasma. This will cause variable and irreproducible results for assays of factors contained in platelets and released by platelet activation. The recommended plasma collection protocol is designed to minimize platelet degranulation. However, since even the best methods for plasma collection may result in some platelet degranulation on occasion, it is recommended that markers for platelet degranulation be determined in samples containing elevated TGF- β 1 levels.**

B. SAMPLE ACTIVATION PROCEDURE

ACTIVATION REAGENT PREPARATION

To activate latent TGF- β 1 to immunoreactive form, prepare the following solutions for acid activation and neutralization. The solutions may be stored in polypropylene bottles at room temperature for up to one month.

Caution: Wear protective clothing and safety glasses during preparation and use of these reagents. Refer to the appropriate SDS before use.

1 N HCl (100 mL) - To 91.67 mL of deionized water, slowly add 8.33 mL of 12 N HCl. Mix well.

1.2 N NaOH/0.5 M HEPES (100 mL) - To 75 mL of deionized water, slowly add 12 mL of 10 N NaOH. Mix well. Add 11.9 g of HEPES. Mix well. Bring final volume to 100 mL with deionized water.

For each new lot of acidification and neutralization reagents, measure the pH of several representative samples after neutralization to ensure that it is within pH 7.2-7.6. Adjust the volume and corresponding dilution factor of the neutralization

reagent as needed.

TGF- β 1 SAMPLE ACTIVATION PROCEDURE

To activate latent TGF- β 1 to immunoreactive TGF- β 1 detectable by the TGF- β 1 immunoassay, follow the activation procedure below. Assay samples is pH 7.2-7.6 after neutralization. **Use polypropylene test tubes.**

Note: Do not activate the kit standards. *The standards contain active recombinant TGF- β 1.*

Cell Culture Supernates	Serum
To 100 μ L of cell culture supernate, add 20 μ L of 1 N HCl.	To 40 μ L serum, add 20 μ L of 1 N HCl.
Mix well.	Mix well.
Incubate 10 minutes at room temperature.	Incubate 10 minutes at room temperature.
Neutralize the acidified sample by adding 20 μ L of 1.2 N NaOH/0.5 M HEPES.	Neutralize the acidified sample by adding 20 μ L of 1.2 N NaOH/0.5 M HEPES.
Mix well.	Mix well.
Prior to the assay, dilute the activated sample with calibrator diluent (1 \times). See <i>the following for suggested dilutions.</i>	Prior to the assay, dilute the activated sample with Calibrator Diluent (1 \times). See <i>the following for suggested dilutions.</i>
The concentration read off the standard curve must be multiplied by the dilution factor, 1.4.	The concentration read off the standard curve must be multiplied by the appropriate dilution factors.

C. SAMPLE PREPARATION

Cell culture supernate samples (human and non-human) tested neat, but may require dilution if high endogenous levels are present (final dilution factor of sample is 1.4 if tested neat). Optimal dilutions should be determined by the end user.

Human serum samples require a 40-fold dilution in Calibrator Diluent (1 \times)[†]. A suggested 40-fold dilution is 10 μ L of activated sample + 390 μ L of Calibrator Diluent

(1×) (final dilution factor of sample is 60). Optimal dilutions should be determined by the end user.

Human urine samples tested neat (final dilution factor of sample is 1.4). Optimal dilutions should be determined by the end user.

Human platelet-poor plasma* samples require a 8-fold dilution in Calibrator Diluent (1×)[†]. A suggested 8-fold dilution is 25 µL of activated sample + 175 µL of Calibrator Diluent (1×) (final dilution factor of sample is 12). Optimal dilutions should be determined by the end user.

Mouse and rat serum/platelet-poor plasma* samples require a 60-fold dilution in Calibrator Diluent (1×)[†]. A suggested 60-fold dilution is 10 µL of activated sample + 590 µL of Calibrator Diluent (1×) (final dilution factor of sample is 90). Optimal dilutions should be determined by the end user.

Porcine serum samples require a 15-fold dilution in Calibrator Diluent (1×). A suggested 15-fold dilution is 10 µL of activated sample + 140 µL of Calibrator Diluent (1×) (final dilution factor of sample is 22.5). Optimal dilutions should be determined by the end user.

Porcine platelet-poor plasma* samples require a 5-fold dilution in Calibrator Diluent (1×). A suggested 5-fold dilution is 40 µL of activated sample + 160 µL of Calibrator Diluent (1×) (final dilution factor of sample is 7.5). Optimal dilutions should be determined by the end user.

Canine serum and platelet-poor plasma* samples require a 40-fold dilution in Calibrator Diluent (1×). A suggested 40-fold dilution is 10 µL of activated sample + 390 µL of Calibrator Diluent (1×) (final dilution factor of sample is 60). Optimal dilutions should be determined by the end user.

* May require different dilutions depending on the extent of platelet contamination.

† See Reagent Preparation section.

D. REAGENT PREPARATION

Bring all reagents to room temperature before use.


Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Calibrator Diluent (1×) - Use deionized or distilled water to prepare Calibrator Diluent (1×).

TGF-1 Standard - Refer to the vial label for reconstitution volume*. Reconstitute the TGF-β1 Standard with Calibrator Diluent (1×). This reconstitution produces a stock solution of 2000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Use polypropylene tubes. Pipette 200 μL of Calibrator Diluent (1×) into each tube. Use the standard stock solution to produce a 2-fold dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted TGF-β1 Standard serves as the high standard (2000 pg/mL). Calibrator Diluent (1×) serves as the zero standard (0 pg/mL).

E. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- For best results, pipette reagents and samples into the center of each well.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.

- When using an automated plate washer, adding a 30 second soak period following the addition of wash buffer, and /or rotating the plate 180 degrees between steps may improve assay precision.
- TMB substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB substrate.

VII. ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.

1. Prepare all reagents, working standards, and activated samples as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 μL of Assay Diluent RD1-21 (*for cell culture supernate/urine samples*) or Assay Diluent RD1-73 (*for serum/platelet-poor plasma samples*) to each well.
4. Add 50 μL of Standard, and activated sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μL) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 μL of TGF- β 1 conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
7. Repeat the aspiration/wash as in step 5.
8. Add 100 μL of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
9. Add 100 μL of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density. Create a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the TGF- β 1 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. RERERENCES

1. Derynck, R. *et al.* (1986) *J. Biol. Chem.* 261:4377.
2. Padgett, R.W. and G.I. Patterson (2001) *Developmental Cell* 1:343.
3. Cox, D.A. and T. Maurer (1997) *Clin. Immunol. Immunopathol.* 83:25.
4. Ruscetti, F.W. and S.H. Bartelmez (2001) *Int. J. Hematol.* 74:18.
5. Gleizes, P-E. *et al.* (1997) *Stem Cells* 15:190.
6. Dubois, C.M. *et al.* (1995) *J. Biol. Chem.* 270:10618.
7. Khalil, N. (1999) *Microbes Infect.* 1:1255.
8. Clark, D.A. and R. Coker (1998) *Int. J. Biochem. Cell Biol.* 30:293.
9. Koli, K. *et al.* (2001) *Microsc. Res. Tech.* 52:354.
10. Mangasser-Stephan, K. and A.M. Gressner (1999) *Cell Tissue Res.* 297:363.
11. Qian, S.W. *et al.* (1990) *Nucleic Acids Res.* 18:3059.
12. Genbank Accession #: AAL87199.
13. Derynck, R. *et al.* (1985) *Nature* 316:701.
14. Manning, A.M. *et al.* (1995) *Gene* 155:307.
15. Kondaiah, P. *et al.* (1988) *J. Biol. Chem.* 263:18313.
16. Scarozza, A.M. *et al.* (1998) *Cytokine* 10:851.
17. Miller, D.A. *et al.* (1989) *Mol. Endocrinol.* 3:1108.
18. Denhez, F. *et al.* (1990) *Growth Factors* 3:139.
19. Derynck, R. and X-H. Feng (1997) *Biochim. Biophys. Acta* 1333:F105.
20. Ten Dijke, P. *et al.* (1996) *Curr. Opin. Cell. Biol.* 8:139.
21. Lawler, S. *et al.* (1994) *Development* 120:165.
22. Susuki, A. *et al.* (1994) *Biochem. Biophys. Res. Commun.* 198:1063.
23. Lopez-Casillas, F. *et al.* (1994) *J. Cell Biol.* 124:557.
24. St. Jacques, S. *et al.* (1994) *Endocrinology* 134:2645.
25. Barbara, N.P. *et al.* (1999) *J. Biol. Chem.* 274:584.
26. Ma, X. *et al.* (2000) *Arterioscler. Thromb. Vasc. Biol.* 20:2546.
27. Lux, A. *et al.* (1999) *J. Biol. Chem.* 274:9984.

28. Stavnezer, J. (1995) *J. Immunol.* 155:1647.
29. Sato, K. *et al.* (2000) *J. Immunol.* 164:2285.
30. Wahl, S.M. *et al.* (2000) *Cytokine Growth Factor Rev.* 11:71.
31. Ashcroft, G.S. and A.B. Roberts (2000) *Cytokine Growth Factor Rev.* 11:125.
32. Matsunaga, S. *et al.* (1999) *Int. J. Oncol.* 14:1063.
33. Mummery, C.L. *et al.* (1999) *Int. J. Dev. Biol.* 43:693.
34. Danielpour, D. *et al.* (1989) *Growth Factors.* 2:61.
35. Danielpour, D. *et al.* (1993) *J. Immunol. Meth.* 158:17.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
	A	B	C	D	E	F	G	H	

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