PRODUCT INFORMATION & MANUAL

Mouse alpha-Fetoprotein/AFP Valukine™ ELISA

Catalog Number: VAL622

For the quantitative determination of natural and recombinant
mouse alpha-Fetoprotein/AFP concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202502.2

TABLE OF CONTENTS

I. BACKGROUND	2
II. OVERVIEW	3
III. ADVANTAGES	4
IV. EXPERIMENT	6
V. KIT COMPONENTS AND STORAGE	7
VI. PREPARATION	9
VII. ASSAY PROCEDURE	11
VIII. REFERENCES	13

I. BACKGROUND

α -Fetoprotein (AFP), a member of the albuminoid superfamily (Albumin, Vitamin D-binding protein, and α -Albumin), is a fetal- and tumor-associated protein well known as a marker for certain cancers and congenital defects (1-4). In humans and rodents, it is a single chain glycoprotein of 66-72 kDa, with 3-5% of the molecular weight resulting from glycosylation (3). Mouse AFP shares 81% and 65% amino acid sequence identity with rat and human AFP, respectively. AFP is synthesized in the fetus primarily by the liver, yolk sac, and tissues of gastrointestinal origin (5-8). During mouse development, AFP is expressed at low but detectable levels in the foregut endoderm, becomes highly activated in hepatoblasts that form the liver bud, and is one of the earliest markers of the hepatocyte lineage (9). Fetal serum levels peak at weeks 10 through 13 and then decline throughout gestation (10). In contrast, maternal serum AFP, derived primarily from fetal circulation, continues to increase into the third trimester before declining (10).

AFP can act as a carrier protein, binding several ligands including steroids, bilirubin, fatty acids, retinoids, and flavonoids (4). In rodents, AFP sequesters endogenous estrogen and blocks the induction of estrogen receptor (ER)-mediated responses (2). Ligand binding may affect the rigidity and conformation of the AFP tertiary structure (2, 4, 11, 12). In addition, AFP has been proposed to have immunosuppressive activity, regulate cell proliferation and apoptosis, initiate intracellular signaling, and contribute to cell invasion (4, 13). Putative cell surface receptors for AFP have been described, and AFP is internalized by many cell types (14). AFP knockout mice are viable, but females are infertile due to an inadequate hormonal environment that blocks normal ovulation (15, 16).

Altered levels of both fetal and maternal AFP have been associated with several congenital abnormalities including hypothyroidism, autoimmune disorders, and heart defects (2, 17-19). In addition, AFP is widely used as a prenatal screening marker for neural tube defects and chromosomal abnormality (2, 20). Low maternal serum AFP levels have been associated with a higher incidence of Down's syndrome, whereas higher levels are associated with spina bifida and anencephaly (21-24). Postnatal production of AFP is associated with several pathologies. AFP is a commonly used marker for hepatocellular cancer, and elevated AFP levels are coincident with several other cancers including hepatoblastoma, germ cell tumors, and certain gastric cancers (25-31).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for mouse AFP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any mouse AFP present is bound by the immobilized antibody. After washing away any unbound substances, a biotinylated detection antibody specific for mouse AFP is pipetted into the wells. After washing away any unbound substances, streptavidin-HRP is pipetted into the wells. Following a wash to remove any unbound reagent, TMB Substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of mouse AFP bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- **FOR RESEARCH USE ONLY NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- This kit is suitable for cell culture supernate and mouse serum.
- The kit should not be used beyond the expiration date on the kit label.
- Do not mix or substitute reagents with those from other lots or sources.
- If samples generate values higher than the highest standard, dilute the samples with Reagent Diluent (1×) and repeat the assay.
- Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	1947.6	504.8	124.5	2086.3	601.6	133.3
Standard Deviation	81.9	20.9	5.6	186.7	38.2	12.8
CV%	4.2	4.1	4.5	8.9	7.2	9.6

B. RECOVERY

The recovery of mouse AFP spiked to different levels throughout the range of the assay in cell culture media was evaluated. The recovery ranged from 106.9 to 129.3% with an average of 118.3%.

The recovery of mouse AFP spiked to different levels throughout the range of the assay in mouse serum was evaluated. The recovery ranged from 87.0 to 113.7% with an average of 101.3%.

C. SENSITIVITY

The minimum detectable dose (MDD) of mouse AFP is typically less than 1.23 pg/mL. The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified NS0-expressed recombinant mouse AFP produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of mouse AFP and diluted with Reagent Diluent (1×) to produce samples with values within the dynamic range of the assay.

Dilution	Average % of Expected	Range (%)
1:2	93.6	89.2-95.9
1:4	92.9	85.1-103.7
1:8	91.9	83.0-105.5
1:16	93.7	85.7-110.2

F. SAMPLE VALUES

Mouse serum - Five mouse serum samples were evaluated for the presence of mouse AFP in this assay. All samples measured ranged from 84.4 to 290.7 ng/mL with an average of 197.9 ng/mL.

G. SPECIFICITY


The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity or interference.

Recombinant mouse	Recombinant human
Leptin	α -Fetoprotein
Prolactin	Natural proteins
Resistin	human α -Fetoprotein
Vitamin D BP	mouse serum albumin

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

pg/mL	OD	Average	Corrected
0	0.080	0.081	-
	0.082		
62.5	0.162	0.165	0.084
	0.168		
125	0.264	0.275	0.194
	0.286		
250	0.493	0.513	0.432
	0.532		
500	0.900	0.945	0.864
	0.989		
1000	1.642	1.700	1.619
	1.757		
2000	2.517	2.585	2.504
	2.653		
4000	3.003	3.029	2.948
	3.055		

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse AFP Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against mouse AFP.	1 plate
Mouse AFP Standard	Recombinant mouse AFP in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	2 vials
Mouse AFP Detection Antibody	Biotinylated AFP antibody, lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Reagent Diluent Concentrate (10×)	A 10× concentrated buffered diluent used to dilute standard, samples, Detection Antibody and HRP.	1 vial
Streptavidin-HRP A (200×)	200× Streptavidin conjugated to horseradish peroxidase .	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant with preservatives.	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution.	1 vial
Stop Solution	2 N sulfuric acid.	1 vial
Plate Sealers	Adhesive strip.	3 strips

B. STORAGE

Unopened Kit	Store at 2-8 °C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Streptavidin-HRP A	May be stored for up to 1 month at 2-8 °C.*
	Wash Buffer (1×)	
	TMB Substrate	
	Stop Solution	
	Standard	Prepare fresh for each assay. Standards may be stored for up to 1 month at -20 °C.*
	Detection Antibody	Aliquot and store for up to 1 month at -20 °C in a manual defrost freezer. *
	Reagent Diluent Concentrate (10×)	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Reagent Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8 °C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ♦ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ♦ Pipettes and pipette tips.
- ♦ Deionized or distilled water.
- ♦ Squirt bottle, manifold dispenser, or automated microplate washer.
- ♦ Horizontal orbital microplate shaker capable of maintaining a speed of 500±50 rpm.
- ♦ 500 mL graduated cylinder.

D. PRECAUTION

- ♦ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ♦ The Stop Solution provided with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Reagent Diluent (1×).

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes at room temperature before centrifugation for 15 minutes at $1000 \times g$. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Samples may require dilution with Reagent Diluent (1×).

B. SAMPLE PREPARATION

Mouse serum samples recommend a 100-fold dilution. For example, add 10 μL of sample into a tube with 90 μL Reagent Diluent (1×) to prepare a 10-fold diluted sample. Mix through and then pipette 10 μL of prepared 10-fold diluted sample into a tube with 90 μL Reagent Diluent (1×) to prepare a final 100-fold diluted sample. Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Note: *Bring all reagents to room temperature before use.*

Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Reagent Diluent (1×) - Use deionized or distilled water to prepare Reagent Diluent (1×).


Detection Antibody (1×) - **Centrifuge briefly before opening. Reconstitution Volume refer to vial label to prepare Detection Antibody (100×).** Allow the Detection Antibody to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Aliquot and store if needed. Dilute to Detection Antibody (1×) with Reagent Diluent (1×). Prepare at least 15 minutes prior to use.

Streptavidin-HRP A (1×) - **Centrifuge briefly before opening.** Dilute to the working concentration specified on the vial label using Reagent Diluent (1×).

Mouse AFP Standard - Centrifuge briefly before opening. Refer to the vial label for the reconstitution volume*. This reconstitution produces a stock solution of 8000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 300 μ L of the Reagent Diluent (1 $\times$) into each tube. Use the stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The 4000 pg/mL standard serves as the high standard. The Reagent Diluent (1 $\times$) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Note: *Bring all reagents and samples to room temperature before use. It is recommended that all samples and standards be assayed in duplicate.*

1. Prepare all reagents and working standards as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 100 μ L of standard and prepared sample per well. Cover with the adhesive strip provided. Incubate for 2 hours at room temperature **on a horizontal orbital microplate shaker set at 500 $\pm$ 50 rpm**. A plate layout is provided for a record of standards and samples assayed. (Samples may require dilution. See Sample Preparation section.)
4. Aspirate each well and wash, repeating the process three times for a total of four washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
5. Add 100 μ L of the Detection Antibody (1 $\times$) in Reagent Diluent (1 $\times$), to each well. Cover with a new adhesive strip and incubate for 2 hours at room temperature **on a horizontal orbital microplate shaker set at 500 $\pm$ 50 rpm**.
6. Repeat the aspiration/wash as in step 4.
7. Add 100 μ L of the working dilution of Streptavidin-HRP A to each well. Cover the plate and incubate for 20 minutes at room temperature **on a horizontal orbital microplate shaker set at 500 $\pm$ 50 rpm. Protect from light**.
8. Repeat the aspiration/wash as in step 4.
9. Add 100 μ L of TMB Substrate to each well. Incubate for 20 minutes at room temperature **on a horizontal orbital microplate shaker set at 500 $\pm$ 50 rpm. Protect from light**.
10. Add 50 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
11. Determine the optical density of each well within 10 minutes, using a microplate

reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

12. CALCULATION OF RESULTS

Average the duplicate readings for each standard and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse AFP concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Gibbs, P.E. et al. (1987) *Biochemistry* 26:1332.
2. Mizejewski, G.J. (2001) *Exp. Biol. Med.* 226:377.
3. Mizejewski, G.J. (2004) *Exp. Biol. Med.* 229:439.
4. Terentiev, A.A. and N.T. Moldogazieva (2013) *Tumor Biol.* 34:2075.
5. Gitlin, D. et al. (1972) *Cancer Res.* 32:979.
6. Andrews, G.K. et al. (1982) *J. Biol. Chem.* 257:5148.
7. Tilghman, S.M. and A. Belayew (1982) *Proc. Natl. Acad. Sci. USA* 79:5254.
8. Jones, E.A. et al. (2001) *J. Anat.* 198:555.
9. Gualdi, R. et al. (1996) *Genes Dev.* 10:1670.
10. Mizejewski, G.J. (2003) *Obstet. Gynecol. Surv.* 58:804.
11. Nunez, E.A. et al. (1989) *J. Nucl. Med. Allied Sci.* 33:18.
12. Uversky, V.N. et al. (1997) *Biochemistry* 36:13638.
13. Zubkova, E. et al. (2013) *Mol. Cell. Biochem.* 379:283.
14. Mizejewski, G.J. (2013) *Tumor Biol.* 34:1317.
15. Gabant, P. et al. (2002) *Proc. Natl. Acad. Sci. USA* 99:12865.
16. De Mees, C. et al. (2006) *Mol. Cell. Biol.* 26:2012.
17. Patni, D. et al. (2012) *J. Pediatr. Endocrinol. Metab.* 25:761.
18. Yamagiwa, S. et al. (2014) *Hepatol Res.* [Epub ahead of print].
19. Jelliffe-Pawlowski, L. et al. (2011) *Prenat. Diagn.* 31:466.
20. Krantz, D.A. et al. (2010) *Clin. Lab. Med.* 30:721.
21. Cuckle, H.S. et al. (1984) *Lancet* 1:926.
22. Donalson, K. et al. (2013) *Prenat. Diagn.* 33:457.
23. Brock, D.J. et al. (1984) *Lancet* 1:767.
24. Canick, J.A. et al. (2003) *Clin. Lab. Med.* 23:385.
25. Abelev, G.I. et al. (1963) *Transplantation* 1:174.
26. Peterson, M.L. et al. (2011) *Semin. Cancer Biol.* 21:21.
27. Alpert, M.E. and R.A. Seeler (1970) *J. Pediatr.* 77:1058.
28. Watanabe, K. (2013) *Int. J. Clin. Oncol.* 18:955.
29. Smith, J.B. and R.T. O'Neill (1971) *Am. J. Med.* 51:767.
30. Billmire, D.F. et al. (2014) *J. Clin. Oncol.* 32:465.
31. Hirajima, S. et al. (2013) *World J. Gastroenterol.* 19:6055.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
	A	B	C	D	E	F	G	H

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