PRODUCT INFORMATION & MANUAL

Mouse/Rat Cystatin C Valukine™ ELISA

Catalog Number: VAL650

For the quantitative determination of natural and recombinant mouse
or rat Cystatin C concentrations

For research use only.
Not for diagnostic or therapeutic procedures.

Bio-Techne China Co. Ltd

P: +86 (21) 52380373 P: 8009881270 F: +86 (21) 52381001
info.cn@bio-techne.com

Please refer to the kit label for expiry date.
Novus kits are guaranteed for 3 months from date of receipt

Version 202411.1

TABLE OF CONTENTS

I.	BACKGROUND	2
II.	OVERVIEW	3
III.	ADVANTAGES	4
IV.	EXPERIMENT	10
V.	KIT COMPONENTS AND STORAGE	11
VI.	PREPARATION	14
VII.	ASSAY PROCEDURE	17
VIII.	REFERENCES	19

I. BACKGROUND

Cystatin C (gene name CST3) is a secreted, extracellular cysteine protease inhibitor that belongs to the cystatin superfamily (1-4). It is a protein of 120 amino acids (aa) and approximately 13 kDa in its non-glycosylated form; a glycosylated form is reported in rat, but not in mouse (2-4). Mouse and rat Cystatin C share 72% aa sequence identity with human Cystatin C and 88% aa sequence identity with each other. Cystatin C is susceptible to endoprotease cleavage producing N-terminally truncated forms (3, 4). Cysteine proteases of the papain family, such as Cathepsins B, H, K, L, and S, are the major targets for Cystatin C (5, 6).

Cystatin C is produced in all tissues and is present in all biological fluids. Cystatin C is freely filtered by the glomeruli. It is then taken up by proximal tubule epithelial cells via megalin-mediated endocytosis and is metabolized so that it does not return to the bloodstream (1, 7-9). Therefore, Cystatin C serum concentration correlates closely to the glomerular filtration rate (GFR). Its measurement in serum or plasma has been proposed as an indicator of drug nephrotoxicity that is less affected by factors such as gender, age, muscle mass, and cirrhosis than creatinine (1, 7, 9). Circulating Cystatin C can, however, be increased during chronic low-level inflammation, in part due to IL-6-mediated increases in Cystatin C production (1). Conversely, the anti-inflammatory cytokines IL-10, IFN- β , and IFN- γ can decrease Cystatin C expression and its circulating levels (10-12).

Cystatin C is involved in several disease processes through its regulation of cysteine protease activity (1). In humans, high circulating Cystatin C in the presence of apparently normal kidney function is an indicator of coronary artery and cardiovascular disease risk (1, 9, 13, 14). In a model of human aortic aneurism, deletion of mouse Cystatin C in ApoE^{-/-} mice promotes inflammation and speeds cathepsin-mediated rupture of the arterial wall tunica elastica (15, 16). Circulating Cystatin C has been reported to influence tumor metastasis. Abnormally low Cystatin C levels allow cathepsin B-mediated degradation of extracellular matrix and promote tumor metastasis, while high Cystatin C levels antagonize TGF- β signaling, slowing cancer invasion and growth (1, 17, 18). Cystatin C is an amyloidogenic protein. In humans, the L68Q variant forms dimers and oligomers more easily than wild type protein under physiological conditions and is the cause for hereditary Cystatin C amyloid angiopathy

(5, 19, 20). Cystatin C also inhibits amyloid- β deposition and protects neuronal cells from toxicity in mouse models of Alzheimer's disease (21-23).

II. OVERVIEW

A. PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for mouse/rat Cystatin C has been pre-coated onto a microplate. Standards, control and samples are pipetted into the wells and any mouse/rat Cystatin C present is bound by the immobilized antibody. After washing away any unbound substances, an enzyme-linked antibody specific for mouse/rat Cystatin C is added to the wells. Following a wash to remove any unbound antibody-enzyme reagent, TMB substrate (Chromogenic agent) is added to the wells and color develops in proportion to the amount of mouse/rat Cystatin C bound in the initial step. The color development is stopped, and the intensity of the color is measured.

B. LIMITATIONS OF THE PROCEDURE

- ◆ **FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**
- ◆ This kit is suitable for cell culture supernates, tissue lysates, mouse/rat serum, mouse/rat plasma and mouse/rat urine.
- ◆ The kit should not be used beyond the expiration date on the kit label.
- ◆ Do not mix or substitute reagents with those from other lots or sources.
- ◆ If samples generate values higher than the highest standard, dilute the samples with Calibrator Diluent (1×) and repeat the assay.
- ◆ Any variation in operator, pipetting technique, washing technique, incubation time or temperature, and kit age can cause variation in binding.

III. ADVANTAGES

A. PRECISION

Intra-assay Precision (Precision within an assay)

Three samples were tested twenty times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays)

Three samples were tested in twenty separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
Mean (pg/mL)	434	894	3065	465	942	3257
Standard Deviation	14.9	29.2	85.9	43.5	51.2	213
CV%	3.4	3.3	2.8	9.4	5.4	6.5

B. RECOVERY

The recovery of mouse/rat Cystatin C spiked into cell culture media was evaluated.

Sample Type	Average % Recovery	Range (%)
Cell culture media (n=4)	107	98-113

C. SENSITIVITY

Forty-one assays were evaluated and the minimum detectable dose (MDD) of mouse/rat Cystatin C ranged from 2.47-12.9 pg/mL. The mean MDD was 3.93 pg/mL.

The MDD was determined by adding two standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

D. CALIBRATION

This immunoassay is calibrated against a highly purified NS0-expressed recombinant mouse Cystatin C produced at R&D Systems.

E. LINEARITY

To assess the linearity of the assay, different samples were containing or spiked with high concentrations of mouse/rat Cystatin C and diluted with Calibrator Diluent (1×) to produce samples with values within the dynamic range of the assay. Samples were diluted prior to assay.

Mouse Sample		Cell culture supernatants (n=4)	Tissue lysates (n=4)	Serum (n=4)	EDTA plasma (n=4)	Heparin plasma (n=4)	Urine (n=4)
1:2	Average % of Expected	98	99	97	97	96	100
	Range (%)	87-101	97-100	96-98	93-100	94-101	99-101
1:4	Average % of Expected	98	97	96	95	95	102
	Range (%)	86-106	94-100	90-100	90-100	89-103	99-107
1:8	Average % of Expected	98	98	96	97	98	104
	Range (%)	88-110	93-103	87-101	92-103	94-108	103-104
1:16	Average % of Expected	99	100	95	101	92	103
	Range (%)	90-111	90-109	89-104	91-110	89-94	98-108

Rat Sample		Cell culture supernates (n=4)	Serum (n=4)	EDTA plasma (n=4)	Heparin plasma (n=4)	Urine (n=4)
1:2	Average % of Expected	92	93	97	91	94
	Range (%)	91-94	81-99	95-98	89-93	90-98
1:4	Average % of Expected	88	91	98	86	90
	Range (%)	84-92	80-100	96-99	83-89	84-95
1:8	Average % of Expected	86	92	98	85	89
	Range (%)	84-88	85-96	95-100	83-86	83-96
1:16	Average % of Expected	80	93	97	85	88
	Range (%)	80-81	87-100	91-101	81-88	83-96

F. SAMPLE VALUES

Mouse/Rat serum/plasma/urine - Mouse and rat samples were evaluated for the presence of mouse and rat Cystatin C in this assay.

Mouse Samples	Mean (ng/mL)	Range (ng/mL)	Standard Deviation (ng/mL)
Serum (n=20)	447	267-638	103
EDTA plasma (n=20)	368	256-515	75.6
Heparin plasma (n=20)	333	254-495	57.6
Urine (n=20)	119	7.45-241	60.4

Rat Samples	Mean (ng/mL)	Range (ng/mL)	Standard Deviation (ng/mL)
Serum (n=20)	2111	1403-2678	360
EDTA plasma (n=20)	1542	963-1972	288
Heparin plasma (n=20)	1535	1002-2054	287
Urine (n=20)	603	239-1383	316

Cell Culture Supernates:

Organs from mice or rats were removed, rinsed in 1× PBS, and kept on ice in 1× PBS. Organs were then cut into 1-2 mm pieces and homogenized using a tissue homogenizer. Cells were seeded into media containing RPMI 1640 supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. Cells were cultured as indicated in the tables below. Aliquots of the cell culture supernates were removed and assayed for levels of mouse/rat Cystatin C.

Mouse Tissue	(ng/mL)
Brain (1 day)	70.1
Heart (3 days)	8.4
Kidney (3 days)	29.0
Liver (3 days)	7.7
Lung (3 days)	12.6
Spleen (3 days)	14.4

Rat Tissue	(ng/mL)
Brain (18 hours)	225
Heart (18 hours)	10.2
Kidney (18 hours)	104
Lung (18 hours)	31.1
Spleen (18 hours)	8.6

3T3-L1 undifferentiated mouse embryonic fibroblast adipose-like cells (2×10^6 cells/T75 flask) were cultured in DMEM supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 100 µg/mL streptomycin sulfate, and incubated for 3 days. The media was removed and 50 mL of fresh media was added and incubated for 4 additional days. An aliquot of the cell culture supernate was removed, assayed for mouse/rat Cystatin C, and measured 166 ng/mL.

3T3-L1 differentiated mouse embryonic fibroblast adipose-like cells (2×10^6 cells/T75 flask) were cultured in DMEM supplemented with 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 100 µg/mL streptomycin sulfate, and incubated for 3 days. The media was removed and 10% fetal bovine serum, 2 mM L-glutamine, 100 U/mL penicillin, 100 µg/mL streptomycin sulfate, 1 µg/mL bovine insulin, 0.5 mM MIX, and 1 µM DEX was added and incubated for 4 additional days. An aliquot of the cell culture supernate was removed, assayed for mouse/rat Cystatin C, and measured 507 ng/mL.

Tissue Lysates - Organs from mice were rinsed with 1 × PBS and homogenized with a tissue homogenizer in 1 × PBS. An equal volume of Cell Lysis Buffer 2 was added and tissues were lysed at room temperature for 30 minutes with gentle agitation. Debris was then removed by centrifugation. An aliquot of each tissue lysate was removed and assayed for levels of mouse Cystatin C.

Mouse Tissue	(ng/mL)
Brain	1254
Heart	210
Kidney	407
Liver	26.7
Lung	169
Spleen	127

G. SPECIFICITY

This assay recognizes natural and recombinant mouse and rat Cystatin C.

The factors listed below were prepared at 80 ng/mL in Calibrator Diluent (1×) and assayed for cross-reactivity. Preparations of the following factors at 80 ng/mL in a mid-range recombinant mouse Cystatin C control were assayed for interference. No significant cross-reactivity or interference was observed.

Recombinant mouse:		Recombinant human:
Cathepsin A	Cathepsin L	Cathepsin O
Cathepsin B	Cathepsin Z	Cathepsin S
Cathepsin C	Cystatin A	Cystatin F
Cathepsin D	Cystatin B	Cystatin S
Cathepsin E	Cystatin E/M	Cystatin SA
Cathepsin H		Cystatin SN


Recombinant human Cystatin C cross-reacts approximately 0.64% in this assay.

Recombinant human Cystatin D cross-reacts approximately 0.11% in this assay.

IV. EXPERIMENT

EXAMPLE STANDARD

The standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

V. KIT COMPONENTS AND STORAGE

A. MATERIALS PROVIDED

Parts	Description	Size
Mouse/Rat Cystatin C Microplate	96 well polystyrene microplate (12 strips of 8 wells) coated with an antibody against mouse/rat Cystatin C.	1 plate
Mouse/Rat Cystatin C Conjugate	An antibody specific for mouse/rat Cystatin C conjugated to horseradish peroxidase.	1 vial
Mouse/Rat Cystatin C Standard	Recombinant mouse/rat Cystatin C in a buffered protein base; lyophilized. Refer to the vial label for reconstitution volume.	1 vial
Mouse/Rat Cystatin C Control	Recombinant mouse/rat Cystatin C in a buffered protein base; lyophilized. The assay value of the control should be within the range specified on the vial label.	1 vial
Assay Diluent RD1W	A buffered protein solution.	1 vial
Calibrator Diluent Concentrate (4×)/RD5-26	A 4× concentrated buffered protein base used to dilute standard and samples.	1 vial
Wash Buffer Concentrate (25×)	A 25× concentrated solution of buffered surfactant.	1 vial
TMB Substrate	TMB ELISA Substrate Solution/TMB Substrate Solution.	2 vials
Stop Solution	Diluted hydrochloric acid.	1 vial
Plate Sealers	Adhesive strip.	3 strips

B. STORAGE

Unopened Kit	Store at 2-8°C. Do not use past kit expiration date.	
Opened/ Reconstituted Reagents	Wash Buffer (1×)	May be stored for up to 1 month at 2-8°C.*
	Assay Diluent RD1W	
	Stop Solution	
	Conjugate	
	TMB Substrate	
	Control	May be stored for up to 1 month at 2-8 °C.*
	Standard	
	Calibrator Diluent Concentrate (4×)/RD5-26	May be stored for up to 1 month at 2-8 °C.* Use and discard diluted Calibrator Diluent (1×). Prepare fresh for each assay.
	Microplate Wells	Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.*

* Provided this is within the expiration date of the kit.

C. OTHER SUPPLIES REQUIRED

- ♦ Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm.
- ♦ Pipettes and pipette tips.
- ♦ Deionized or distilled water.
- ♦ Squirt bottle, manifold dispenser, or automated microplate washer.
- ♦ 100 mL and 500 mL graduated cylinder.
- ♦ Test tubes for dilution of standards and samples.
- ♦ If using tissue lysate samples, the following is also required: Cell Lysis Buffer 2 (R&D Systems, Catalog # 895347).

D. PRECAUTION

- ♦ Cystatin C is detectable in saliva. Take precautionary measures to prevent contamination of kit reagents while running this assay.
- ♦ Some components in this kit contain a preservative which may cause an allergic skin reaction. Avoid breathing mist.
- ♦ The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face, and clothing protection when using this material.

VI. PREPARATION

A. SAMPLE COLLECTION AND STORAGE

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Cell Culture Supernates - Remove particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Tissue Lysates - Cells must be lysed prior to assay as directed in the Sample Values section. Samples may require dilution with Calibrator Diluent (1×).

Serum - Allow blood samples to clot for 2 hours at room temperature before centrifuging for 20 minutes at $2000 \times g$. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Plasma - Collect plasma using EDTA or heparin as an anticoagulant. Centrifuge for 20 minutes at $2000 \times g$ within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Samples may require dilution with Calibrator Diluent (1×).

Note: Citrate plasma has not been validated for use in this assay. Grossly lipemic samples are not suitable for use in this assay.

Urine - Collect urine using a metabolic cage. Remove any particulates by centrifugation and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles. Centrifuge again before assaying to remove any additional precipitates that may appear after storage. Samples may require dilution with Calibrator Diluent (1×).

B. SAMPLE PREPARATION

Mouse serum and plasma samples, and rat urine samples recommend a 200-fold dilution into Calibrator Diluent (1×). A suggested 200-fold dilution is 10 μ L of sample + 90 μ L of Calibrator Diluent (1×). Complete the 200-fold dilution by adding 10 μ L of the diluted sample + 190 μ L of Calibrator Diluent (1×). Optimal dilutions should be determined by the end user.

Rat serum and plasma samples recommend a 400-fold dilution into Calibrator Diluent

(1×). A suggested 400-fold dilution is 10 µL of sample + 90 µL of Calibrator Diluent (1×). Complete the 400-fold dilution by adding 10 µL of the diluted sample + 390 µL of Calibrator Diluent (1×). Optimal dilutions should be determined by the end user.

Mouse urine samples recommend a 40-fold dilution into Calibrator Diluent (1×). A suggested 40-fold dilution is 10 µL of sample + 390 µL of Calibrator Diluent (1×). Optimal dilutions should be determined by the end user.

C. REAGENT PREPARATION

Bring all reagents to room temperature before use.

Note: *Cystatin C is detectable in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.*

Mouse/Rat Cystatin C Control - Reconstitute the control with 1.0 mL of deionized or distilled water. Mix thoroughly. Assay the control undiluted.


Wash Buffer (1×) - If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Dilute 20 mL of Wash Buffer Concentrate (25×) into deionized or distilled water to prepare 500 mL of Wash Buffer (1×).

Calibrator Diluent (1×) - Use deionized or distilled water to prepare Calibrator Diluent (1×).

Mouse/Rat Cystatin C Standard- Refer to the vial label for the reconstitution volume* Reconstitute the Mouse/Rat Cystatin C Standard with Calibrator Diluent (1×). This reconstitution produces a stock solution of 8000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions.

*If you have any question, please seek help from our Technical Support.

Pipette 200 µL of Calibrator Diluent (1×) into each tube. Use the standard stock solution to produce a dilution series (below). Mix each tube thoroughly before the next transfer. The undiluted Mouse/Rat Cystatin C Standard (8000 pg/mL) serves as the high standard. The Calibrator Diluent (1×) serves as the zero standard (0 pg/mL).

D. TECHNICAL HINTS

- When mixing or reconstituting protein solutions, always avoid foaming.
- To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.
- It is recommended that the samples be pipetted within 15 minutes.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- TMB Substrate should remain colorless until added to the plate. Keep TMB Substrate protected from light. TMB Substrate should change from colorless to gradations of blue.
- Stop Solution should be added to the plate in the same order as the TMB Substrate. The color developed in the wells will turn from blue to yellow upon addition of the Stop Solution. Wells that are green in color indicate that the Stop Solution has not mixed thoroughly with the TMB Substrate.

VII. ASSAY PROCEDURE

Bring all reagents and samples to room temperature before use. It is recommended that all samples, control and standards be assayed in duplicate.

Note: *Cystatin C is detectable in saliva. It is recommended that a face mask and gloves be used to protect kit reagents from contamination.*

1. Prepare all reagents, working standards, control and samples as directed in the previous sections.
2. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal.
3. Add 50 μ L of Assay Diluent RD1W to each well.
4. Add 50 μ L of standard, control and prepared sample per well. Cover with the adhesive strip provided. **Incubate for 2 hours at room temperature.** A plate layout is provided for a record of standards and samples assayed.
5. Aspirate each well and wash, repeating the process four times for a total of five washes. Wash by filling each well with Wash Buffer (400 μ L) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
6. Add 100 μ L of Mouse/Rat Cystatin C Conjugate to each well. Cover with a new adhesive strip. **Incubate for 2 hours at room temperature.**
7. Repeat the aspiration/wash as in step 5.
8. Add 100 μ L of TMB Substrate to each well. **Incubate for 30 minutes at room temperature. Protect from light.**
9. Add 100 μ L of Stop Solution to each well. Gently tap the plate to ensure thorough mixing.
10. Determine the optical density of each well within 10 minutes, using a microplate reader set to 450 nm. If wavelength correction is available, set to 540 nm or 570 nm. If wavelength correction is not available, subtract readings at 540 nm or 570 nm from the readings at 450 nm. This subtraction will correct for optical imperfections in the plate. Readings made directly at 450 nm without correction may be higher and less accurate.

11. CALCULATION OF RESULTS

Average the duplicate readings for each standard, control and sample and subtract the average zero standard optical density (O.D.). Create a standard curve by reducing the data using computer software capable of generating a four-parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse/rat Cystatin C concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.

If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

VIII. REFERENCES

1. Seronie-Vivien, S. *et al.* (2008) Clin. Chem. Lab. Med. 46:1664.
2. Solem, M. *et al.* (1990) Biochem. Biophys. Res. Commun. 172:945.
3. Esnard, A. *et al.* (1988) FEBS Lett. 236:475.
4. Cole, T. *et al.* (1989) Eur. J. Biochem. 186:34.
5. Janowski, R. *et al.* (2001) Nat. Struct. Biol. 8:316.
6. Hall, A. *et al.* (1995) J. Bio. Chem. 270:5115.
7. Laterza, O.F. *et al.* (2002) Clin. Chem. 48:699.
8. Kaseda, R. *et al.* (2007) Biochem. Biophys. Res. Commun. 357:1130.
9. Taglieri, N. *et al.* (2009) Clin. Chem. 55:1932.
10. Xu, Y. *et al.* (2011) J. Immunol. 186:3666.
11. Staun-Ram, E. and A. Miller (2011) J. Neuroimmunol. 232:200.
12. Frendeus, K.H. *et al.* (2009) Int. J. Biochem. Cell Biol. 41:2262.
13. Kiyosue, A. *et al.* (2010) Circ. J. 74:2441.
14. Xie, L. *et al.* (2010) Cardiovasc. Res. 87:628.
15. Sukhova, G.K. *et al.* (2005) Circ. Res. 96:368.
16. Schulte, S. *et al.* (2010) Am. J. Pathol. 177:456.
17. Sokol, J.P. and W.P. Schiemann (2004) Mol. Cancer Res. 2:183.
18. Yu, W. *et al.* (2010) PLoS ONE 5:e13973.
19. Abrahamson, M. *et al.* (1992) Hum. Genet. 89:377.
20. Nagai, A. *et al.* (2008) Front. Biosci. 13:3470.
21. Mi, W. *et al.* (2007) Nat. Genet. 39:1440.
22. Sun, B. *et al.* (2008) Neuron 60:247.
23. Tizon, B. *et al.* (2010) J. Alzheimers Dis. 19:885.

PLATE LAYOUT

Use this plate layout to record standards and samples assayed.

	A	B	C	D	E	F	G	H
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								

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