

RNAscope® 2.5 HD Duplex Detection Kit (GREEN/RED) Quick Guide

For FFPE Tissues

Introduction

This quick guide is intended for advanced users who are familiar with the procedures in the *Formalin-Fixed Paraffin-Embedded (FFPE) Sample Preparation Pretreatment Guide User Manual, Part 1* (Catalog No. 322452-USM) and *RNAscope® 2.5 HD Duplex Detection Kit User Manual, Part 2* (Catalog No. 322500-USM). Refer to the user manual for safety guidelines. For every chemical, read the Safety Data Sheet (SDS) and follow handling instructions. Wear appropriate protective eyewear, clothing, and gloves. For the latest services and support information, go to: www.acdbio.com/support.


Part 1 Prepare and Pretreat Samples

Workflow Steps	
<p>PREPARE FFPE SECTIONS</p>	<ol style="list-style-type: none"> 1. Immediately place dissected tissue sample in fresh 10% NBF for 16–32 HRS at ROOM TEMPERATURE (RT). 2. Dehydrate, embed in paraffin, and cut the sample into 5 +/- 1 µm sections. Mount sections on Superfrost® Plus slides. <hr/> <p>OPTIONAL STOPPING POINT (1). Use sectioned tissue within 3 months. Store sections with desiccants at RT.</p>
<p>PREPARE SLIDES ~1.5 HOURS</p> <p>Bake Slides ↓ Deparaffinize FFPE Sections</p>	<p>Bake Slides</p> <ol style="list-style-type: none"> 1. Bake slides in a dry oven for 1 HR at 60°C. <hr/> <p>OPTIONAL STOPPING POINT (2). Use sectioned tissue within 1 week. Store sections with dessicants at RT.</p> <p>Deparaffinize FFPE Sections</p> <ol style="list-style-type: none"> 1. In a fume hood: <ul style="list-style-type: none"> • Fill two Tissue-Tek® Clearing Agent dishes with ~200 mL fresh xylene. • Fill two Tissue-Tek® Staining dishes with ~200 mL fresh 100% EtOH. 2. Place slides in a Tissue-Tek® Slide Rack in xylene 2 x 5 MIN. 3. Incubate slides in 100% EtOH 2 x 1 MIN. 4. Remove slides from the rack. Air dry slides for 5 MIN at RT. <hr/> <p>OPTIONAL STOPPING POINT (3). Air dry overnight at RT (must use within 24 hrs) or proceed directly to the next step.</p>
<p>PRETREAT SAMPLES ~1–2 HOURS</p> <p>Prepare Oven and Reagents ↓ Apply RNAscope® Hydrogen Peroxide ↓</p>	<p>Prepare Oven and Reagents (30 MIN at 40°C)</p> <ol style="list-style-type: none"> 1. Set HybEZ™ oven to 40°C and warm HybEZ™ Humidity Control Tray containing wet Humidifying Paper for 30 MIN before use. Keep tray warm during the assay. 2. Prepare 700 mL fresh 1X Target Retrieval in a beaker. Cover with foil, bring to a mild boil, and maintain. Do not boil more than 30 MIN before use. <p>Apply Hydrogen Peroxide (10 MIN at RT)</p> <ol style="list-style-type: none"> 1. Add ~5–8 drops of Hydrogen Peroxide to each section for 10 MIN at RT. 2. Place slides into a Tissue-Tek® Slide Rack submerged in distilled water. 3. Wash slides in the distilled water by moving the rack up and down 3–5 times and repeat with fresh distilled water.

<p>↓ Perform Target Retrieval ↓ Create Hydrophobic Barrier ↓ Apply Protease Plus</p>	<p>Perform RNAscope® Target Retrieval</p> <ol style="list-style-type: none"> 1. With a pair of forceps <i>very slowly</i> submerge the slide rack into the boiling 1X Target Retrieval solution. Refer to Appendix A of the <i>Formalin-Fixed Paraffin-Embedded (FFPE) Sample Preparation Pretreatment Guide User Manual, Part 1</i> (Cat. No. 322452) for specific pretreatment time, depending on your tissue type. 2. <i>Immediately</i> transfer hot slide rack to a staining dish containing distilled water. 3. Wash slides in the distilled water by moving the rack up and down 3–5 times and repeat with fresh distilled water. 4. Wash slides in fresh 100% EtOH by moving the rack up and down 3–5 times and air dry. <p>Create Hydrophobic Barrier</p> <ol style="list-style-type: none"> 1. Draw 2–4 times around tissue using the Immedge™ hydrophobic barrier pen. Dry completely ~2 MIN or OVERNIGHT at RT. <p>Apply Protease Plus</p> <ol style="list-style-type: none"> 1. Place slides in the HybEZ™ Slide Rack, and add ~5 drops of Protease Plus to each section. 2. Place the HybEZ™ Slide Rack in the prewarmed HybEZ™ Humidity Control Tray. Seal tray and insert back into the HybEZ™ Oven. Incubate at 40°C for 30 MIN. <p>Note: If needed, prepare RNAscope® 2.5 assay materials during this step.</p> <ol style="list-style-type: none"> 3. Wash slides in the distilled water by moving the rack up and down 3–5 times and repeat with fresh distilled water.
--	---

Part 2: RNAscope® 2.5 Duplex Assay

Workflow Steps	
<p>PREPARE THE MATERIALS ~10–30 MIN</p>	<ol style="list-style-type: none"> 1. Prepare 3 L of 1X Wash Buffer by adding 2.94 L distilled water and 1 bottle (60 mL) of 50X Wash Buffer to a large carboy. Mix well. Warm 50X Wash Buffer up to 40°C for 10–20 MIN before making 1X Wash Buffer. 2. Prepare 50% Hematoxylin and 0.02% Ammonia water. 3. Equilibrate remove Amps 1–10 from the refrigerator and keep at RT <p>Prepare Probes</p> <ol style="list-style-type: none"> 4. Warm probes for 10 MIN at 40°C, then cool to RT. 5. Briefly spin the C2 probe. 6. Mix 1:50 ratio of C2 probe to C1 probe by pipetting 1 volume of C2 probe to 50 volumes of C1 probe into a tube. Invert the tube several times. 7. Mixed probes can be stored at 4°C for up to 6 months.
<p>RUN THE ASSAY ~7 Hours</p> <p>Hybridize Probe ↓ Hybridize Amp 1 ↓ Hybridize Amp 2 ↓ Hybridize Amp 3 ↓ Hybridize Amp 4</p>	<p>Hybridize Probe (2 HRS at 40°C)</p> <ol style="list-style-type: none"> 1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops probe to each section. 2. Insert the sealed tray containing HybEZ™ Slide Rack back into the HybEZ™ Oven for 2 HRS at 40°C. Remove slide rack. 3. Wash slides in 1X Wash Buffer for 2 MIN at RT. Repeat with fresh 1X Wash Buffer. <p>Hybridize Amp 1 (30 MIN at 40°C)</p> <ol style="list-style-type: none"> 1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 1 to each section. 2. Insert the sealed tray containing HybEZ™ Slide Rack into the HybEZ™ Oven for 30 MIN at 40°C. Remove slide rack. 3. Wash slides in 1X Wash Buffer for 2 MIN at RT. Repeat with fresh 1X Wash Buffer.



Hybridize Amp 2 (15 MIN at 40°C)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 2 to each section.
2. Insert the sealed tray containing HybEZ™ Slide Rack into the HybEZ™ Oven for **15 MIN** at **40°C**. Remove slide rack.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 3 (30 MIN at 40°C)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 3 to each section.
2. Insert the sealed tray containing HybEZ™ Slide Rack into the HybEZ™ Oven for **30 MIN** at **40°C**. Remove slide rack.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 4 (15 MIN at 40°C)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 4 to each section.
2. Insert the sealed tray containing HybEZ™ Slide Rack into the HybEZ™ Oven for **15 MIN** at **40°C**. Remove slide rack, but do *not* place the tray back into the oven.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 5 (30 MIN at RT)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 5 to each section.
2. Incubate the sealed tray containing HybEZ™ Slide Rack for **30 MIN** at **RT**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 6 (15 MIN at RT)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 6 to each section.
2. Incubate the sealed tray containing HybEZ™ Slide Rack for **15 MIN** at **RT**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Detect the Signal (10 MIN at RT)

1. Briefly spin **RED-B** and mix 1 volume of **RED-B** to 60 volumes of **RED-A** (must use within **3–5 MIN**). For examples add 2.5 µL of **Red-B** to 150 µL of **Red-A** per section.
2. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and pipette ~120 µL of **RED** solution onto each tissue section.
3. Incubate sealed tray containing HybEZ™ Slide Rack for **10 MIN** at **RT**.
4. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 7 (15 MIN at 40°C)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 7 to each section.
2. Incubate the sealed tray containing HybEZ™ Slide Rack for **15 MIN** at **40°C**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 8 (30 MIN at 40°C)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 8 to each section.
2. Incubate the sealed tray containing HybEZ™ Slide Rack for **30 MIN** at **40°C**.
3. Wash slides in 1X Wash Buffer for **2 MIN** at **RT**. Repeat with fresh buffer.

Hybridize Amp 9 (30 MIN at RT)

1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 9 to each section.
2. Incubate the sealed tray containing HybEZ™ Slide Rack for **30 MIN** at **RT**.

	<p>3. Wash slides in 1X Wash Buffer for 2 MIN at RT. Repeat with fresh buffer.</p> <p>Hybridize Amp 10 (15 MIN at RT)</p> <ol style="list-style-type: none"> 1. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and add ~4 drops Amp 10 to each section. 2. Incubate the sealed tray containing HybEZ™ Slide Rack for 15 MIN at RT. 3. Wash slides in 1X Wash Buffer for 2 MIN at RT. Repeat with fresh buffer. <p>Detect the Signal (10 MIN at RT)</p> <ol style="list-style-type: none"> 1. Briefly spin GREEN-B and mix 1 volume of Green-B to 50 volumes of GREEN-A (must use within 3–5 MIN). For examples add 3 µL of GREEN-B to 150 µL of GREEN-A per section. 2. Remove excess liquid from slides, place in the HybEZ™ Slide Rack, and pipette ~120 µL of GREEN solution onto each tissue section. 3. Incubate sealed tray containing HybEZ™ Slide Rack for 10 MIN at RT. 4. Remove solution from slides and wash 5 MIN with 1X wash buffer 5. Rinse quickly with distilled water. <p>Counterstain the Slides (30 Sec at RT)</p> <ol style="list-style-type: none"> 1. Place slides in 50% Hematoxylin I for 30 Sec at RT. Wash 3–5 times in tap water and repeat with fresh tap water. 2. Wash slides 10 SEC in 0.02% Ammonia water, and then wash 3–5 times in tap water. <p>Mount the Slides</p> <ol style="list-style-type: none"> 1. Dry slides in a 60°C dry oven for 15 MIN. 2. Dip the slides into FRESH pure xylene and immediately place 1–2 drops of EcoMount on the slide before the xylene dries. Place a coverslip over the section. 3. Air dry for 30 MIN.
EVALUATE THE RESULTS	Examine tissue sections under a standard bright field microscope at 20–40X magnification

Troubleshooting

For troubleshooting information, please contact technical support at support.acd@bio-techne.com.

For Research Use Only. Not for diagnostic use.

NOTICE TO PURCHASER: PLEASE REFER TO THE RNASCOPE® 2.5 ASSAY-USER MANUAL FOR LIMITED USE LABEL LICENSE OR DISCLAIMER INFORMATION. Advanced Cell Diagnostics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This manual is subject to change without notice. Although this manual has been prepared with every precaution to ensure accuracy, Advanced Cell Diagnostics, Inc. assumes no liability for any errors, omissions, or for any damages resulting from the use of this information.

© 2019 Advanced Cell Diagnostics. All rights reserved. RNAScope® and HybEZ™ are trademarks of Advanced Cell Diagnostics, Inc. All other trademarks belong to their respective owners.

Headquarters

7707 Gateway Blvd Suite 200, Newark, CA 94560 Phone 1-510-576-8800 Toll Free 1-877-576-3636

For support, email support.acd@bio-techne.com

www.acdbio.com

