

Justice Laboratory Software

23 Church Street
PO Box 1227
Denville, New Jersey 07834
Tel. (973) 586-8551
Fax (973) 627-0947
www.ChromPerfect.com

chromperfect Custom Software

Revision Notes, Proteinsimple

...

Subject points from an extracted document 09/15/2016

6.0.12 Release September 15, 2016

- Completed modifications for OEM customers
- Improved Chromperfect browser to include both Raw and Bnd files

6.0.12 ProteinSimple Release Notes 12/12/11 – 9/15/16

- No features from ver6.04 removed
- Bug fix to correctly calculate RMS noise on the baseline
- New Menu bar included, original menu still available to users
- IDE update to Visual Studio 2013
- Targeted .Net Framework 4.0
- Support for Windows 8.1 and Windows 10 added
- Bug fix for printing eGrams while in compare mode in Analysis
- Auto-name conventions improved

- Operator field made available
- RMS noise value added to the Summary applet
- DataGridView replaces spreadsheet control from Grape City
- eGram to Excel included
- Bug fix for zero offset in pI plots
- Evaluation of update sent J. Wu

7.4.0 Proteinsimple Release Notes 2/23/2018

- All revisions through 6.0.12 included
- Targeted .Net Framework 4.5
- Evaluation update sent to Rick Le Bold via link.

7.4.2 Proteinsimple Release Notes 5/18/2018

- All revisions through 7.4.0 included
- Added Method printout to be selected by Tab Form information, (ie Events)
- Added high res PNG capability
- Added high res PNG in Compare mode from Analysis