

# Quantist<sup>™</sup> Release Notes

P/N 805-0004

Revision A, July 2022

Copyright © 2022 Bio-Techne. All rights reserved.

---

## Version 1.0.0

### Change History

<b>Change Type</b>	<b>Notes</b>
--------------------	--------------

Feature	Initial Release
---------	-----------------

### Known Issues

- Well Assignments View can crash application if assignments are performed quickly back-to-back (BTQ-569)
- Loading Luminex CSV File with standard dilutions all containing the same concentration will cause application to become unresponsive (BTQ-572)

brands:

**RD**SYSTEMS

**NOVUS**  
BIOLOGICALS

**TOCRIS**

protein  
simple

**ACD**<sup>™</sup>

**exosome**d<sub>x</sub>