

# Product Datasheet

## Alcian Blue (pH 2.5) Stain Kit H-3501-NB

Unit Size: 1 Kit

Store at room temperature.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H-3501-NB](http://www.novusbio.com/H-3501-NB)

Updated 2/14/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H-3501-NB](http://www.novusbio.com/reviews/destination/H-3501-NB)




**H-3501-NB****Alcian Blue (pH 2.5) Stain Kit**

<b>Product Information</b>	
<b>Unit Size</b>	1 Kit
<b>Concentration</b>	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
<b>Storage</b>	Store at room temperature.
<b>Product Description</b>	
<b>Description</b>	<p>The Alcian Blue (pH 2.5) Stain Kit is intended for use in the histological visualization of sulfated and carboxylated acid mucopolysaccharides and sulfated and carboxylated sialomucins (glycoproteins).</p> <p>Acidic Sulfated Mucosubstances: Blue  Hyaluronic Acid: Blue  Sialomucins: Blue  Nuclei: Red</p> <p>Control Tissue: Small Intestine, Appendix, Colon</p>
<b>Kit Components</b>	Alcian Blue Solution (pH 2.5) - 30 mL, Nuclear Fast Red (Enhanced) - 30 mL, Acetic Acid Solution (3%) - 2 x 30 mL
<b>Notes</b>	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.

## Images

Alcian Blue (pH 2.5) Stain Kit [H-3501-NB] - Colon tissue stained with Alcian Blue (pH 2.5) Stain Kit.



Alcian Blue (pH 2.5) Stain Kit [H-3501-NB]



Alcian Blue (pH 2.5) Stain Kit [H-3501-NB]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H-3501-NB](http://www.novusbio.com/reviews/submit/H-3501-NB)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

