Product Datasheet

Histone H2AX [p Ser139] Antibody [mFluor Violet 610 SE] NB100-384MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB100-384MFV610

Updated 11/8/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NB100-384MFV610



NB100-384MFV610

Histone H2AX [p Ser139] Antibody [mFluor Violet 610 SE]

Product Information		
Unit Size	0.1 ml	
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.	
Storage	Store at 4C in the dark.	
Clonality	Polyclonal	
Preservative	0.05% Sodium Azide	
Isotype	IgG	
Conjugate	mFluor Violet 610 SE	
Purity	Immunogen affinity purified	
Buffer	50mM Sodium Borate	
Product Description		
Host	Rabbit	
Gene ID	3014	
Gene Symbol	H2AX	
Species	Human, Mouse, Rat, Canine	
Reactivity Notes	Rat reactivity reported in scientific literature (PMID: 27102221), Canine reactivity reported in scientific literature (PMID: 23365434).	
Marker	DNA Double-strand break marker	
Specificity/Sensitivity	The epitope maps to a region surrounding phosphorylated serine 139 of human histone H2AX.	
Immunogen	This Histone H2AX [p Ser139] Antibody was developed against to a region surrounding phosphorylated serine 139 of human histone H2AX [Swiss-Prot entry P16104] (GeneID 3014).	
Product Application Details		
Applications	Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Chromatin Immunoprecipitation (ChIP)	
Recommended Dilutions	Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Chromatin Immunoprecipitation (ChIP)	
Application Notes	Optimal dilution of this antibody should be experimentally determined.	

Images

Histone H2AX [p Ser139] Antibody [mFluor Violet 610 SE] [NB100-384MFV610] - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.

	mFluor™∖	/iolet 610
	LASER (nm)	FILTER
^{mFluor[™] Violet 6¹⁰}	Violet (405)	605/30
	EXCITATION MAX (nm)	EMISSION MAX (nm)
CAUTION -Research Lise Only	421	613

www.novusbio.com





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Products Related to NB100-384MFV610

NBP1-42464	Recombinant Human Histone H2AX His Protein
DYC2288-2	Histone H2AX [p Ser139] [Biotin]
NB100-304	53BP1 Antibody - BSA Free
NB100-143	Nbs1 Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-384MFV610

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

