Product Datasheet

FCAR/CD89 Antibody (MIP8a) [mFluor Violet 610 SE] NB100-63266MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

G

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB100-63266MFV610

Updated 7/11/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NB100-63266MFV610



NB100-63266MFV610

FCAR/CD89 Antibody (MIP8a) [mFluor Violet 610 SE]

-
0.1 ml
Please see the vial label for concentration. If unlisted please contact technical services.
Store at 4C in the dark.
Monoclonal
MIP8a
0.05% Sodium Azide
lgG1
mFluor Violet 610 SE
Protein A or G purified
50mM Sodium Borate
Mouse
2204
FCAR
Human
NB100-63266 recognizes the human CD89 cell surface antigen, a 50-75kD cell surface glycoprotein that is also known as the IgA receptor (Fc alpha R). CD89 is expressed by peripheral blood monocytes and neutrophils. Clone MIP8a blocks binding of IgA to the Fc alpha R, and also inhibits neutrophil phagocytosis of IgA complexes.
Recombinant soluble human Fc alpha R
mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Western Blot, ELISA, Flow Cytometry, Immunoprecipitation

Applications	Western Blot, ELISA, Flow Cytometry, Immunoprecipitation
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunoprecipitation
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

FCAR/CD89 Antibody (MIP8a) [mFluor Violet 610 SE] [NB100-63266MFV610] - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.

mFluor [™] Violet 61 MFluor [™] Violet 610 Excitation Max (mm) EMISSION MAX (m)			
hFluor" Violet (405) 605/30		mFluor™∖	/iolet 610
runt. Alloier .	-	LASER (nm)	FILTER
EXCITATION MAX (nm) EMISSION MAX (n	'luor'' Violet 6 ¹⁰	Violet (405)	605/30
		EXCITATION MAX (nm)	EMISSION MAX (nm)
421 613		421	613

www.novusbio.com





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: technical@novusbio.com Orders: orders@novusbio.com General: novus@novusbio.com

Products Related to NB100-63266MFV610

210-TA-005TNF-alpha [Unconjugated]3939-FA-050FCAR/CD89 [Unconjugated]D6050IL-6 [HRP]7954-GM-010/CFGM-CSF [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-63266MFV610

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

