Product Datasheet

Integrin alpha 2b/CD41 Antibody (IVA30) [mFluor Violet 500 SE] NB500-496MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB500-496MFV500

Updated 9/20/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB500-496MFV500



NB500-496MFV500

Integrin alpha 2b/CD41 Antibody (IVA30) [mFluor Violet 500 SE]

Product Information		
Unit Size	0.1 ml	
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.	
Storage	Store at 4C in the dark.	
Clonality	Monoclonal	
Clone	IVA30	
Preservative	0.05% Sodium Azide	
Isotype	IgG2	
Conjugate	mFluor Violet 500 SE	
Purity	Protein A purified	
Buffer	50mM Sodium Borate	
Product Description		
Host	Mouse	
Gene ID	3674	
Gene Symbol	ITGA2B	
Species	Human, Bovine	
Specificity/Sensitivity	The antibody IVA30 reacts with CD41/CD61 complex (bovine), one of the earliest markers of the megakaryocytic lineage. CD41/CD61 complex acts as the receptor for fibrinogen, fibronectin, Von Willebrand factor, vitronectin, thrombin and mediates platelets aggregation. CD41/CD61 is expressed only by platelets and megakaryocytes; the complex may be absent or strongly reduced in Glanzmann's thrombasthenia (GT).	
Immunogen	Bovine trombocytes.	
Product Application Details		
Applications	ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, CyTOF-ready	
Recommended Dilutions	Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, CyTOF-ready	
Application Notes	Optimal dilution of this antibody should be experimentally determined.	

Images

Integrin alpha 2b/CD41 Antibody (IVA30) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.

	mFluor™∖	/iolet 500
	LASER (nm)	FILTER
^{INFluor^w Violet ⁵⁰⁰}	Violet (405)	515/20
	EXCITATION MAX (nm)	EMISSION MAX (nm)
ENOUS Notechne trunk	410	501
CAUTION - Research Use Onl		

www.novusbio.com





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Products Related to NB500-496MFV500

n

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB500-496MFV500

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

