## **Product Datasheet**

### Clusterin/APOJ Antibody (Hs-3) [mFluor Violet 450 SE] NB500-552MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

G

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB500-552MFV450

Updated 9/20/2023 v.20.1

# Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NB500-552MFV450



#### NB500-552MFV450

Clusterin/APOJ Antibody (Hs-3) [mFluor Violet 450 SE]

Product Information			
Unit Size	0.1 ml		
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.		
Storage	Store at 4C in the dark.		
Clonality	Monoclonal		
Clone	Hs-3		
Preservative	0.05% Sodium Azide		
Isotype	lgG1		
Conjugate	mFluor Violet 450 SE		
Purity	Protein A purified		
Buffer	50mM Sodium Borate		
Product Description			
Host	Mouse		
Gene ID	1191		
Gene Symbol	CLU		
Species	Human, Bovine (Negative), Canine (Negative), Feline (Negative), Porcine (Negative)		
Specificity/Sensitivity	Apolipoprotein J (Hs-3)		
Immunogen	Freshly ejaculated human sperms were washed in PBS and extracted in 3% acetic acid, 10% glycerol, 30 mM benzaminidine. The acid extract was dialyzed against 0.2% acetic acid and subsequently used for immunization.		
Product Application Details			
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin		
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin		
Application Notes	Optimal dilution of this antibody should be experimentally determined.		

#### Images

Clusterin/APOJ Antibody (Hs-3) [mFluor Violet 450 SE] [NB500-552MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.

nd		mFluor™ Violet 450		
		LASER (nm)	FILTER	
	mFluor" Violet 450	Violet (405)	450/45	
		EXCITATION MAX (nm)	EMISSION MAX (nm)	
		406	445	
	CAUTOR - Research Use Only			





#### Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

#### Products Related to NB500-552MFV450

NBP1-90128PEP	Clusterin/APOJ Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
DCLU00	Clusterin/APOJ [HRP]
D6050	IL-6 [HRP]

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB500-552MFV450

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

