Product Datasheet

Clusterin/APOJ Antibody (Hs-3) [mFluor Violet 610 SE] NB500-552MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

G

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB500-552MFV610

Updated 9/20/2023 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NB500-552MFV610



NB500-552MFV610

Clusterin/APOJ Antibody (Hs-3) [mFluor Violet 610 SE]

| Product Information | | | |
|-----------------------------|---|--|--|
| Unit Size | 0.1 ml | | |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. | | |
| Storage | Store at 4C in the dark. | | |
| Clonality | Monoclonal | | |
| Clone | Hs-3 | | |
| Preservative | 0.05% Sodium Azide | | |
| Isotype | IgG1 | | |
| Conjugate | mFluor Violet 610 SE | | |
| Purity | Protein A purified | | |
| Buffer | 50mM Sodium Borate | | |
| Product Description | | | |
| Host | Mouse | | |
| Gene ID | 1191 | | |
| Gene Symbol | CLU | | |
| Species | Human, Bovine (Negative), Canine (Negative), Feline (Negative), Porcine (Negative) | | |
| Specificity/Sensitivity | Apolipoprotein J (Hs-3) | | |
| Immunogen | Freshly ejaculated human sperms were washed in PBS and extracted in 3% acetic acid, 10% glycerol, 30 mM benzaminidine. The acid extract was dialyzed against 0.2% acetic acid and subsequently used for immunization. | | |
| Product Application Details | | | |
| Applications | Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin | | |
| Recommended Dilutions | Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin | | |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. | | |
| | | | |

Images

Clusterin/APOJ Antibody (Hs-3) [mFluor Violet 610 SE] [NB500-552MFV610] - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.

| nd | | mFluor™ Violet 610 | | |
|----|-----------------------------|---------------------|-------------------|--|
| | and the second | LASER (nm) | FILTER | |
| | mFluor" Violet 610 | Violet (405) | 605/30 | |
| | | EXCITATION MAX (nm) | EMISSION MAX (nm) | |
| | | 421 | 613 | |
| | CAUTION -Research Lise Only | | | |





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Products Related to NB500-552MFV610

| NBP1-90128PEP | Clusterin/APOJ Recombinant Protein Antigen |
|---------------|--|
| 210-TA-005 | TNF-alpha [Unconjugated] |
| DCLU00 | Clusterin/APOJ [HRP] |
| D6050 | IL-6 [HRP] |
| | |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB500-552MFV610

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

