## **Product Datasheet**

### AP2 alpha Antibody (A6/2/2) [mFluor Violet 500 SE] NB600-202MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB600-202MFV500

Updated 7/24/2023 v.20.1

# Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NB600-202MFV500



#### NB600-202MFV500

AP2 alpha Antibody (A6/2/2) [mFluor Violet 500 SE]

Product Information		
Unit Size	0.1 ml	
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.	
Storage	Store at 4C in the dark.	
Clonality	Monoclonal	
Clone	A6/2/2	
Preservative	0.05% Sodium Azide	
Isotype	IgG1	
Conjugate	mFluor Violet 500 SE	
Purity	Protein A purified	
Buffer	50mM Sodium Borate	
Product Description		
Host	Mouse	
Gene ID	7020	
Gene Symbol	TFAP2A	
Species	Human	
Specificity/Sensitivity	This antibody detects both the alpha and beta forms of AP2.	
Immunogen	Peptide to C-Terminal of AP2 alpha coupled to KLH. PNSHTDNNAKSSDKEEKHRK	
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.	
Product Application Details		
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Radioimmunoassay, CyTOF-ready	
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Radioimmunoassay, CyTOF-ready	
Application Notes	Optimal dilution of this antibody should be experimentally determined.	

#### Images

AP2 alpha Antibody (A6/2/2) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.

	mFluor™∖	/iolet 500
	LASER (nm)	FILTER
mFluor" Violet 500	Violet (405)	515/20
	EXCITATION MAX (nm)	EMISSION MAX (nm)
	410	501
CAUTOR -Research Use On		

www.novusbio.com





#### **Novus Biologicals USA**

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

#### Products Related to NB600-202MFV500

H00007020-Q01-10ugRecombinant Human AP2 alpha GST (N-Term) Protein210-TA-005TNF-alpha [Unconjugated]H00007020-P01-10ugRecombinant Human AP2 alpha GST (N-Term) Protein1129-ER-050ErbB2/Her2 [Unconjugated]

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB600-202MFV500

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

