

Product Datasheet

HPI3 Antibody (HIC1-7H10) [mFluor Violet 450 SE] NBP1-18947MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-18947MFV450

Updated 9/20/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-18947MFV450



NBP1-18947MFV450

HPi3 Antibody (HIC1-7H10) [mFluor Violet 450 SE]

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	HIC1-7H10
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	mFluor Violet 450 SE
Purity	Tissue culture supernatant
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Species	Human, Mouse
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 26515645).
Marker	Endocrine pan-islet cell Marker
Specificity/Sensitivity	Specific for multiple endocrine cell types.
Immunogen	Human pancreatic enriched islet cells containing low levels of exocrine and ductal cells.
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.

Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Frozen
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

HPi3 Antibody (HIC1-7H10) [mFluor Violet 450 SE] [NBP1-18947MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.



mFluor™ Violet 450

LASER (nm)	FILTER
Violet (405)	450/45
EXCITATION MAX (nm)	EMISSION MAX (nm)
406	445



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-18947MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

