Product Datasheet

Uterus Tissue Slides (Leiomyosarcoma) NBP2-30327

Unit Size: 5 Slides

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-30327

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications
Submit a review at www.novusbio.com/reviews/destination/NBP2-30327



NBP2-30327

Lysate Tissue Condition

Uterus Tissue Slides (Leiomyosarcoma)

Oterus Tissue Sildes (Leiomyosarcoma)	
Product Information	
Unit Size	5 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C. Do not freeze.
Reconstitution Instructions	These slides are paraffin coated to prevent sample oxidization, it is recommended that slides are first de-paraffinized by baking at 60 degrees C for 1 hour in a vertical orientation prior to performing antigen retrieval procedures.
Product Description	
Description	Layout: 1X1 Diameter: 4 Thickness:Single Tissue Slides Pos:A01 SI No:1 Age:27 Sex:F Organ:Uterus Pathology:Leiomyosarcoma Species:Human
Species	Human
Specificity/Sensitivity	Uterus Leiomyosarcoma Format:Formalin Fixed
Lysate Type	Tissue

Leiomyosarcoma





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP2-30327

NBP2-78095 Human Uterus Tissue MicroArray (Tumor)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-30327

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

