

# Product Datasheet

## Myeloid Cell Marker Antibody (BM-1) [mFluor Violet 610 SE]

**NBP2-34562MFV610**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-34562MFV610](http://www.novusbio.com/NBP2-34562MFV610)

Updated 3/28/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-34562MFV610](http://www.novusbio.com/reviews/destination/NBP2-34562MFV610)



**NBP2-34562MFV610**

Myeloid Cell Marker Antibody (BM-1) [mFluor Violet 610 SE]

| <b>Product Information</b>         |  |
|------------------------------------|--|
| <b>Unit Size</b>                   | 0.1 ml   |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| <b>Storage</b>                     | Store at 4C in the dark.   |
| <b>Clonality</b>                   | Monoclonal   |
| <b>Clone</b>                       | BM-1   |
| <b>Preservative</b>                | 0.05% Sodium Azide   |
| <b>Isotype</b>                     | IgG1 Kappa   |
| <b>Conjugate</b>                   | mFluor Violet 610 SE   |
| <b>Purity</b>                      | Protein A or G purified  |
| <b>Buffer</b>                      | 50mM Sodium Borate   |
| <b>Product Description</b>         |  |
| <b>Host</b>                        | Mouse  |
| <b>Species</b>                     | Human  |
| <b>Marker</b>                      | Macrophage / Granulocyte Marker  |
| <b>Specificity/Sensitivity</b>     | Recognizes 183kDa protein with DNA-binding characteristics, which is identified as a myeloid specific antigen. It reacts with myeloid precursor cells and granulocytes in bone marrow. Its antigen appears to be restricted to M2 and M3 acute myelogenous leukemia (AML) subtypes. Markers of myeloid cells are useful in the identification of different levels of cellular differentiation. This monoclonal antibody reacts with early precursor and mature forms of human myeloid cells. It is useful in the identification of myelogenous leukemias, distinguishing granulocytic sarcomas from lymphoid malignancies and also in the study of differentiation and transformation of human myeloid cells. The biological function of this antigen is not clear, although it has been proposed that it may play a role in the differentiation of myeloid cells. |
| <b>Immunogen</b>                   | Human peripheral blood mononuclear cells   |
| <b>Notes</b>                       | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.  |
| <b>Product Application Details</b> |  |
| <b>Applications</b>                | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready  |
| <b>Recommended Dilutions</b>       | Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, CyTOF-ready  |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.   |



## Images

Myeloid Cell Marker Antibody (BM-1) [mFluor Violet 610 SE] [NBP2-34562MFV610] - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.



mFluor™ Violet 610

| LASER (nm)   | FILTER |
|--------------|--------|
| Violet (405) | 605/30 |

| EXCITATION MAX (nm) | EMISSION MAX (nm) |
|---------------------|-------------------|
| 421                 | 613               |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-34562MFV610](http://www.novusbio.com/reviews/submit/NBP2-34562MFV610)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

