

Product Datasheet

Human Breast Tumor Tissue MicroArray (Tumor) NBP2-42080

Unit Size: 5 Slides

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-42080

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-42080



NBP2-42080

Human Breast Tumor Tissue MicroArray (Tumor)

Product Information	
Unit Size	5 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C. Do not freeze.

Product Description	
Description	Please see online datasheet for well details: www.novusbio.com/NBP2-42080
Species	Human
Lysate Type	Tissue
Lysate Tissue	Breast
Lysate Tissue Condition	Tumor

Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin
Application Notes	Each individual tissue spot should be >90% core retention and total core lost from 1 slide should be no more than 5-10%. Please bake slides at 60C for 30 minutes before use.

Complete product information can be found online at www.novusbio.com/NBP2-42080.
Please contact technical service with any questions regarding the use of this product.

Publications

Niraj Jatin Patel, Sahaja Thippiani, Jasmine Jathan, Gauri Gaur, Janhavi Y. Sawant, Jay M. Pandya, Eva Sapi
Evidence for the presence of *Borrelia burgdorferi* in invasive breast cancer tissues European Journal of Microbiology & Immunology 2024-03-07 [PMID: 38451280]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Tissue Micro Arrays are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-42080

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

