## **Product Datasheet**

### CEACAM5/CD66e Antibody (C66/1292) [Alexa Fluor® 750] NBP2-54412AF750

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-54412AF750

Updated 10/23/2024 v.20.1

# Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-54412AF750



#### NBP2-54412AF750

Immunogen

| CEACAM5/CD66e Antibody (C66 | 5/1292) [Alexa Fluor® 750]   |
|-----------------------------|--|
| Product Information         |  |
| Unit Size                   | 0.1 ml   |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Store at 4C in the dark.   |
| Clonality                   | Monoclonal   |
| Clone                       | C66/1292   |
| Preservative                | 0.05% Sodium Azide   |
| Isotype                     | IgG2b  |
| Conjugate                   | Alexa Fluor 750  |
| Purity                      | Protein A or G purified  |
| Buffer                      | 50mM Sodium Borate   |
| Product Description         |  |
| Host                        | Mouse  |
| Gene ID                     | 1048   |
| Gene Symbol                 | CEACAM5  |
| Species                     | Human  |
| Specificity/Sensitivity     | This antibody recognizes proteins of 80-200kDa, identified as different members<br>of CEA family. CEA is synthesized during development in the fetal gut and is re-<br>expressed in increased amounts in intestinal carcinomas and several other<br>tumors. This monoclonal antibody does not react with nonspecific cross-reacting<br>antigen (NCA) and with human polymorphonuclear leucocytes. It shows no<br>reaction with a variety of normal tissues and is suitable for staining of<br>formalin/paraffin tissues. CEA is not found in benign glands, stroma, or malignant<br>prostatic cells. Antibody to CEA is useful in detecting early foci of gastric<br>carcinoma and in distinguishing pulmonary adenocarcinomas (60-70% are CEA<br>+) from pleural mesotheliomas (rarely or weakly CEA+). Anti-CEA positivity is<br>seen in adenocarcinomas from the lung, colon, stomach, esophagus, pancreas, |

www.novusbio.com



gallbladder, urachus, salivary gland, ovary, and endocervix.

Recombinant full-length human CEACAM5/CD66e protein (Uniprot: P06731)

|                              | Alexa Fluor (R) products are provided under an intellectual property license from<br>Life Technologies Corporation. The purchase of this product conveys to the<br>buyer the non-transferable right to use the purchased product and components<br>of the product only in research conducted by the buyer (whether the buyer is an<br>academic or for-profit entity). The sale of this product is expressly conditioned on<br>the buyer not using the product or its components, or any materials made using<br>the product or its components, in any activity to generate revenue, which may<br>include, but is not limited to use of the product or its components: (i) in<br>manufacturing; (ii) to provide a service, information, or data in return for<br>payment; (iii) for therapeutic, diagnostic or prophylactic purposes; or (iv) for<br>resale, regardless of whether they are resold for use in research. For information<br>on purchasing a license to this product for purposes other than as described<br>above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad,<br>CA 92008 USA or outlicensing@lifetech.com. This conjugate is made on<br>demand. Actual recovery may vary from the stated volume of this product. The<br>volume will be greater than or equal to the unit size stated on the datasheet. |
|------------------------------|---|
| Product Application Details  |   |
| Applications                 | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, CyTOF-ready  |
| <b>Recommended Dilutions</b> | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, CyTOF-ready  |
| Application Notes            | Optimal dilution of this antibody should be experimentally determined.  |
|                              |   |

Notes





#### Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

#### Products Related to NBP2-54412AF750

| IC0041S       | Mouse IgG2b Isotype Control (133303) [Alexa Fluor® 750] |
|---------------|---|
| NBP1-85742PEP | CEACAM5/CD66e Recombinant Protein Antigen               |
| DKK300        | Kallikrein 3/PSA [HRP]                                  |
| AVI10449-050  | CEACAM5/CD66e [Biotin]                                  |

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-54412AF750

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

