

Product Datasheet

C5b-9 Antibody (2A1)

NBP2-77377-0.1mg

Unit Size: 0.1 mg

Store at 4°C.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-77377

Updated 11/6/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-77377



NBP2-77377-0.1mg

C5b-9 Antibody (2A1)

Product Information

| | |
|----------------------|---|
| Unit Size | 0.1 mg |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4°C. |
| Clonality | Monoclonal |
| Clone | 2A1 |
| Preservative | No Preservative |
| Isotype | IgG1 |
| Purity | Protein G purified |
| Buffer | PBS |

Product Description

| | |
|------------------|-----------|
| Host | Mouse |
| Species | Rat |
| Immunogen | Rat C5b-9 |

Product Application Details

| | |
|------------------------------|---|
| Applications | Western Blot, Flow Cytometry, Immunoassay, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Immunoassay |
| Application Notes | Dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50. IA: plates were coated with 30 ug/ml in 50 mM carbonatebuffer, pH 10.6 for 16 h at 4 degrees C. WB: A non-reduced sample treatment and SDS-Page was used. The band size is >200 kDa. IHC-P: Tissue sections fixed in formalin were pretreated with protease type XXIV for 10 minutes at 37 degrees C before incubation. IHC-F: Tissue sections were fixed in acetone for 10 minutes at room temperature before incubation. |

Publications

Holste KG, Ye F, Koduri S et al. Attenuation of ventriculomegaly and iron overload after intraventricular hemorrhage by membrane attack complex inhibition Journal of neurosurgery 2023-11-10 [PMID: 37948699]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-77377

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

