# **Product Datasheet**

# Bone Marrow Tissue Slides (Diabetes)- Frozen NBP2-77712

Unit Size: 5 Slides

Store at -70C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-77712

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP2-77712



#### NBP2-77712

Bone Marrow Tissue Slides (Diabetes)- Frozen

<b>Product Information</b>	
Unit Size	5 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at -70C. Avoid freeze-thaw cycles.


Storage	Store at -70C. Avoid freeze-thaw cycles.	
Product Description		
Description	Each slide contains a single human adult tissue section, 5-10 um thickness that is mounted on a positively charged glass slide and fixed by cold acetone. The slides included in this package are adjacent/serial sections of bone marrow tissue from one donor. Tissue was snap frozen in liquid nitrogen immediately after excision and embedded in OCT. This product can be used for both immunohistochemistry and in-situ hybridization, especially for antibodies that don't recognize paraffin embedded tissue. At least one of the tissue slides from each lot was stained with H & E to ensure the quality (not included in the package).  Samples are IRB-approved from consented donors. Documentation on tissues' clinical histories may be available upon request. Donor information is also available upon request. Please contact nb-technical@bio-techne.com for any questions or requests.  Please contact nb-custom@bio-techne.com for all custom requests related to this product.	
Species	Human	
Notes	Donor information available upon request	
Lysate Type	Tissue	
Lysate Tissue	Bone/Cartilage	
Lysate Tissue Condition	Diabetes	

# **Product Application Details**

Applications	Immunohistochemistry, Immunohistochemistry-Frozen, In-situ Hybridization
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Frozen, In-situ Hybridization

## **Images**

Hematoxylin & Eosin Stain: Bone Marrow Tissue Slides (Diabetes) [NBP2-77712] - Tissue: Human Bone Marrow, Pathology: Diabetes







## Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### **Products Related to NBP2-77712**

NB820-59283

Human Bone Marrow Whole Tissue Lysate (Adult Whole Diabetes)

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-77712

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

