# **Product Datasheet**

# Uterus Tissue Slides (Endometrial Adenocarcinoma)-Paraffin NBP2-77941

Unit Size: 5 Slides

Store at 4C.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP2-77941

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-77941


#### NBP2-77941

Uterus Tissue Slides (Endometrial Adenocarcinoma)- Paraffin

	· · · · · · · · · · · · · · · · · · ·
Product Information	
Unit Size	5 Slides
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C.
Product Description	
Description	Each slide contains a single tissue section with 5 um thickness that is mounted on a positively charged glass slide. The slides included in this package are adjacent/serial sections of uterus tissue from one donor. Tissue was fixed in formalin immediately after excision and embedded in paraffin. This product can be used for both immunohistochemistry and in-situ hybridization. At least one of the tissue slides from each lot was stained with H&E to ensure the quality (not included in the package). Samples are IRB-approved from consented donors. Documentation on tissues' clinical histories may be available upon request. Donor information is also available upon request. Please contact nb-technical@bio-techne.com for any questions or requests. For FFPE DNA and RNA isolation, we may be able to provide tissue section slices in tubes. Please contact nb-custom@bio-techne.com for all custom requests related to this product.
Species	Human
Notes	Donor information available upon request
Lysate Type	Tissue
Lysate Tissue	Uterus
Lysate Tissue Condition	Endometrial Adenocarcinoma
Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, Dual RNAscope ISH-IHC
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin, In-situ Hybridization, Dual RNAscope ISH-IHC
Application Notes	Please bake slides at 60C for 30 minutes before use.

#### Images

Hematoxylin & Eosin Stain: Uterus Tissue Slides (Adenocarcinoma) [NBP2-77941] - Tissue: Human Uterus, Pathology: Endometrium Adenocarcinoma







## Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

#### **Bio-Techne Canada**

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

#### Products Related to NBP2-77941

NBP2-78095

Human Uterus Tissue MicroArray (Tumor)

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Slides are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-77941

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

