

Product Datasheet

Tangential Flow Filtration (TFF)-Large EV Isolation Adaptor NBP3-41045

Unit Size: 2 Units

Store at room temperature.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-41045

Updated 11/7/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-41045



NBP3-41045**Tangential Flow Filtration (TFF)-Large EV Isolation Adaptor**

Product Information	
Unit Size	2 Units
Concentration	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
Storage	Store at room temperature.
Product Description	
Description	<p>Large EV Isolation Adaptors for Tangential Flow Filtration (TFF)-Large EV Isolation (Catalog # NBP3-26838)</p> <p>Tangential flow filter-EV are a filter cartridge of hollow polyethersulfone fibers used for purification and concentration of extracellular vesicles. The filter allows a rapid, reproducible and scalable purification of EVs. It can be used on the lab bench for purifying small volume of samples or connected with a mechanical system for purifying larger volumes. During the purification, the sample is simultaneously concentrated.</p> <p>Luer lock connectors for TFF-EVs large filter. Connect the adaptors with the tangential flow filter (TFF) both sides and apply biomatrix into the TFF with large syringes or extend with pipes for peristaltic pump.</p>
Suitable Sample Type	EV containing matrices (urine, cell culture supernatant, plasma, serum, plant extracts)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-41045

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

