# **Product Datasheet**

# Profilin 1 Antibody (CL3524) - Azide and BSA Free NBP3-43957

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NBP3-43957

Updated 3/11/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NBP3-43957



#### NBP3-43957

Profilin 1 Antibody (CL3524) - Azide and BSA Free


Product Information	
Unit Size	100 ul
Concentration	LYOPH mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL3524
Reconstitution Instructions	Centrifuge the vial of lyophilized antibody at 12,000 x g for 20 seconds. Reconstitute by adding sterile, distilled water to achieve a final antibody concentration of 1mg/ml.
Isotype	lgG2a
Purity	Protein A purified
Buffer	Lyophilized from a 0.2 um filtered solution in PBS with Trehalose
Product Description	
Host	Mouse

•	
Host	Mouse
Gene ID	5216
Gene Symbol	PFN1
Species	Human
Immunogen	This antibody was generated using a recombinant protein sequence of P07737, with the exact immunogen sequence remaining proprietary.

<b>Product Application Details</b>	
Applications	Western Blot, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

# **Images**

Staining of human fallopian tube shows strong cytoplasmic positivity in glandular cells.



Page 2 of 4 v.20.1 Updated 3/11/2025 Staining of human liver shows strong cytoplasmic positivity in leukocytes. Staining of human testis shows very weak cytoplasmic positivity in cells in seminiferous ducts. Analysis in U-251MG cells transfected with control siRNA, target specific siRNA probe #1 and #2, using Anti-Profilin 1 antibody. Remaining relative intensity is presented. Loading control: Anti-GAPDH. Staining of human tonsil shows strong cytoplasmic positivity in germinal and non-germinal center cells.



Analysis in human cell line HeLa, human cell line HEK 293, human cell line A-431 and human cell line HepG2.





## Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

## Bio-Techne Canada 21 Canmotor Ave

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

#### **Bio-Techne Ltd**

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

#### **General Contact Information**

www.novusbio.com Technical Support: nb-technical@biotechne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

#### Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-43957

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

