

# ImmEdge Hydrophobic Barrier Pen

Catalog #310018

## Description

The ImmEdge Hydrophobic Barrier Pen provides a stable, water-repellent barrier that keeps reagents localized on the tissue specimens, preventing reagent mixing when multiple sections are mounted on the same slide.

## Applications

The ImmEdge Hydrophobic Barrier Pen is intended for the following laboratory applications:

- In situ hybridization
- Immunohistochemistry
- Immunocytochemistry
- Immunofluorescence

## Contents and Storage

The components of the ImmEdge Hydrophobic Barrier Pen and quantities supplied are listed in the following table. The product has a shelf life until the date of expiry indicated on the label. It is important to follow the product storage conditions as specified on the label so that the product maintains its reliability.

Product	Cat. No.	Component	Quantity	Storage
ImmEdge Hydrophobic Barrier Pen	310018	Hydrophobic barrier pen	1	RT

## Precautions

Handle biological specimens and other materials according to local laboratory guidelines and regulations. Avoid cross-contamination of samples and reagents to ensure quality results. Use equipment and reagents as specified by the manufacturer. Wear appropriate protective eyewear, clothing, and gloves. Read the Safety Data Sheets (SDS) and follow the handling instructions available at [https://www.bio-techne.com/p/in-situ-hybridization-ancillaries/immedge-hydrophobic-barrier-pen\\_310018](https://www.bio-techne.com/p/in-situ-hybridization-ancillaries/immedge-hydrophobic-barrier-pen_310018).

## Contact Information:

Advanced Cell Diagnostics  
 7707 Gateway Blvd. Ste #200, Newark CA 94560 USA  
 510.576.8800 1.877.576.3636  
[support.ACD@bio-techne.com](mailto:support.ACD@bio-techne.com)