

Data Sheet

Product Name: Paraffin Tissue Array

Catalog No.: NBP2-78114

Lot No.: A605038

Species: Human Mouse Rat Monkey (Rh) Guinea Pig Porcine
 Bovine Hamster Dog Monkey (Cy) Rabbit Plant

Tissue Type: Normal Adult Fetal Tumor Disease Cell line

Tissue Array Diagram:

	1	2	3	4	5	6	7	8	9	10	11
A	O	O	O	O	O	O	O	O	O	O	O
B	O	O	O	O	O	O	O	O	O	O	O
C	O	O	O	O	O	O	O	O	O	O	O
D	O	O	O	O	O	O	O	O	O	O	O
E	O	O	O	O	O	O	O	O	O	O	O
F	O	O	O	O	O	O	O	O	O	O	O

Components:

1. 5 Tissue Array Slides
2. Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____



Data Sheet

Product Name: Paraffin Tissue Array

Catalog No.: NBP2-78114

Lot No.: A605038

Donor Information

Position	Age	Sex	Organ	Pathological Diagnosis	Tumor History	Tumor Size (cm)	Differentiation	TNM	Tissue ID
A1	29	F	Positive Control	Human Normal Placenta	N/A	N/A	N/A		HT00031
A2			Negative Control	Paraffin					
A3	40	F	Ovary	Adenocarcinoma, mucous	5 months	7x5x5	Poorly	T3N2M1	HT00756
A4	60	F	Ovary	Adenocarcinoma	1 month	5x5x4	Poorly	T3N1M1	HT00757
A5	50	F	Ovary	Adenocarcinoma	1 month	5x5x4.5	Moderately	T3N0M0	HT00758
A6	63	F	Ovary	Cystadenocarcinoma, serous	1 month	6x5x5	Moderately	T3N1M0	HT00759
A7	55	F	Ovary	Clear Cell Adenocarcinoma	3 months+	diameter 13	Poorly	T3N0M0	HT00760
A8	39	F	Ovary	Matched Normal Tissue					HT00761
A9	39	F	Ovary	Serous papillary cystadenoma	2 months	15x14x9	Poorly	T4N5M1	HT00762
A10	51	F	Ovary	Matched Normal Tissue					HT00763
A11	51	F	Ovary	Thecoma	N/A	N/A	Moderately	TxNxMx	HT00764
B1	66	M	Bladder	Transitional Cell Carcinoma	1 month+	4x3x3	Moderately	T2N0M0	HT00765
B2	35	M	Bladder	Transitional Cell Carcinoma	6 years	N/A	Moderately	TxNxMx	HT00766
B3	61	M	Bladder	Transitional Cell Carcinoma	15 days	4x3.5x1.5	Poorly	T2N0M0	HT00767
B4	77	M	Bladder	Transitional Cell Carcinoma, Papillary	2 months	3x3x1.5	Well	T2N0M0	HT00768
B5	75	M	Bladder	Transitional Cell Carcinoma	2 years	6x5x2.5	Poorly	T3N0M0	HT00769
B6	69	M	Bladder	Transitional Cell Carcinoma, Papillary	2 months	8x8x1.5	Moderately	T3N0M0	HT00770
B7	69	M	Bladder	Transitional Cell Carcinoma, papillary	3 years	3x3x2	Moderately	T2N0M0	HT00771
B8	58	F	Bladder	Matched Normal Tissue					HT00772
B9	58	F	Bladder	Transitional cell carcinoma grade	3 months	4x4x2.5	Moderately	T2N0M0	HT00773
B10	20	M	Bladder	Matched Normal Tissue					HT00774
B11	20	M	Bladder	Transitional Cell Carcinoma	1 month	3x3x1.5	Well	T2N0M0	HT00775
C1	48	F	Breast	Invasive ductal carcinoam	1 year	Diameter 4	Moderately	T2N0M0	HT00776
C2	44	F	Breast	Invasive Lobular Carcinoma	4 months	6x4x4	Well	T2N17M1	HT00777
C3	82	F	Breast	Invasive ductal carcinoam	4 years	5x6	Moderately	T3N4M1	HT00778
C4	41	F	Breast	Invasive ductal carcinoam	2 months	1x1x0.8	Moderately	T2N1M0	HT00779
C5	78	F	Breast	Invasive ductal carcinoam	1 year	3x2	Moderately	T2N1M0	HT00780
C6	26	F	Breast	Fibroadenoma	2 months	1.8x1.5x0.7	N/A	T2N0M0	HT00781
C7	58	F	Breast	Invasive ductal carcinoma	2 months	6x5x4	Moderately	T3N3M0	HT00782
C8	56	F	Breast	Matched Normal Tissue					HT00783
C9	56	F	Breast	Invasive ductal carcinoma	1 year	5x2x3	Moderately	T3N0M0	HT00784
C10	36	F	Breast	Matched Normal Tissue					HT00785
C11	36	F	Breast	Invasive ductal carcinoma	2 months	2x2	Well	T2N0M0	HT00786
D1	50	M	Kidney	Clear Cell Carcinoma	40 days	3.5x3.5x3.2	Well	T2N0M0	HT00787

D2	70	M	Kidney	Clear Cell Carcinoma	7 days	5.5x5	Poorly	T3N0M0	HT00788
D3	64	M	Kidney	Clear Cell Carcinoma	20 days	4x3x3	Moderately	T2N0M0	HT00789
D4	48	M	Kidney	Clear Cell Carcinoma	30 days	3x3x2.5	Moderately	T2N0M0	HT00790
D5	50	M	Kidney	Clear Cell Carcinoma	20 days	4x3x3	Moderately-Poorly	T2N0M0	HT00791
D6	50	F	Kidney	Clear Cell Carcinoma	20 days	6.5x5x5	Moderately	T3N0M0	HT00792
D7	70	M	Kidney	Clear Cell Carcinoma	15 days	6x5x5	Moderately	T3N0M0	HT00793
D8	60	M	Kidney	Matched Normal Tissue					HT00794
D9	60	M	Kidney	Transational Cell Carcinoma	1 month	6x4x2.5	Moderately	T3N0M0	HT00795
D10	33	M	Kidney	Matched Normal Tissue					HT00796
D11	33	M	Kidney	Renal cell carcinoma	2 months	2.5x2x2	Moderately	T2N0M0	HT00797
E1	70	F	Uterus	Endometrium Adenocarcinoma, papillary	20 days+	3x2x1	Moderately	T2N0M0	HT00798
E2	49	F	Uterus	Endometrium Adenocarcinoma	2 years+	10x8x1.5	Moderately	T3N0M0	HT00799
E3	61	F	Uterus	Adenocarcinoma	2 months+	Diameter 0.3	Well	T1N2M0	HT00800
E4	25	F	Uterus	Squamous Cell Carcinoma	4 months+	5x2.5	Poorly	T3N0M0	HT00801
E5	36	F	Uterus	Adenomyoma of uterus	7 years	Diameter 7.5	N/A		HT00802
E6	42	F	Uterus	Adenoma	1 year+	Diameter 4	N/A		HT00803
E7	66	F	Uterus	Endometrium Adenocarcinoma	3 months+	3x0.4x0.2	Well	T2N0M0	HT00804
E8	52	F	Uterus	Matched Normal Tissue					HT00805
E9	52	F	Uterus	Leiomyoma	3 years	5x5x4	N/A		HT00806
E10	50	F	Uterus	Matched Normal Tissue					HT00807
E11	50	F	Uterus	Leiomyoma	6 months	3x2x2	N/A		HT00808
F1	58	M	Lung	Squamous Cell Carcinoma	3 months	5x4x4	Moderately	T3N2M0	HT00809
F2	46	M	Lung	Squamous Cell Carcinoma	45 days	6x5x5	Moderately	T3N2M0	HT00810
F3	63	M	Lung	Squamous Cell Carcinoma	40 days	3x3x2	Moderately	T2N3M0	HT00811
F4	62	M	Lung	Squamous Cell Carcinoma	60 days	5x4x4	Moderately-Poorly	T3N1M0	HT00812
F5	60	F	Lung	Adenocarcinoma	24 days	4x3.5x3.2	Poorly	T2N4M0	HT00813
F6	60	M	Lung	Squamous Cell Carcinoma	42 days	3.7x3x3	Moderately	T2N4M0	HT00814
F7	62	F	Lung	squamous cell carcinoma	20 days	3x2.5x2.2	Poorly	T2N1M0	HT00815
F8	50	F	Lung	Matched Normal Tissue					HT00816
F9	50	F	Lung	Adenocarcinoma	1 month	3x3x2.5	Moderately-Poorly	T2N5M1	HT00817
F10	63	M	Lung	Matched Normal Tissue					HT00818
F11	63	M	Lung	Bronchio alveolar carcinoma	15 days	7x5x5	Moderately	T3N4M1	HT00819

T - Primary tumor

Tx - Primary tumor cannot be assessed

T0 - No evidence of primary tumor

Tis - Carcinoma in situ; intraepithelial or invasion of lamina propria

T1 - Tumor invades submucosa

T2 - Tumor invades muscularis propria

T3 - Tumor invades through muscularis propria into subserosa or into non-peritonealized pericolic or perirectal tissues.

T4 - Tumor directly invades other organs or structures and/or perforate visceral peritoneum

N - Regional lymph nodes

Nx - Regional lymph nodes cannot be assessed

N0 - No regional lymph node metastasis

N1 - Metastasis in 1 to 3 regional lymph nodes

N2 - Metastasis in 4 or more regional lymph nodes

M - Distant metastasis

Mx - Distant metastasis cannot be assessed

M0 - No distant metastasis

M1 - Distant metastasis