

# Product Datasheet

## IL-4R alpha Antibody [Janelia Fluor® 669] AF-230-PBJF669

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF-230-PBJF669](http://www.novusbio.com/AF-230-PBJF669)

Updated 7/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF-230-PBJF669](http://www.novusbio.com/reviews/destination/AF-230-PBJF669)



**AF-230-PBJF669**

IL-4R alpha Antibody [Janelia Fluor® 669]

| Product Information         |                                                                                                                                                                                                                                                                                                               |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                        |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                   |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                      |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                    |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                            |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                           |
| Conjugate                   | Janelia Fluor 669                                                                                                                                                                                                                                                                                             |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                                                                                                                                     |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                            |
| Product Description         |                                                                                                                                                                                                                                                                                                               |
| Host                        | Goat                                                                                                                                                                                                                                                                                                          |
| Gene ID                     | 3566                                                                                                                                                                                                                                                                                                          |
| Gene Symbol                 | IL4R                                                                                                                                                                                                                                                                                                          |
| Species                     | Human                                                                                                                                                                                                                                                                                                         |
| Specificity/Sensitivity     | Detects human IL-4 R $\alpha$ in direct ELISAs and Western blots. In direct ELISAs and Western blots, less than 5% cross-reactivity with recombinant human (rh) IL-1 sRI, rhIL-1 sRII, rhIL-2 sR $\alpha$ , rhIL-2 sR $\beta$ , rhIL-2 sR $\gamma$ , rhIL-3 sR $\alpha$ , and rhIL-5 sR $\alpha$ is observed. |
| Immunogen                   | <i>S. frugiperda</i> insect ovarian cell line Sf 21-derived recombinant human IL-4 R $\alpha$                                                                                                                                                                                                                 |
| Notes                       | Sold under license from the Howard Hughes Medical Institute, Janelia Research Campus.                                                                                                                                                                                                                         |
| Product Application Details |                                                                                                                                                                                                                                                                                                               |
| Applications                | Western Blot, Flow Cytometry, CyTOF-ready                                                                                                                                                                                                                                                                     |
| Recommended Dilutions       | Western Blot, Flow Cytometry, CyTOF-ready                                                                                                                                                                                                                                                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                        |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to AF-230-PBJF669**

---

|                 |                                           |
|-----------------|-------------------------------------------|
| NBP2-61210-50ug | Recombinant Human IL-4R alpha His Protein |
| 210-TA-005      | TNF-alpha [Unconjugated]                  |
| 230-4RB-025     | IL-4R alpha [Unconjugated]                |
| M6000B-1        | IL-6 [HRP]                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF-230-PBJF669](http://www.novusbio.com/reviews/submit/AF-230-PBJF669)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

