

Product Datasheet

TWEAK/TNFSF12 Antibody [DyLight 650] AF1237W

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF1237W

Updated 7/13/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF1237W



AF1237W**TWEAK/TNFSF12 Antibody [DyLight 650]**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 650
Purity	Antigen Affinity-purified
Buffer	50mM Sodium Borate
Product Description	
Host	Goat
Gene ID	8742
Gene Symbol	Tnfsf12
Species	Mouse
Specificity/Sensitivity	Detects TWEAK/TNFSF12 in direct ELISAs and Western blots. In direct ELISAs, approximately 30% cross-reactivity with recombinant human TWEAK is observed and less than 1% cross-reactivity with recombinant mouse (rm) BAFF, rmFas Ligand, rmOX40 Ligand, rmTNF- α , and rmTRAIL is observed.
Immunogen	<i>E. coli</i> -derived recombinant mouse TWEAK/TNFSF12 Arg105-His249 Accession # O54907
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, Neutralization
Recommended Dilutions	Western Blot, Neutralization
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF1237W

NB410-28088C-500ul	Goat IgG Isotype Control [DyLight 650]
NBP2-46811PEP	TWEAK/TNFSF12 Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
1090-TW-025/CF	TWEAK/TNFSF12 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF1237W

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

