

# Product Datasheet

## Cystatin C Antibody [HRP] AF1238H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1238H](http://www.novusbio.com/AF1238H)

Updated 7/7/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1238H](http://www.novusbio.com/reviews/destination/AF1238H)



**AF1238H**

Cystatin C Antibody [HRP]

| Product Information         |                                                                                                                                                                        |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                                               |
| Clonality                   | Polyclonal                                                                                                                                                             |
| Preservative                | No Preservative                                                                                                                                                        |
| Isotype                     | IgG                                                                                                                                                                    |
| Conjugate                   | HRP                                                                                                                                                                    |
| Purity                      | Antigen Affinity-purified                                                                                                                                              |
| Buffer                      | PBS                                                                                                                                                                    |
| Product Description         |                                                                                                                                                                        |
| Host                        | Goat                                                                                                                                                                   |
| Gene ID                     | 1471                                                                                                                                                                   |
| Gene Symbol                 | CST3                                                                                                                                                                   |
| Species                     | Mouse                                                                                                                                                                  |
| Specificity/Sensitivity     | Detects mouse and human Cystatin C in direct ELISAs. In direct ELISAs, less than 10% cross-reactivity with recombinant human Cystatin D is observed.                   |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant mouse Cystatin C (R & D Systems, Catalog # 1238-PI)<br>Ala21-Ala140<br>Accession # P21460                              |
| Product Application Details |                                                                                                                                                                        |
| Applications                | Western Blot, Immunohistochemistry, Immunoprecipitation, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair) |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunoprecipitation, ELISA Standard (Matched Pair), ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair) |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to AF1238H**

---

|                |                                      |
|----------------|--------------------------------------|
| AB-108-CH      | Goat IgG Isotype Control [HRP]       |
| NBP2-61381-1mg | Recombinant Human Cystatin C Protein |
| 7954-GM-010/CF | GM-CSF [Unconjugated]                |
| MSCTC0         | Cystatin C [HRP]                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1238H](http://www.novusbio.com/reviews/submit/AF1238H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

