

# Product Datasheet

## Semaphorin 7A Antibody [DyLight 405] AF1835E

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1835E](http://www.novusbio.com/AF1835E)

Updated 7/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1835E](http://www.novusbio.com/reviews/destination/AF1835E)



**AF1835E****Semaphorin 7A Antibody [DyLight 405]**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                                                         |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                  |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                             |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                                                                                                                                                              |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                      |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                     |
| <b>Conjugate</b>                   | DyLight 405                                                                                                                                                                                                                                                                                                                             |
| <b>Purity</b>                      | Antigen Affinity-purified                                                                                                                                                                                                                                                                                                               |
| <b>Buffer</b>                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                      |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                                                         |
| <b>Host</b>                        | Goat                                                                                                                                                                                                                                                                                                                                    |
| <b>Gene ID</b>                     | 8482                                                                                                                                                                                                                                                                                                                                    |
| <b>Gene Symbol</b>                 | SEMA7A                                                                                                                                                                                                                                                                                                                                  |
| <b>Species</b>                     | Mouse                                                                                                                                                                                                                                                                                                                                   |
| <b>Specificity/Sensitivity</b>     | Detects mouse Semaphorin 7A in direct ELISAs and Western blots. In these formats, less than 1% cross-reactivity with recombinant mouse (rm) Semaphorin 3C, rmSemaphorin 3E, rmSemaphorin 3F, rmSemaphorin 6A, rmSemaphorin 6B, rmSemaphorin 6C, recombinant human (rh) Semaphorin 3A, rhSemaphorin 3B, and rhSemaphorin 6D is observed. |
| <b>Immunogen</b>                   | Mouse myeloma cell line NS0-derived recombinant mouse Semaphorin 7A<br>Gln45-Ala646<br>Accession # Q9QUR8                                                                                                                                                                                                                               |
| <b>Notes</b>                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                       |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                                                         |
| <b>Applications</b>                | Western Blot                                                                                                                                                                                                                                                                                                                            |
| <b>Recommended Dilutions</b>       | Western Blot                                                                                                                                                                                                                                                                                                                            |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to AF1835E**

---

|                 |                                                            |
|-----------------|------------------------------------------------------------|
| NB410-28088V    | Goat IgG Isotype Control [DyLight 405]                     |
| NBP2-61977-10ug | Recombinant Human Semaphorin 7A Isoform 1 hlgG-His Protein |
| 233-FB-025      | FGF basic/FGF2/bFGF [Unconjugated]                         |
| 2068-S7-050     | Semaphorin 7A [Unconjugated]                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1835E](http://www.novusbio.com/reviews/submit/AF1835E)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

