

# Product Datasheet

## JNK2 Antibody [mFluor Violet 610 SE] AF1846MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1846MFV610](http://www.novusbio.com/AF1846MFV610)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1846MFV610](http://www.novusbio.com/reviews/destination/AF1846MFV610)



**AF1846MFV610**

JNK2 Antibody [mFluor Violet 610 SE]

| Product Information         |                                                                                                                                                                                                                                                                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                       |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                          |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                        |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                               |
| Conjugate                   | mFluor Violet 610 SE                                                                                                                                                                                                                                                                                                              |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                                                                                                                                                         |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                |
| Product Description         |                                                                                                                                                                                                                                                                                                                                   |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                                            |
| Gene ID                     | 5601                                                                                                                                                                                                                                                                                                                              |
| Gene Symbol                 | MAPK9                                                                                                                                                                                                                                                                                                                             |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                 |
| Specificity/Sensitivity     | Detects human, mouse, and rat JNK2 in Western blots. Recognizes p54 JNK (isoforms 1 and/or 4, both 424 aa) and p46 JNK (isoforms 2 and/or 3, both 382 aa) in Western blots. Because the peptide immunogen corresponds to a region of JNK2 with low homology to other JNKs, the antibody does not detect recombinant JNK1 or JNK3. |
| Immunogen                   | Human, Mouse, and Rat JNK2 synthetic peptide<br>PGIKFEELFPDWIFPSESERD                                                                                                                                                                                                                                                             |
| Notes                       | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.                                                                                             |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                   |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF1846MFV610**

---

|                  |                                                   |
|------------------|---------------------------------------------------|
| NBP2-24891MFV610 | Rabbit IgG Isotype Control [mFluor Violet 610 SE] |
| NBP1-98883-0.1mg | Recombinant Human JNK2 His Protein                |
| 210-TA-005       | TNF-alpha [Unconjugated]                          |
| DYC2236-2        | JNK2 [Biotin]                                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1846MFV610](http://www.novusbio.com/reviews/submit/AF1846MFV610)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

