

# Product Datasheet

## IL-2 Antibody [DyLight 350] AF1890D

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1890D](http://www.novusbio.com/AF1890D)

Updated 7/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1890D](http://www.novusbio.com/reviews/destination/AF1890D)



**AF1890D**

IL-2 Antibody [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                          |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Host                        | Goat                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Gene ID                     | 3558                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Gene Symbol                 | IL2                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Species                     | Feline                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Specificity/Sensitivity     | Detects feline IL-2 in ELISAs and Western blots. In sandwich immunoassays, approximately 3% cross-reactivity with recombinant canine IL-2 is observed, approximately 1% cross reactivity with recombinant human IL-2 and recombinant porcine IL-2 is observed, and less than 0.2% cross-reactivity with recombinant mouse IL-2, recombinant equine IL-2, recombinant bovine IL-2, recombinant rat IL-2, and recombinant cotton rat IL-2 is observed. |
| Immunogen                   | <i>E. coli</i> -derived recombinant feline IL-2; Ala21-Thr154 (Cys146Ser)<br>Accession # Q07885                                                                                                                                                                                                                                                                                                                                                      |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                                                                                                                                    |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Applications                | Western Blot, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), ELISA Standard (Matched Pair), Immunocytochemistry, Neutralization                                                                                                                                                                                                                                                                                     |
| Recommended Dilutions       | Western Blot, ELISA Standard (Matched Pair), Neutralization, ELISA Capture (Matched Antibody Pair), ELISA Detection (Matched Antibody Pair), Immunocytochemistry                                                                                                                                                                                                                                                                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                               |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to AF1890D**

---

|                 |                                        |
|-----------------|----------------------------------------|
| NB410-28088UV   | Goat IgG Isotype Control [DyLight 350] |
| NBP2-35054-10ug | Recombinant Human IL-2 Protein         |
| 210-TA-005      | TNF-alpha [Unconjugated]               |
| 202-IL-010      | IL-2 [Unconjugated]                    |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1890D](http://www.novusbio.com/reviews/submit/AF1890D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

