

# Product Datasheet

## ICAM-5 Antibody [DyLight 650] AF1950W

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1950W](http://www.novusbio.com/AF1950W)

Updated 7/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1950W](http://www.novusbio.com/reviews/destination/AF1950W)



**AF1950W**

ICAM-5 Antibody [DyLight 650]

| Product Information         |                                                                                                                                                                                                                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                             |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                        |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                           |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                         |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                 |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                |
| Conjugate                   | DyLight 650                                                                                                                                                                                                                                                                        |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                                                                                                          |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                 |
| Product Description         |                                                                                                                                                                                                                                                                                    |
| Host                        | Goat                                                                                                                                                                                                                                                                               |
| Gene ID                     | 7087                                                                                                                                                                                                                                                                               |
| Gene Symbol                 | ICAM5                                                                                                                                                                                                                                                                              |
| Species                     | Human                                                                                                                                                                                                                                                                              |
| Specificity/Sensitivity     | Detects human ICAM-5 in direct ELISAs and Western blots. In direct ELISAs and Western blots, approximately 10% cross-reactivity with recombinant mouse ICAM-5 is observed and less than 2% cross-reactivity with recombinant human (rh) ICAM-1, rhICAM-2 and rhICAM-3 is observed. |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human ICAM-5<br>Ala28-Glu570<br>Accession # Q9UMF0                                                                                                                                                                                 |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                  |
| Product Application Details |                                                                                                                                                                                                                                                                                    |
| Applications                | Western Blot                                                                                                                                                                                                                                                                       |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                                                                       |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                             |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to AF1950W**

---

|                    |                                           |
|--------------------|-------------------------------------------|
| NB410-28088C-500ul | Goat IgG Isotype Control [DyLight 650]    |
| NBP3-05774-50ug    | Recombinant Human ICAM-5 hIgG-His Protein |
| 7268-CT-100        | CTLA-4 [Unconjugated]                     |
| 1173-M5-050        | ICAM-5 [Unconjugated]                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1950W](http://www.novusbio.com/reviews/submit/AF1950W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

