

# Product Datasheet

## IL-12/IL-23 p40 Antibody [DyLight 488] AF1969K

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1969K](http://www.novusbio.com/AF1969K)

Updated 7/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1969K](http://www.novusbio.com/reviews/destination/AF1969K)



**AF1969K**

IL-12/IL-23 p40 Antibody [DyLight 488]

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 ml  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.               |
| Storage                     | Store at 4C in the dark.  |
| Clonality                   | Polyclonal  |
| Preservative                | 0.05% Sodium Azide  |
| Isotype                     | IgG   |
| Conjugate                   | DyLight 488   |
| Purity                      | Antigen Affinity-purified   |
| Buffer                      | 50mM Sodium Borate  |
| Product Description         |   |
| Host                        | Goat  |
| Gene ID                     | 3593  |
| Gene Symbol                 | IL12B   |
| Species                     | Canine  |
| Specificity/Sensitivity     | Detects canine IL-12/IL-23 p40 in direct ELISAs and Western blots.  |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant canine IL-12/IL-23 p40 Ile23-Ser329<br>Accession # Q28268 |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                         |
| Product Application Details |   |
| Applications                | Western Blot, Neutralization  |
| Recommended Dilutions       | Western Blot, Neutralization  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF1969K**

---

|              |  |
|--------------|--|
| NB410-28088G | Goat IgG Isotype Control [DyLight 488] |
| 210-TA-005   | TNF-alpha [Unconjugated]               |
| DY499-05     | IL-12/IL-23 p40 [Biotin]               |
| M6000B-1     | IL-6 [HRP]                             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1969K](http://www.novusbio.com/reviews/submit/AF1969K)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

