

# Product Datasheet

## NOV/CCN3 Antibody [DyLight 488] AF1976K

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1976K](http://www.novusbio.com/AF1976K)

Updated 7/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1976K](http://www.novusbio.com/reviews/destination/AF1976K)



**AF1976K**

NOV/CCN3 Antibody [DyLight 488]

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Conjugate</b>     | DyLight 488                                                                                 |
| <b>Purity</b>        | Antigen Affinity-purified                                                                   |
| <b>Buffer</b>        | 50mM Sodium Borate                                                                          |

**Product Description**

|                                |                                                                                                                                                                                                                                     |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                    | Goat                                                                                                                                                                                                                                |
| <b>Gene ID</b>                 | 4856                                                                                                                                                                                                                                |
| <b>Gene Symbol</b>             | NOV                                                                                                                                                                                                                                 |
| <b>Species</b>                 | Mouse                                                                                                                                                                                                                               |
| <b>Specificity/Sensitivity</b> | Detects mouse NOV/CCN3 in direct ELISAs and Western blots. In Western blots, approximately 20% cross-reactivity with recombinant human NOV is observed and less than 1% cross-reactivity with recombinant mouse WISP-1 is observed. |
| <b>Immunogen</b>               | Mouse myeloma cell line NS0-derived recombinant mouse NOV/CCN3 Ser26-Ile354<br>Accession # Q64299                                                                                                                                   |
| <b>Notes</b>                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                   |

**Product Application Details**

|                              |                                                                        |
|------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot                                                           |
| <b>Recommended Dilutions</b> | Western Blot                                                           |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF1976K**

---

|                |                                        |
|----------------|----------------------------------------|
| NB410-28088G   | Goat IgG Isotype Control [DyLight 488] |
| NBP2-35084-5ug | Recombinant Human NOV/CCN3 Protein     |
| NBP2-07996     | HIST1H4F Overexpression Lysate         |
| 1640-NV-050    | NOV/CCN3 [Unconjugated]                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1976K](http://www.novusbio.com/reviews/submit/AF1976K)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

