

# Product Datasheet

## Reg3B Antibody [DyLight 350] AF1996D

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF1996D](http://www.novusbio.com/AF1996D)

Updated 7/13/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF1996D](http://www.novusbio.com/reviews/destination/AF1996D)



**AF1996D**

Reg3B Antibody [DyLight 350]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 350
Purity	Antigen Affinity-purified
Buffer	50mM Sodium Borate
Product Description	
Host	Goat
Gene ID	18489
Gene Symbol	Reg3b
Species	Rat
Specificity/Sensitivity	Detects rat Reg3B in direct ELISAs and Western blots. In direct ELISAs, approximately 75% cross-reactivity with recombinant mouse (rm) Reg3B, and approximately 20% cross-reactivity with rmReg3G is observed, and less than 2% cross-reactivity with recombinant human (rh) Reg3A, rmReg1, rmReg2, rhReg4, rhReg1A, rmReg3A, and rmReg3D is observed. In Western blots, approximately 40% cross-reactivity with rmReg3B is observed.
Immunogen	<i>E. coli</i> -derived recombinant rat Reg3B Glu27-Gly175 Accession # P25031
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF1996D**

---

NB410-28088UV	Goat IgG Isotype Control [DyLight 350]
DKK300	Kallikrein 3/PSA [HRP]
5110-RG-050	Reg3B [Unconjugated]
NBP1-19773	Serpin E1/PAI-1 Antibody - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF1996D](http://www.novusbio.com/reviews/submit/AF1996D)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

