

Product Datasheet

TNF RI/TNFRSF1A Antibody [mFluor Violet 450 SE] AF225MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF225MFV450

Updated 7/7/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF225MFV450



AF225MFV450

TNF RI/TNFRSF1A Antibody [mFluor Violet 450 SE]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	mFluor Violet 450 SE
Purity	Antigen Affinity-purified
Buffer	50mM Sodium Borate
Product Description	
Host	Goat
Gene ID	7132
Gene Symbol	TNFRSF1A
Species	Human
Specificity/Sensitivity	Detects human sTNF RI/TNFRSF1A in direct ELISAs and Western blots. In direct ELISAs, less than 20% cross-reactivity with recombinant mouse TNF sRI and recombinant rat TNF sRI is observed. This anti-human TNF RI antibody can interact with the cell surface 55 kDa human TNF RI and exhibits TNF agonist activities on the A549 human lung carcinoma cell line. This antibody also exhibits TNF agonist activities on the L $\square$ 929 mouse fibroblast cell line, suggesting that it has cross-species activity with mouse cell surface TNF RI.
Immunogen	<i>E. coli</i> -derived recombinant human sTNF RI/TNFRSF1A Accession # NP_001056
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, Flow Cytometry, Agonist Activity, CyTOF-ready, Immunocytochemistry
Recommended Dilutions	Western Blot, Flow Cytometry, Agonist Activity, Immunocytochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF225MFV450

NB410-28088MFV450	Goat IgG Isotype Control [mFluor Violet 450 SE]
NBP2-35274-5ug	Recombinant Human TNF RI/TNFRSF1A Protein
210-TA-005	TNF-alpha [Unconjugated]
DRT100	TNF RI/TNFRSF1A [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF225MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

