

# Product Datasheet

## TFPI-2 Antibody [PE] AF2545P

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF2545P](http://www.novusbio.com/AF2545P)

Updated 5/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF2545P](http://www.novusbio.com/reviews/destination/AF2545P)



**AF2545P**

TFPI-2 Antibody [PE]

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C in the dark.                                                                    |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.05% Sodium Azide                                                                          |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Conjugate</b>     | PE                                                                                          |
| <b>Purity</b>        | Antigen Affinity-purified                                                                   |
| <b>Buffer</b>        | PBS                                                                                         |

**Product Description**

|                                |                                                                                                                                                                                                                                                                              |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                    | Goat                                                                                                                                                                                                                                                                         |
| <b>Gene ID</b>                 | 7980                                                                                                                                                                                                                                                                         |
| <b>Gene Symbol</b>             | TFPI2                                                                                                                                                                                                                                                                        |
| <b>Species</b>                 | Human                                                                                                                                                                                                                                                                        |
| <b>Specificity/Sensitivity</b> | Detects human TFPI $\square$ 2 in direct ELISAs and Western blots. In direct ELISAs, approximately 10% cross $\square$ reactivity with recombinant mouse TFPI-2 is observed and less than 5% cross $\square$ reactivity with recombinant human TFPI $\square$ 1 is observed. |
| <b>Immunogen</b>               | Mouse myeloma cell line NS0-derived recombinant human TFPI-2<br>Asp23-Lys213<br>Accession # P48307                                                                                                                                                                           |

**Product Application Details**

|                              |                                                                        |
|------------------------------|------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Immunoprecipitation, Neutralization                      |
| <b>Recommended Dilutions</b> | Western Blot, Immunoprecipitation, Neutralization                      |
| <b>Application Notes</b>     | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to AF2545P**

---

|               |                                    |
|---------------|------------------------------------|
| AB-108-CP     | Goat IgG Isotype Control [PE]      |
| NBP2-33552PEP | TFPI-2 Recombinant Protein Antigen |
| 2914-HT-100MG | Holo-Transferrin [Unconjugated]    |
| 2545-PI-010   | TFPI-2 [Unconjugated]              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF2545P](http://www.novusbio.com/reviews/submit/AF2545P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

