

# Product Datasheet

## Vasohibin Antibody [DyLight 680] AF3319Y

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF3319Y](http://www.novusbio.com/AF3319Y)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF3319Y](http://www.novusbio.com/reviews/destination/AF3319Y)



**AF3319Y**

Vasohibin Antibody [DyLight 680]

| Product Information         |                                                                                                             |
|-----------------------------|-------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                 |
| Storage                     | Store at 4C in the dark.                                                                                    |
| Clonality                   | Polyclonal                                                                                                  |
| Preservative                | 0.05% Sodium Azide                                                                                          |
| Isotype                     | IgG                                                                                                         |
| Conjugate                   | DyLight 680                                                                                                 |
| Purity                      | Antigen Affinity-purified                                                                                   |
| Buffer                      | 50mM Sodium Borate                                                                                          |
| Product Description         |                                                                                                             |
| Host                        | Goat                                                                                                        |
| Gene ID                     | 22846                                                                                                       |
| Gene Symbol                 | VASH1                                                                                                       |
| Species                     | Human                                                                                                       |
| Specificity/Sensitivity     | Detects human Vasohibin Isoform 1 in direct ELISAs and Western blots.                                       |
| Immunogen                   | Mouse myeloma cell line NS0-derived recombinant human Vasohibin Isoform 1 Pro2-Val365<br>Accession # Q7L8A9 |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                           |
| Product Application Details |                                                                                                             |
| Applications                | Western Blot                                                                                                |
| Recommended Dilutions       | Western Blot                                                                                                |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF3319Y**

---

|                    |                                                  |
|--------------------|--------------------------------------------------|
| NB410-28088FR      | Goat IgG Isotype Control [DyLight 680]           |
| H00022846-P01-10ug | Recombinant Human Vasohibin GST (N-Term) Protein |
| DVE00              | VEGF [HRP]                                       |
| NBL1-17697         | Vasohibin Overexpression Lysate                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF3319Y](http://www.novusbio.com/reviews/submit/AF3319Y)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

