

# Product Datasheet

## TSC22/TSC22D1 Antibody [DyLight 650] AF3779W

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF3779W](http://www.novusbio.com/AF3779W)

Updated 5/8/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF3779W](http://www.novusbio.com/reviews/destination/AF3779W)



**AF3779W**

TSC22/TSC22D1 Antibody [DyLight 650]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Polyclonal                                                                                  |
| Preservative                | 0.05% Sodium Azide                                                                          |
| Isotype                     | IgG                                                                                         |
| Conjugate                   | DyLight 650                                                                                 |
| Purity                      | Antigen Affinity-purified                                                                   |
| Buffer                      | 50mM Sodium Borate                                                                          |
| Product Description         |                                                                                             |
| Host                        | Goat                                                                                        |
| Gene ID                     | 8848                                                                                        |
| Gene Symbol                 | TSC22D1                                                                                     |
| Species                     | Mouse, Rat                                                                                  |
| Specificity/Sensitivity     | Detects mouse and rat TSC22 in direct ELISAs and Western blots.                             |
| Immunogen                   | <i>E. coli</i> -derived recombinant mouse TSC22<br>Met1-Ala143<br>Accession # P62500.2      |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.           |
| Product Application Details |                                                                                             |
| Applications                | Western Blot                                                                                |
| Recommended Dilutions       | Western Blot                                                                                |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF3779W**

|                    |                                                      |
|--------------------|------------------------------------------------------|
| NB410-28088C-500ul | Goat IgG Isotype Control [DyLight 650]               |
| H00008848-P02-10ug | Recombinant Human TSC22/TSC22D1 GST (N-Term) Protein |
| DKK300             | Kallikrein 3/PSA [HRP]                               |
| H00008848-P01-10ug | Recombinant Human TSC22/TSC22D1 GST (N-Term) Protein |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF3779W](http://www.novusbio.com/reviews/submit/AF3779W)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

