

Product Datasheet

LAMP-1/CD107a Antibody [FITC] AF4320F

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF4320F

Updated 5/8/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF4320F



AF4320F**LAMP-1/CD107a Antibody [FITC]**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	FITC
Purity	Antigen Affinity-purified
Buffer	PBS
Product Description	
Host	Goat
Gene ID	3916
Gene Symbol	LAMP1
Species	Mouse
Specificity/Sensitivity	Detects mouse LAMP1/CD107a Luminal Domain in direct ELISAs and Western blots. In direct ELISAs, less than 10% cross-reactivity with recombinant human LAMP1 is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse LAMP1/CD107a Leu25-Asn370 Accession # P11438
Product Application Details	
Applications	Western Blot, CyTOF-ready, Immunocytochemistry, Intracellular Staining by Flow Cytometry
Recommended Dilutions	Western Blot, Intracellular Staining by Flow Cytometry, Immunocytochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF4320F

NB410-28088F-500ul	Goat IgG Isotype Control [FITC]
NBP2-59780PEP	LAMP-1/CD107a Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
4320-LM-050	LAMP-1/CD107a

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF4320F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

