

Product Datasheet

TSC1 Antibody [DyLight 550] AF4379L

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/AF4379L

Updated 5/8/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/AF4379L



AF4379L

TSC1 Antibody [DyLight 550]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 550
Purity	Antigen Affinity-purified
Buffer	50mM Sodium Borate
Product Description	
Host	Sheep
Gene ID	7248
Gene Symbol	TSC1
Species	Human, Mouse, Rat
Specificity/Sensitivity	Detects human, mouse, and rat TSC1 in direct ELISAs and Western blots.
Immunogen	<i>E. coli</i> -derived recombinant human TSC1 Asp156-Thr300 Accession # Q92574
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, CyTOF-ready, Intracellular Staining by Flow Cytometry
Recommended Dilutions	Western Blot, Intracellular Staining by Flow Cytometry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to AF4379L

H00007248-Q01-10ug	Recombinant Human TSC1 GST (N-Term) Protein
NBL1-17350	TSC1 Overexpression Lysate
AF887	AKT [p Ser473] Antibody [Unconjugated] - Pan Specific
NB200-157	4EBP1 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/AF4379L

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

