

# Product Datasheet

## SIRP gamma/CD172g Antibody [DyLight 594] AF4486M

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF4486M](http://www.novusbio.com/AF4486M)

Updated 5/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF4486M](http://www.novusbio.com/reviews/destination/AF4486M)



**AF4486M**

SIRP gamma/CD172g Antibody [DyLight 594]

| Product Information         |                                                                                                                                                                                                                                        |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                               |
| Clonality                   | Polyclonal                                                                                                                                                                                                                             |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                     |
| Isotype                     | IgG                                                                                                                                                                                                                                    |
| Conjugate                   | DyLight 594                                                                                                                                                                                                                            |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                                                              |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                     |
| Product Description         |                                                                                                                                                                                                                                        |
| Host                        | Sheep                                                                                                                                                                                                                                  |
| Gene ID                     | 55423                                                                                                                                                                                                                                  |
| Gene Symbol                 | SIRPG                                                                                                                                                                                                                                  |
| Species                     | Human                                                                                                                                                                                                                                  |
| Specificity/Sensitivity     | Detects human SIRP $\gamma$ /CD172g in direct ELISAs and Western blots. In direct ELISAs, less than 5% cross-reactivity with recombinant human (rh) SIRP $\alpha$ , rhSIRP $\beta$ 1, and recombinant mouse SIRP $\alpha$ is observed. |
| Immunogen                   | Chinese hamster ovary cell line CHO-derived recombinant human SIRP $\gamma$ /CD172g Val64-Ser364, (Val263Ala) and (Ser286Leu) Accession # Q9P1W8                                                                                       |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                      |
| Product Application Details |                                                                                                                                                                                                                                        |
| Applications                | Western Blot, Flow Cytometry, CyTOF-ready                                                                                                                                                                                              |
| Recommended Dilutions       | Western Blot, Flow Cytometry, CyTOF-ready                                                                                                                                                                                              |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF4486M**

---

|             |                                         |
|-------------|-----------------------------------------|
| 7268-CT-100 | CTLA-4 [Unconjugated]                   |
| 9999-SB-050 | SIRP gamma/CD172g [Unconjugated]        |
| NBP2-22203  | ERK1 Antibody (1E5)                     |
| NB100-309   | ATM Antibody (2C1) - Azide and BSA Free |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF4486M](http://www.novusbio.com/reviews/submit/AF4486M)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

