

# Product Datasheet

## VEGF Antibody [FITC] AF564F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF564F](http://www.novusbio.com/AF564F)

Updated 7/14/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF564F](http://www.novusbio.com/reviews/destination/AF564F)



**AF564F****VEGF Antibody [FITC]**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                                                                                             |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                 |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                    |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                          |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Conjugate</b>                   | FITC                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>                      | Antigen Affinity-purified                                                                                                                                                                                                                                                                                                                                                   |
| <b>Buffer</b>                      | PBS                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Host</b>                        | Goat                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Gene ID</b>                     | 7422                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Gene Symbol</b>                 | VEGFA                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Species</b>                     | Rat                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Specificity/Sensitivity</b>     | Detects rat VEGF in direct ELISAs and Western blots. In direct ELISAs, approximately 20% cross-reactivity with recombinant human (rh) VEGF165 and rhVEGF121 is observed and less than 2% cross-reactivity with rhVEGF-B, recombinant mouse (rm) VEGF-B, rhVEGF-C, rhVEGF-D, and rmVEGF-D is observed. In Western blots, detection of recombinant mouse VEGF165 is observed. |
| <b>Immunogen</b>                   | Mouse myeloma cell line NS0-derived recombinant rat VEGF164 Ala27-Arg190<br>Accession # AAL07526                                                                                                                                                                                                                                                                            |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Applications</b>                | Western Blot, Immunohistochemistry, Neutralization                                                                                                                                                                                                                                                                                                                          |
| <b>Recommended Dilutions</b>       | Western Blot, Immunohistochemistry, Neutralization                                                                                                                                                                                                                                                                                                                          |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF564F**

---

|                    |                                 |
|--------------------|---------------------------------|
| NB410-28088F-500ul | Goat IgG Isotype Control [FITC] |
| NB100-2381PEP      | VEGF Antibody Blocking Peptide  |
| 210-TA-005         | TNF-alpha [Unconjugated]        |
| DVE00              | VEGF [HRP]                      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF564F](http://www.novusbio.com/reviews/submit/AF564F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

