

# Product Datasheet

## **Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 Antibody [mFluor Violet 450 SE] AF6664MFV450**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/AF6664MFV450](http://www.novusbio.com/AF6664MFV450)

Updated 5/10/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/AF6664MFV450](http://www.novusbio.com/reviews/destination/AF6664MFV450)



**AF6664MFV450**

Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 Antibody [mFluor Violet 450 SE]

| Product Information         |                                                                                                                                                                                                                                       |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                           |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                              |
| Clonality                   | Polyclonal                                                                                                                                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                    |
| Isotype                     | IgG                                                                                                                                                                                                                                   |
| Conjugate                   | mFluor Violet 450 SE                                                                                                                                                                                                                  |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                                                             |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                    |
| Product Description         |                                                                                                                                                                                                                                       |
| Host                        | Sheep                                                                                                                                                                                                                                 |
| Gene ID                     | 11041                                                                                                                                                                                                                                 |
| Gene Symbol                 | B4GAT1                                                                                                                                                                                                                                |
| Species                     | Human, Mouse                                                                                                                                                                                                                          |
| Specificity/Sensitivity     | Detects human and mouse $\beta$ -1,4-Glucuronyltransferase 1/B4GAT1 in Western blot. In direct ELISAs, less than 1% cross-reactivity with recombinant human (rh) B4GalT1 and rhB3GNT 2 is observed.                                   |
| Immunogen                   | Chinese hamster ovary cell line CHO-derived recombinant human $\beta$ -1,4-Glucuronyltransferase 1/B4GAT1<br>Asp43-Cys415<br>Accession # O43505                                                                                       |
| Notes                       | mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Product Application Details |                                                                                                                                                                                                                                       |
| Applications                | Western Blot                                                                                                                                                                                                                          |
| Recommended Dilutions       | Western Blot                                                                                                                                                                                                                          |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to AF6664MFV450**

---

|               |                                                                               |
|---------------|-------------------------------------------------------------------------------|
| NBP1-88952PEP | Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 Recombinant Protein Antigen |
| MEP00B        | Erythropoietin/EPO [HRP]                                                      |
| 6664-GT-050   | Beta-1,3-N-Acetylglucosaminyltransferase 1/B3GNT1 [Unconjugated]              |
| NB100-524     | NOD2 Antibody (2D9) - BSA Free                                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/AF6664MFV450](http://www.novusbio.com/reviews/submit/AF6664MFV450)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

